














LÂMPADAS ESPECIAIS PARA EQUIPAMENTO HOSPITALAR

CÓDIGO DO ARTIGO	DESCRIÇÃO	FIGURA
ML8010 Zeiss SLIT LAMP 12V 50W LT03104, 12V 50W. INAMI microscópio Lâmpada de fenda. Wattage: 50W, 50 Watt, 50 Watts Voltage: 12V, 12 Volt, 12 Volts Glass: T-3 Tubular Clear	12v-50w, 12v 4, 2A, 12v-4,2A Zeiss SLIT LAMP 12V 50W P44S. Microlites 12-050-52M, Inami L0169-H1, 12v 50w Base: Special E14/26x18 Flange Average Life: 50 Hours Filament: C-BAR6 Lumens: 1400, Gas: Krypton	
ML8011 71-61-86 Baush and Lomb Slit Lamp bulb for model OSL- 2000, OSL-3000. Mentor Slit Lamp bulb for model G10, GH12, S10, SH12. Burton SL-825, Mentor S-10, G-10, SH-10, GH-12, 22-4010, Woodlyn HR-1, HR-3, HR-5	ML8011 Woodlyn Slit Lamp bulb 12V- 50w, 12V 4.2 A, 12V 50 W European Screw with Special Flange base. Length: 66mm. 12V 50w. 45070	
ML8012 6v- 20w, 6v 20w P44S	ML8012 6v-20w P44S 6v 20w P44S	
ML8014 12v-20w, 12v 20w P44S	ML8014 12v-20w P44S 12v 20w P44S	
ML8051 Lamp for Veatch E200 & E500. Replacement light bulb for Veatch E200 and E500 slit lamps 12V 30W Bulb and silver base attachment sold included.	Bulb For Veatch E200 & E500 Slit Lamp, Slit Lamp. SLVE20-500. 12V 30W, 12v-30w. Replacement light bulb for Veatch E200 and E500 slit lamps. 12V-30W Bulb and silver base attachment sold included.	
ML9901 SLMA03-500, Bulb, Slit Lamp Halogen H-S Style. 6 Volts 20 Watts. 6V 20W, 6V-20W	Slit Lamp Halogen , Fits Haag-Streit 900BM (Incandescent style) Topcon SL-3D, Topcon SL-3E (Incandescent style), Marco IIB, Marco IV, Marco V, Marco VG, Burton L1000 6 Volts 20 Watts Height 60 mm	
ML9022 1019, 6W 20W Marco IIB- Photo/III/V/VG/VG-Photo Halogen 6v/20w.1019	1019, 6W-20W, Slit lamp bulb for the following slit lamps: Marco IIB, IIB-R, IIB-Photo, III, V, V-R, V- Photo, VG, VG-Photo	
ML7999 Lamps Zeiss Model FF450 + IR. Zeiss FF 450 Plus IR. Zeiss FF-450 Plus IR. Camera to image with a Merge/OIS WS5000 system. Carl Zeiss FF 450 Not available any more	Lamps Zeiss Model FF450 + IR. Carl Zeiss FF 450 IR FUNDUS CAMERA / RETINAL CAMERA Not available any more	
ML7997 Zeiss Fixation Bulb for Zeiss FF450 IR. Medical lamps, Medical bulbs. Led Lamp 301350-9052-000, Light sources: Medical - Microscope	ML4180 3,5V Spare Bulb Retinoscope Zumax Lamp 3,5V Retinoscópico Zumax Zumax. lamps R4180-LED Bulb	






<p>ML2284 Haag-Streit HS940-750 6V 4,5A, 940-750 (OP 2510) (OP2510) 6V 27W TAKAGI. (Takagi 330-PF)(Takagi MT- 325UD 8708) 6v-4,35A P40S Haag Streit Goldman Perimeter 6v 4,35A P40S</p>		
<p>ML2202 TAKAGI, HAAG STREIT 6V 4,5A 900-930, 6V 27W P44S (OP2366S) (SM900-930) 3000353, (OR 1900) 366 900930, MT 303, MT-303, MT 304, ZITEM-129</p> <p>ML2203 TAKAGI, HAAG STREIT 6V 4,5A P44S HS900-930 , SL-250 HY 366 (MT-303), Nidek SL250, SL-450, Nidek SL450, Eiko 41318</p>	<p>ML7794 Topcon SL-3F, 6V 4,5A (40310-25600) SL-6E, SL-3F. Lamp for Shin Nippon SL-101. Burton Medical P4892 • Marco 1008 • Marco 1008-G • Topcon 40310-25600 • Ushio 8000373 SM-40310-25600. Shin-Nippon / Rexxam SL-101 / SL- 102 / SL-80. Topcon Models SL-3 / SL-3D / SL-5D / SL-6E Dr. Fisher 00843096, 6V 27W</p>	 <p>ML2203 ML2202</p>
<p>ML9221 Lamp 9v 4,5A SPEC. P44S Lamp 9v/4,5A SPEC. P44S Lamp 9v-4,5A SPEC. P44S</p>		
<p>ML9224 ROYAL 4710 10V 4A Lamp 9v-4,5A (ROYAL 4710) SPEC. P39S, Lamp. TOC 10V-4A 10V/4A, 64x35mm, Diâmetro da base 39mm, Com 6 furos de 6,3mm</p>	<p>ROYAL 4710, 10V 4A P39S, Lamp 9v-4,5A (ROYAL 4710) SPEC. P39S, Lamp. TOC 10V-4A 10V 4A Lamp Measures 64x35mm, Base diameter 39mm, With 6 holes diameter 6,3mm each</p>	
<p>ML9102 Haag Streit BM/BQ/BX-900 Slit Lamp Main Bulb 7,5v 38w, 7.5v-38w [22-0120-0960]</p>	<p>0900-0960 (Haag-Streit Slit Lamp 900 BM, BQ, BP, BX Halogen 7,5V-38W Haag-Streit 900 BM/BQ/BP/BX Halogen 7,5v/38w Discontinued, See ML2202</p>	
<p>ML5066 DCA 104512 , 21v 150w 21v-150w</p>		
<p>ML6609 SYLVANIA DEF bulb A1/211 21.5V 150W Lamp REF: DEFSY 21.5V-150W AVG 10HRS</p>		
<p>ML4711 Lamp. A1/211 Atlas L1, 21,5V 150W (A1-211 ATLAS, Lamp. A1/211 Atlas L1 (DEF/21.5V/150W ANSI) 21,5V-150W</p>		
<p>ML6610 DFC 125V 150W G17Q ANSI DFC 125 v-150 w. KODAK BROWNIE 8-A-15 Projector. DNF 125V 150W G17Q ANSI</p>		




<p>ML8001 TOPCON 40340-20700 / JC 6V 20W, 6V-20W, CL70536 SL 3G Ref. 4474020400 SL3G 4474020400 6v 20w TOPCON 40340-20700 6v 20w. Lamp 5628/3</p>	<p>6V 20W TOPCON 40340-20700 / JC 6V 20W - Medical Lamp. USHIO SM40300-20700, Hosobuchi-HF9102, 7507/H, Topcon 40340-20700, Topcon SL2E, Topcon SL-2E, Topcon SL7F, SL 7F 6v-20w, Hikari M-01005. Microscópio Sbsisa. Topcon SL-D301</p>	
<p>ML8002 Topcon Slit SL-7F, Topcon Slitlamp 6v-20w Falange Keeler 1030-L-2050; Keeler Slit Lamp 1030-L-2050; Topcon Slit Lamp SL-1E; Topcon SL-7F; (Topcon SL1 - SL2 - SL3 - SL4) (M-01005). HF9102, Topcon 447032530 (Topcon SL-1 - SL-2 - SL-3 - SL-4)</p>	<p>6V 20W Topcon 4474020400 6v 20w Slit Lamp. CARLEY 1778 8Z D7, TOPCON-44680-25700, CARLEY CL1778, USHIO SM40300-20700, Hosobuchi-HF9102, TOPCON SLIT-7F. Topcon 40340-20400, Topcon SL2E, Topcon SL-2E, Topcon SL7F, SL 7F 6v 20w, Hikari M-01005. 6V/20W. Topcon Slit SL 7F, Topcon SL-D301, Topcon SL-1 /Topcon SL-2 /Topcon SL-3 /Topcon SL4</p>	 <p>Base: 25x25mm (ranhura 4,5x1,5mm)</p>
<p>ML9002 Huvitz slit 6V 20W Huvitz CCP-3100P. Replacement bulb for Huvitz CCP-3100 auto projector</p>	<p>6V 20W M-01208 6V 20W Huvitz CCP 3100 Chart Projection Lamp. 6V 20W lamp for Huvitz CCP 3100. 6V 20W lamp for Huvitz CCP-3100</p>	
<p>ML4291 Huvitz slit 12V 30W Halogen bulb. Huvitz slit lamp HS-7000 12V 30W. Huvitz HS 5000/Huvitz HS 5500/Huvitz HS 7000 Slit Lamp Bulb</p>	<p>12V 30W 12V 6A 30W Huvitz HS-7000 Lamp 12V 6A 30W halogen bulb Huvitz slit lamp. Huvitz HS-7000 Huvitz HS-5000/ Huvitz HS-5500/ Huvitz HS7000 Slit Lamp Bulb</p>	
<p>ML8003 Topcon SL-D7/8Z 12v 30w slit lamp special base. Topcon SL-D7/8Z Slit Lamp Halogen. LT03067. Topcon SL-8Z, SL-D8Z, SL-D7. 4468025700 (12V 30W)</p>	<p>12V 30W. SL-D7/8Z/D8Z 12V 30W Topcon 12V30W 42440-20030 Medical Dental Ophthalmic. LT03067, Topcon: 40340-20770. SL-D701 Topcon SL-D7, Topcon SL-8Z</p>	 <p>Base: 25x20mm (ranhura 4,5x1,5mm)</p>
<p>ML8008 TOPCON-42412-20400. TOPCON 42412-20400 ACP-8R. CL1780. 5628/5 12V-50W BASE. ACP-8R, ACP-8 12V 50W BASE. ACP-8R Topcon ACP-8</p>	<p>12V 50W. TOPCON 42412-20400LS. TP 42412-20400LS OPHTHALMIC KEYWAY AND MOUNTING HOLE IN BRASS BASE Replacement Bulbs/Lamps for TOPCON "ACP8,42412-20400"; TP 42412-20400LS;OPHTHALMIC;KEYWAY AND MOUNTING HOLE IN BRASS</p>	
<p>ML8009 Topcon ACP-8 and ACP-1000 Chart Projector ACP 8 12v 50w Halogen lamp (ACP-1000). LT03068. Topcon 42412-20400 ACP-8 Chart Projector ACP-1000 12V 50W 4,17 A Flansch CODE: LT03068. 12V/50W</p>	<p>12v 50w Topcon ACP-8, Auto Chart Projector, Replacement Light Bulb 12v-50w Special base, G6.35pins. Cross Reference: Topcon slit lamp code:LT03068 volt:12 Watt:50V Base: Special, Pins:2-pins,G6.35 Life: 100 hrs. (ACP 8, ACP 1000). Focus Pine 12V 50W Special - Topcon 42412-2040 ACP-8</p>	
<p>ML8119 Chart Projector ACP-1000 12v 50w Halogen lamp (ACP-1000). LT03068, 42412-20400. Chart Projector ACP-1000 12V 50W Flansch</p>	<p>12v 50w Topcon ACP-1000, Auto Chart Projector, Replacement Light Bulb 12v-50w Special base, G6.35pins. Cross Reference: Topcon slit lamp code:LT03068 volt:12 Watt:50V Base: Special, Pins:2-pins,G6.35 Life:100 hrs</p>	
<p>ML8004 Topcon SL-12V, Lâmpada Fenda SL 12 V, Lâmpada Fenda SL12V, Lâmpada Fenda SL 12V, SL/12V Lamp SL12V</p>	<p>Topcon SL12V, Topcon SL 12V, Topcon SL-12V, Lâmpada de Fenda Topcon SL 12V Lamp Topcon SL 12 V Lamp SL 12 V, Lamp SL12 V</p>	






<p>ML7422 Hosobuchi 0-3120 6V 1,5A BA15 (6v 9W BA15) (37,5x18,5mm) Hosobuchi 1702</p>	<p>Hosobuchi 0-3120 6V 1,5A (37,5x18,5mm) With connecting wires, 3mm holes</p>	
<p>XBO R 100W/45C XBO R 101W/45C XBO R 180W/45C Osram XBO R 100W/45C OFR Osram XBO R 101W/45C OFR Osram XBO R 180W/45C OFR Pentax EPM-1000 Pentax EPM1000</p>	<p>(XBO R 180W/45C for Zeiss Opmi Sensera) (We have also with 4 pole plug) ML4405 Osram XBO R 100W/45C OFR Pentax EPK-100P/EPK-1000 With 5 Pole Plug</p>	
<p>ML4404 XBO R 100W/45C OFR ozone free xenon bulb for endoscopic light sources With 4 pole Plug Com 4 fichas de ligação Pentax EPM-1000, Pentax EPM1000 Osram XBO R 100W/45C</p>	<p>ML4527 Roche LightCycler 480 Xenon Lamp Cat. No. 04 686 136 001. Xenon Lamp 04 686 136 001. Light Cycler 480 Xenon Lamp Cat. No. 04686136 001. Osram Xenon Reflector Lamp 100 watt ozone free xenon bulb for endoscopic light sources</p>	 <p style="text-align: center;">ML4404 ML4527</p>
<p>Osram HBO R 103 W/45 Mercury Vapor Light Bulb, 100 watt 23 volt Short Arc Double Ended HBO Mercury HID Osram Light Bulb. Sylvania 69311 (HBO R 103 W/45) Average Life: 300 hours. (SL69311). HBO-R-103- W-45-23V/69311</p>	<p>A lâmpada HBO R 103W destina- se especialmente a aplicações para condução de luz e possui um reflector elíptico que permite a ligação directa da luz num condutor. HBOR103W45 23v 100w, 23v-100w</p>	
<p>OSRAM XBO R 180W/45C OSRAM 180/45C XBO-R 14v 180W. The XBO R 180 W/45C OFR is a 180 watt, ozone free quartz Xenon Short Arc Lamp with focussing 44 degree beam angle. Typical lifetime is 500 hours. This lamp has lead wires and connector on the 67mm diameter x 90mm long dichroic coated elliptical reflector with a 45mm working distance.</p>	<p>ML9997 XBO R 180W/45C OFR. XBOR180W/45C, XBO-R180W/45C OFR, Product Warning: Xenon Short Arc lamps emit harmful ULTRAVIOLET RADIATION, are under HIGH PRESSURE, require a HIGH VOLTAGE ignition, and can GENERATE OZONE (unless designated "OFR".) For Zeiss OPMI LUMERA 700 and RESCAN 700 (XBO R 180W/45C for Zeiss Opmi Sensera)</p>	
<p>OSRAM XBO R 181W/45C Xenon lamp 181w Wire Leads</p>	<p>ML9988 XBO R 181W/45C OFR XENON SHORT ARC 14VDC 69184 OSRAM OZONE- FREE COOLING REQUIRED 500 HOUR</p>	
<p>Osram XBO R 300W/60C XBO R 300W/60C OFR, LS7601 300 Watt Replacement Bulb with Mounting Bracket. HID/XBO/300W/XBOR300-60COFR. OSXBOR300W60COFR, OS XBO R300W/60COFR #69167, Sylvania 69167. 20V 300W, Average Lifetime: 1000 hours For Microscope Zeiss Opmi Pentero. XBOR300</p>	<p>Osram XBO R 300W/60C OFR 81524, 611129, 857749 04613569167, 4008321137722 4050300611129, 4050300857749 CL70079 CARLEY-LAMPS HID/XBO/300W/XBOR300-60COFR Osram XBO-R-300W/60C/OFR HID/XBOR/300W/XBOR300-60COFR X-R-300WATT-60-CABLE-OFR. 300 watt 20 volt T17. Sylvania 69167 - XBO R 300 W 60C OFR Healthcare Medical Scientific Light Bulb</p>	
<p>ML7133 OSRAM HXP R 120W/45C HXP R 120W/45C VIS OSRAM HXP R 120W/45C VIS, Product description: HXP R 120W/45C VIS, HXP-R 120W-45C VIS, EAN/ Product: 4050300882772, Reflector diameter, max. 67 mm,</p>	<p>OSRAM HXP R 120W/45C VIS Technical - Electrical Data Type of current AC Construction voltage 85V Wattage 120W Lamp voltage 60...90 V - initial voltage range, Construction current 1.4 A Lamp current 1.6 A Rated wattage 120W, Diameter 67.0 mm Length with base excl. base pins/connect 77 mm Cable length 200 mm Lifespan 2000 h. OSRAM HXP R 120W/45X VIS</p>	



<p>STRYKER QUANTUM 3000 Osram HTI 250W/32 846598/32V/250W, 846598/250W, 846598, 846604/250W 8488 220130300 OSRAM 0220130300 STRYKER Osram HTI 250W/32</p>	<p>HID/MH-HTI/270W/MH250-CONNECTOR- MR21DSPECIAL-32MMWD-5600K M56-038-011 STRYKER MH250-CONNECTOR-MR21DSPECIAL- 32MMWD-5600K-STRYKER-HTI QUANTUM-3000/STRYKER INTERLIGHT QUANTUM-3000 STRYKER</p>	
<p>STRYKER QUANTUM 400 QUANTUM 400 OSHTI250W22S5, STRYKER QUANTUM 4000 FIBER OPTIC LIGHT SOURCE STRYKER 220- 180-300. SURGICAL / QUANTUM 5000, STRYKER With 4 pole Plug</p>	<p>OSRAM HTI 250W/22C STRYKER QUANTUM 400 Quantum 400 HTI 250-22C, Description: 270W MH 45V NO CONNECTOR 54080 OSRAM SYLVANIA SHORT ARC LAMP HTI250W/22 This is a Bulb Matrix Value Replacement. (ALB-UPC: 04613554080). 45v 270w With 4 pole Plug</p>	
<p>STRYKER QUANTUM 4000 STRYKER QUANTUM 5000 FIBER OPTIC LIGHT SOURCE STRYKER 220-180-300 STRYKER QUANTUM 5000 STRYKER QUANTUM 400 Osram HTI 250/22C LAMPADA DE FONTE LUZ STRYKER DO TROLLEY DE LAPAROSCOPIA Refª 220-180-300 Q5000 BULB</p>	<p>STRYKER QUANTUM 4000 Description: 270W MH 45V NO CONNECTOR 54080 OSRAM SYLVANIA SHORT ARC LAMP HTI250W/22 This is a BulbMatrix Value Replacement. (ALB-UPC: 04613554080). Osram HTI 250/22C 45v 270w. C/suporte (Completo) MODULE OSHTI250W22S5</p>	
<p>Osram HTI 250W/32C Osram HTI250/32C Light Source Storz Light Source Richard Wolf 5140 Endoscope Light Source With 5 pole Plug</p>	<p>OSRAM HTI 250W/22C Stryker Q5000 endoscopy lamp METAL HALIDE 45V 6A Endoscope Light Source With 5 pole Plug</p>	
<p>Osram HTI 400W/24 Osram HTI 403W/24 Osram HTI 404W/24</p>	<p>55V 24W. Wattage: 400W, 400 Watt Watts Voltage: 55V, 55 Volt Volts Amperage: 7.3A, 300Hz Square-Wave AC, Base: Special Plug. Bulb Type: PDR-67 Electrodes: Tungsten Rods with 4.5mm</p>	
<p>ML5502 Conmed / Linvatec light source LS-7500 Only Bulb ML5511 Conmed / Linvatec light source LS-7500 Complete Module</p>	<p>Lâmpada para Fonte de Luz Lamp For Conmed LS-7500, Conmed LS 7500 Conmed LS7500 Chip 100115D</p>	
<p>ML7098A 300W, Pentero, 1277-2200 Complete Module Zeiss Reference 1277-2200 (no Exchange) Microscope Zeiss Pentero OPMI Pentero Multivision neu 390</p>	<p>ML7098 300W, Pentero, 1294-6580 Complete Module Zeiss Ref1294-658, ZE1294658 (Exchange) Microscope Zeiss Pentero OPMI Pentero Multivision neu 390</p>	
<p>ML8496 Replacement lamp for all models of Carl Zeiss HXP 200 fluorescence microscope illuminators. For models:435716 - Illuminator HXP 200 C</p>	<p>Carl Zeiss HXP 200 fluorescence microscope illuminators Lamp HT-B8916 for Zeiss HXP200 and 435716</p>	
<p>ML9491 Lamp Osram P-VIP 240/0.8 E20.9n, 76V 240W FaP2.5. P-VIP 240/0.8 E 301 Osram P-VIP 240/0.8 E20.9n Osram Projector lamp</p>	<p>76,0 V 240 W Soquete FaP2.5. 4008321806390, 4052899247765, 4052899514935, 4052899555280, 4052899599307, 4062172040914, 4062172354714, 4062172354851, 4062172374408, 4062172393874, 4062172393812, 4062172393775, 4062172393898, 4062172395366, 4062172396622</p>	



<p>ML5808 Lamps Philips RP502P Philips 35W</p>	<p>Lamps Philips RP502P 35W, Xenon lamp 35W used in ATMOS, Philips RP 50/2-P Philips Philips RP50-2P Philips</p>	
<p>ML2882 3255 Narva 25.1606/22, 24v 100w BA20D (Narva 25.1606722) Narva 24v 100w 25.1606/22 .</p>	<p>3255 Narva 251606/22 24v-100w BA20D .Diameter 60mm x Long 85mm. Life: 100h Antique Lamp. Vintage Lamp</p>	
<p>ML8829 Ohmeda Giraffe socket 6600-1176-400 Incubadora Ohmeda Giraffe</p>	<p>ML1991 Agglutinoscope Merz+dade AG, Lamp 229V 25W E14/25 (230V 25W E14, 20x11, 10x301501S, 142011649, T20x115mm clear ROHS. aglutinoscópio Merz+dade AG lamp</p>	
<p>ML2887 Narva 25.1606/22 24v 100w BA20D Narva 25.1606722)</p>	<p>Ref. 8.102.77 Compleat box with green glass filter and lamp Narva 25.1606722</p>	
<p>ML7406 Vidros azuis com 28cm de Diâmetro ML7404 Vidros Brancos Martelados Com 28cm de diâmetro</p>	<p>ML7406 Blue Glass (Stop ray Red) Diameter 280mm for lamps Hanaulux ML7404 White Glass, Diameter 280mm For Lamps Hanaulux</p>	
<p>ML7402 Vidro Verde (Stop ray Red) c/diâmetro de 16,5cm e altura de 20cm para foco SIGLA</p>	<p>Green Glass (Stop ray Red) 165x200mm for SIGLA</p>	
<p>ML4480 P-129352-157 (P129352157) For Amasco Examiner 10 25V 155W BA15D 24v 150w BA15D Frosted</p>	<p>P-129352-157 25V 155W (FOSCA) FROSTED FOCO MARCA AMASCO EXAMINER 10, 25V-155W 24v 150w, 24V-150W BA15D Frosted</p>	
<p>ML4481 Lamps for Burton Visionary 24v 40w BA15D, 25v 35w BA15D. 844409004418 HIKARI, JC-24V-40W-BA15D-M01045</p>		
<p>ML4482 24V 50W BA15D, 24V 50W BA15D H/24V/50W/BA15D/T3.5/31MMLCL</p>		
<p>ML4483 24v 150w BA15D Lamp for Surgical DHT 24V/150W/BA15D/130540 NARVA</p>		
<p>ML4484 24v 100w BA15D, 24v 100w Ba15d, 24v-100w BA15D, 5429/2B</p>		
<p>ML3378 ESR 120v 100w BA15D ESR 120v-100w BA15D</p>		
<p>ML3379 ESP 120v 150w BA15D, ESP 120v-150w BA15D, 15246, 44384 EXAMINATION-LIGHT/49455 CASTLE/MDT</p>		
<p>ML3380 ESS 120v 250w BA15D ESS 120v-250w BA15D</p>		



<p>ML7488 F4T5/BL-B 7G Tube Black Light bulb. Fluorescent , Luz Negra 6 INCH 4W, F4T5/BLB, HR4T5BLB</p>	<p>ML7489P (Philips) Philips TUV 4Watt G4T5 ultraviolet Tube / Bulb, TUVG4T5, Germicidal lamps G4 T5 (TUV4W) 16x149mm Philips 63872427 Philips 301219T ML7489 (Osram) Osram Puritec Germicidal Lamp HNS 4W G5 G4T5/OF. 16x149mm</p>	 <p>ML7489 ML7489P</p>
<p>ML5112 Fuji Medical Part Number 123N0007D lamps Fuji photo film tube 123N0007D</p>	<p>Fuji Medical 123N0007D Fuji (length 480mm diameter 15mm) for device digitalisator Fuji photo film Model: CR-IR 346RU S/N: 36526210B replacement erase lamps for Fuji computed radiography readers Original Medical Fuji 123N0007D</p>	
<p>ML7454 G36T5/4P/SE 39W UV Output (watts): 12W G36T54P/SE, G36T54PSE Germicidal lamps use a special glass that transmits 254nm energy that inactivates bacteria, viruses, molds, and cysts by preventing replication</p>	<p>Fluorescent Lamp, T5 Lamp Shape, 4- Pin (G10q) Base, Length 34 In., 24 Watts, Lamp Designation G36T5/4P/SE, Clear, UVC source, 5/8 In. Dia., Average Life 9K hr., Cool Lamp Tone, Rapid, Standards Meets 2012 Lighting Legislation Energy Standards. Only shipped in quantities of 24.</p>	
<p>ML4717 Lamp TUV PL-L 95W/4P HO 1CT/25 TUV PL-L 95W/4P HO ref. 88829740</p>	<p>Lamp TUV PL-L 95W/4P HO 1CT/25, Lamp 88829740 Emitem radiação UV de onda curta com um pico de 353.7 nm(UV-C) para acção germicida; -O vidro da lâmpada filtra a linha de formação de ozônio de 185 nm</p>	
<p>Waldmann F85/100W-PUVA F85/100W PUVA, F85/100W- UV21, F85/100W-TL01. PUVA 100, 1000, 3001, 3002, 3003, 6001, 6002, 701, UV 1000KL, UV 1000K, UV 3002K, UV 3003K, UV 7001K, UV 8001K, UV 7002, UV 5040 KL, UV 5001 KL, UV 1000 KL</p>	<p>Indications: -Psoriasis, atopic dermatitis -Vitiligo, Lichen ruber planus -Mycosis Fungoides (MF) -Scleroderma, acne vulgaris -Uraemic pruritus -Polymorphous light dermatosis -Others Phototherapy of Psoriasis</p>	
<p>TL100W/01, TL 100W 01 Dermatological phototherapy lamp dim:1760MMx38MM medical lamp nm:300 nm:315 phototherapy lamp Therapy lamp ultra-voilet-B/01. Emission: between 305 to 315 nm with a peak of 311 nm</p>	<p>Lâmpada dermatologia fototerapia Dim: 1760MMx38MM Medical lamp, phototherapy lamp therapy lampultra-voilet-B/01. (F/72INCH/T12/R17D/100W/300N) TL100W 01. 007873/122V/100W/R17D, 007873 (639745),8711500614254 PHILIPS F/69INCH/T12/R17D/100W/310NM/122V F/70INCH/T12/R17D/100W/UVB-01 F/72INCH/T12/R17D/100W/01 F/72INCH/T12/R17D/100W/01-UVB</p>	 <p>R17d recessed double contact</p> <p>T12 Tubular 12 / 8 inch 1-1/2 inch 38mm diameter T38mm</p>
<p>Electronic Starters N61347-2-1,4-125W Waldmann OKO6 - the most popular starter to work in temperatures from -25 for C to 100 of C, the fluorescent 4-125W, OKO2 - for serial connections, for temperatures of -25 with C to 100 of C, for 4-22W fluorescent lamps and compact fluorescent lamps and wheeled the dim power range. Electronic.</p>	<p>Arranadores electrónicos FRIGOSTART - to work in temperaturachod -40 for C to 100 of C, for 4-125W fluorescent lamps (used in freezers) TURBO - for applications where an instant start fluorescent lamps (-5 for C to 100 of C) for fluorescent 4-65W SOLARSTART - to tanning salons (4-180W). Source: company materials S12 Sunstarter, 115w - 140w Tubes Cleo Body Tone Starte 100-180W Cosmostarter/E (15w-225w Sunbeds)</p>	



<p>Waldmann UVA7001K, UV7001K Phototherapy for psoriasis Ultraviolet light therapy for atopic dermatitis. 97612638 UV lamps for device UV 7001K, UVA 7001K, UVA 7001K, UV 7001K, UV 7001 K, UV7001 K Lamp. dim:1760MMx38MM Phototherapy lamp, therapy lamp ultra-voilet-B/01. TL100W/01, F85/100W-UV01 (F85/100W-TL01) UVA GOLD 100W (TL85/100W/09) (WALDMANN F85/100W-PUVA), TL100W-10R. When replacing the fluorescent lamp it is recommended to also replace the starter Waldmann electronic ballast</p>	<p>Combined therapy: 27 pcs of F85/100W-PUVA +13 pcs of F85/100W-TL01 or 27 pcs of F85/100W-PUVA +13 pcs of F85/100W-UV21 or 27 pcs of F85/100W-PUVA +13 pcs of TL100W-UV01 or 20 pcs of F85/100W-PUVA +20 pcs of F85/100W-TL01 or 20 pcs of F85/100W-PUVA +20 pcs of F85/100W-UV6 Waldmann Starter (for 100W tubes) Waldmann ballast (for 100W tubes)</p>	
<p>Waldmann UVA7001K, UV7001K Ultraviolet light therapy for atopic dermatitis. 97612638 UV lamps for device UV7001K, UVA 7001K, UVA 7001K, UV 7001K, UV 7001 K, UV7001 K Lamp. dim:1760MMx38MM Phototherapy lamp, therapy lamp ultra-voilet-B/01. TL100W/01, F85/100W-UV01. UVA GOLD 100W (TL85/100W/09) F85/100W-01, (WALDMANN F85/100W-PUVA), TL100W-10R.) Therapy-System. When replacing the fluorescent lamp it is recommended to also replace the starter Waldmann electronic ballast</p>	<p>Monotherapy: 40 pcs of F85/100W-PUVA or 40 pcs of F85/100W-TL01 or 40 pcs of F85/100W-UV6 or 40 pcs of F85/100W-TL10R or 40 pcs of F85/100W-UV21 or 40 pcs of F85/100W-UV01 or 40 pcs of TL 100W-01 Waldmann Starter (for 100W tubes) Waldmann ballast (for 100W tubes)</p>	
<p>Waldmann UV 7002 UV 7002, UVA 7002, UVA7002 - Therapy System. Radiation (UVA): Lamp waldmann: F79/120W-PUVA 320-410nm. Radiation (UVB): Lamp waldmann: F79/120W-TL01 310-315nm. When replacing the fluorescent lamp it is recommended to also replace the starter Waldmann electronic ballast Waldmann UV 7002. Waldmann Therapy-System</p>	<p>Combined therapy: 21 pcs F79/120W- PUVA + 21 pcs F79/120W -TL01 Monotherapy: 42 pcs F79/120W - PUVA or 42 pcs F79/120W - TL01 Waldmann Starter (for 120W tubes) Waldmann ballast (for 120W tubes)</p>	
<p>Waldmann UV 5040 KL (UV 5040KL) UV5040KL UV-B (UV6) 285-350nm UV-B (UV21) 285-350nm Waldmann UV 1000 KL, UV1000KL, UV 1000KL, UV1000 KL Monotherapy: 26 F85/100W-PUVA 26 F85/100W-UV6 26 F85/100W-UV01(TL01) Combined therapy: 20 F85/100W-PUVA + 6 F85/100W 16 F85/100W-PUVA + 10 F85/100W- UV01</p>	<p>UV 5040 KL Combined therapy: 26 F85/100W-PUVA + 16 F85/100w-UV21 or 24F85/100-PUVA + 16 F75/100W-TL01 Monotherapy: 40 F85/100W-PUVA or 40 F85/100W-UV6 or 40 F85/100W-TL01 Waldmann Starter (for 100W tubes)</p>	



<p>ML1094 Electronic ballast for the UVA 7001K Waldmann 550082070-00063079 For Lamps Waldmann UV 100W. 230x40x30mm</p>		
<p>F20T12/BB 31745 Luz Azul TL20w/52. TL20w/03T (ref.9511780661). Fototerapia Blue Vickers ref. 80030 BL4. AMEDA ZUG mod. FTL. (TL20W/03T) TL 20W/52 Philips LÂMPADA TL 20W/52 RS G13 Phototherapy of Hyperbilirubine</p>	<p>Philips F20T12/BB 31745 Luz Azul TL20w/52. TL20w/03T (ref.9511780661). Fototerapia Blue Vickers ref. 80030 BL4. AMEDA ZUG mod. FTL. TL20W/01RS. UV-B Medical. (TL20W/03T) Phototherapy of Psoriasis. TL 20W/52 SLV. F20T12/BB 20 Watt T12 Fluorescent Bulb Bilirubin Blue. Philips TL20/52 AIR-Shilds mod. PT.533-1E Phototherapy of Bilirubin</p>	 
<p>TL20W/01 RS. UV-B Medical. This is the most efficacious waveband for the treatment of psoriasis, and is therefore much more focused. What's more, because the overall dosage of this narrowband radiation can be closely controlled, these lamps are suitable for home therapy.</p>	<p>Lamps have a very narrow waveband emission: between 305 to 315 nm with a peak of 311 nm. The UVB Narrowband (01) low pressure fluorescent Phototherapy lamps are proven to be the most effective lamps for treatment of skin diseases. Can be used for treatment of: Psoriasis and Parapsoriasis Vitiligo Atopic Dermatititis Mycosis fungoides Other Skin diseases Prescription Required!!!</p>	
<p>PHILIPS TLK 40W/03 40W T12 24 LONG 47V VHO SUPER ACTINIC 308072 LIGHT EMISSIONS TECHNOLOGY TLK 40W/03 Philips BL TL-K 40W/10-R With Reflector</p>	<p>Super Actinic PHILIPS TLK 40W/03 For treatment of Hyperbilirubine. Philips Actinic BL TL-K 40W/10-R UV-A with reflector. Philips TLK 40W/10R UV-A</p>	 
<p>ML7719, F20T12 C-50. F20T12-C50, TL20W/33RS (F20T12/D). AMEDA ZUG mod. FTL. F20T12/D / ECO "Daylight". Average Life Hours 9000, , Length 24 In, Pre Heat, 1 1/2 In (38mm) Diameter</p>	<p>For treatment of hyperbilirubinie F20T12/C50; F20T12/C50/ECO, Fluorescent Lamp, T12 Lamp Shape, 20 Watts, F20T12/C50/ECO, Base Medium BiPin (G13), Color Temperature 5000, Lamp Tone Enhanced Daylight, CRI 90, Rated. Philips TL-M RS 20W 33-640 - 59cm</p>	
<p>ML9229 UV Lamp UV083245 UV LAMP K-14 5000 HOURS Brand: Olympus Manufacturer code: KAT14L(UV083245) LPC0005,</p>	<p>KAT14L(UV083245) UV LAMP K-14 5000 HOURS Original Lamp Lâmpada UV 083245 Lamp UV Lamp UV 083245</p>	
<p>ML1490 Sylvania F14w/133 G13 Coll White 4300K. 5410288000398. Dimensões Prod. Comprimento (mm) 375.4 Diâmetro 26</p>	<p>Used for many applications such as frige lighting and under cabinet lighting. As used in various juke boxes. 14w T8 Fluorescent tube. Manufacturer EAN: 5410288000398 F14W/133-T8, Lamp 14W, 220V, G13, 0,395A. 361mm</p>	
<p>GE F32T8/GO/ECOCVG Fluorescent Shatter- Resistant Gold UV Blocking Light Bulb Case of 36 length 48.000 in (1219.2 mm) Diameter (DIA): 1.000 in (25.4 mm)</p>	<p>GE F40/GO/CVG CovRguard Fluorescent T12 Shatter- Resistant Gold UV Blocking Light Bulb. The exclusive gold sleeve blocks all wave lengths below 520 nm. Case of 30</p>	
<p>ML4744 Philips brand T5 germicidal fluorescent lamp. 40 Watts and 34" with single pin base for use in germicidal sterilization Compatible Ballast - ICN2P60SC351. (Please always confirm compatibility with original equipment manufacturer</p>	<p>Philips Lighting 292672 TUV36T5/SP (G36T5/SP) Philips 40W 34in T5 TUV. Philips TUV36T5/SP 40W. Bulb has applicaiton in both professional water and air disinfection units.</p>	











<p>ML3124 MEDESTIME-Bilicrystal System Duo 2 PLS11w/25000K/2G7 DZS 11/75W DZC 11/75W, BC101/8 Only Medical Lamps</p>	<p>MEDESTIME LAMPS MEDESTIME PLS 11w/25000K/2G7 MEDESTIME- Bilicrystal G11W-4P/BLUE DZS 11/75W DZC 11/75W, BC101/8 For Bilicrystal phototherapy devices. BiliCrystal® lights tubes for intensive treatment of hyperbilirubinemia. Lâmpadas Medestime Lâmpadas Medestime Lamp Medestime</p>	
<p>ML6102 Dulux L Blue 18 W/71 2G11 Osram Dulux L Blue 4Pin. Dulux L 18w/71 Dulux L18w/71. Phototherapy Ardo (Ameda) Phototherapy Bilibed</p>	<p>Danger 2M21010. RADIUM Ralutec Long 18w/71 2G11. Fototerapia Luz Azul Only Medical Lamps</p>	
<p>ML7731 Electronic Ballast PC 1x18-24 W Basic Lamp TC-L, TC-F L145mmx H 21mm x W 22mm Dulux L Blue 18 W/71</p>		
<p>ML4494 WEDECO KATADYN 8013927 Lamp U.V. Wedeco Katdyn K-36L. WEDECO KATADYN UV Lamp, 36W. 8013927. WEDECO KATADYN UV Lamp K-36L Part No.8013927 germicidal lamp</p>	<p>WEDECO KATADYN UV Lamp, 36W.8013927 (Lâmpada para equipamento de esterilização) Lamp for sterilization equipment Lamp UV. Wedeco Katdyn K36L 295x25mm</p>	
<p>ML1877 Sylvania 21625 - F15T8/BLB Fluorescent Tube Black Light</p>	<p>Product Features Color/Finish: Black Light Blue Average Lifetime: 7500 hours F15T8/BLB 15 W MINIATURE FLUORESCENT UV LAMP 18 INCH G13 BASE BLACK LIGHT</p>	
<p>ML4415 Sankyo Denki G8T5E / 8W Lamp UV-B G8T5E / 8W - Sankyo Denki 8W T5 11.3IN GERMICIDAL SANKYO UVB UV-B 306 NM, Length: 290mm, Diameter: 15.5mm, Life: 3000 Hours</p>	<p>ML4415 USHIO G8T5E / 8W Lamp UV-B G8T5E / 8W - USHIO 8W T5 11.3IN GERMICIDAL USHIO UVB UV-B 306 NM, USHIO 3000318 Length: 290mm, Diameter: 15.5mm, Life: 3000 Hours</p>	
<p>ML9006 Sankyo Denki G6T5 6W Germicidal UV-C Tube G.5 XG1068/42V/6W/G5, TUV-6W-T5- G5-D16-X-226/3. G6T5/OF RG3 6W Puritec HNS G5 6W. G6T5/OF RG3 TUV-6W-T5-G5-D16-X-226/3 PHILIPS, XG-1068 O RBITEC XG1068/6W , 34-0013-01-UVP, XG1068/42V/6W/G5 ORBITEC XG1068/42V/6W ORBITEC, TUV- 6T5-6W-250-PK PHILIPS, LT6W-UV- FILTERED-C NARVA, 4DA16, LT6W/UV-C , T5/GERMICID/6W/G5, 020010, 23375, 7856, 8596.76,</p>	<p>ML2079 SANKYO DENKI (JAPAN) 20W G13 588,5mm Germicidal Lamp G20T10, Germicidal lamps containing ultraviolet (UV) light, can be used to kill germs or sterilize air, sterilize water or liquids and also sterilize surfaces for various materials. Sankyo Denki G20T10 20W. Sankyo-Denki G20T10 UVC 20w G13, length 588,5mm, diameter -32,5mm, UV 7,5w</p>	 <p>GERMICIDAL LAMP G20T10(GL20) 10 pcs CAUTION RADIATION ULTRAVIOLET Never look into the burning lamp. Protect your skin and eyes against the ultraviolet rays. SANKYO DENKI CO., LTD. JAPAN</p>
<p>ML9030 Sankyo Denki 30W G30T8 Ultraviolet 99V UVC T8. 30 watt 99 volt T8 Purification Germicidal UVC Fluorescent Lamp</p>	<p>ML9016 Sankyo Denki G16T5 16W Germicidal UV-C Tube G.5 XG1068/6W, 34-0013-01-UVP, XG1068/42V/6W/G5, TUV-6W- T5-G5-D16-X-226/3</p>	 <p>ML9030</p> <p>ML9016</p>




<p>Keeler 1011-P-7034, 1011-P-5055</p> 	<p>Keeler 1011-P-7042, 1011-P-5223</p> 	<p>Keeler 1011-P-7050, 1011-P-5047</p> 
<p>Keeler 1012-P-7000, 1012-P-5006</p> 	<p>Keeler 1012-P-7001, 1012-P-5014</p> 	<p>Keeler 1012-P-7002, 1012-P-5102 1030-P-7195</p> 
<p>Keeler 1012-P-7003, 1012-P-5110</p> 	<p>Keeler 1012-P-7004</p> 	<p>Keeler 1012-P-7006 Keeler 1012-P-5241</p> 
<p>Keeler 1013-P-7001, 1013-P-5091</p> 	<p>Keeler 1013-P-7009</p> 	<p>Keeler 1024-P-7156, KEELER-1024P7156</p> 
<p>Keeler 1299-P-7075 (6V 10W)</p> 	<p>Look at the bottom of your bulb to identify if your instrument is 2.8V (Dry cell battery) or 3.6V (Rechargeable battery).</p>  <p>Keeler</p>	<p>Riester 10413 Riester 10421</p>  <p>Riester 10472 Riester 10473</p> 
<p>Riester 10488 Vacuum Lamps 2.7 V</p>  <p>Riester 10489 Vacuum Lamps 2.5 V</p> 	<p>Riester 10600 Riester 10601</p> 	<p>Riester 10602 Riester 10605 Riester 10606</p> 
<p>Riester 10607</p> 	<p>Riester 10608</p> 	<p>Riester 10609 Riester 12321, LED 3,5 V</p> 
<p>Riester 10610</p> 	<p>Riester 10611</p> 	<p>Riester 10615. Riester 10627 (LED 3,5V)</p> 
<p>Riester 10620</p> 	<p>Riester 10625 Riester 11429</p> 	<p>Riester 11320 Riester 11429</p> 










<p>ML4488 6v 30w (6V-5A) P27D S.4179</p>	<p>I/6V/30W/P27D/T5.5/FC8/100H I/6V/30W/P27D/T5.5/FC8/CL/100H S-0.4179, S-4179</p>	
<p>ML4489 6v 33w P28S 6V 5,5A ML4489 6V/5.5A</p>	<p>I/6V/5.5A/P28S1/S8/FC8/80H I/6V/33W/P28S1/S8/FC8/80H/2756H I/6V/33W/P28S1/S8/FC8/CL/80H/ I/6V/33W/P28S1/S8/FC8/CL/80H</p>	
<p>ML4490 6v 15w P25D 6V 2.5A P25D Lamp OP2130 6v 15w P25D</p>	<p>77903-6V-15W-P25D NIKON 06-015-12N CANIMPEX 06-015-12N MICROLITES 6VOLT-15WATT-P25D-T8 NIKON 41100 EIKO, 41100 UPC Eiko 77903 Eiko 6v 15w P25D Lamp OP2130 6v 15w P25D</p>	
<p>ML4452 70464 RODENSTOCK 70464 12V 25W P25S RODENSTOCK 70464 70464 12V-25W P25S</p>	<p>70464 RODENSTOCK, 12-025-80 CANIMPEX, 12-025-80 MICROLITES, 70464/12V/25W/P25S RODENSTOCK, 125547, 125559/12V/25W/P25S 685050 IMPERIAL-OPTICAL, 4300410-10-CE RODENSTOCK, 6300300040 RODENSTOCK LT70464/12V/25W/P25S,</p>	
<p>ML4456 2756/506 6v 15w 70506 6v-15w. RODENSTOCK 70506, 5190300160. Osram 70506. 2756/506. O-70506 6V 15W SCBAY W/TRIANGL FLANGE. 125545, 125558, LT 70506,</p>	<p>210164 RODENSTOCK 467540Z MICROLAMP, R-02000 RODENSTOCK R02000 RODENSTOCK RD-506-82 RD70506/6V/15W/P25S RD70506/6V/15W</p>	
<p>ML7209 3893NK 2129 6v 15w P25D 99-72910, 04877334010 UPC 031293410948 EIKO 4030524400 TOPCON</p>	<p>3893/NK 2129 6v-15w P25D 77910 NIKON, 77910 UPC 77910 USHIO, 78097 NIKON 125530-LT77910 125530/6V/15W/SPECIAL 125530/LT77910, 125530, 212915/6V/15W BAILEY 8000242 USHIO. 2129 Nikon</p>	
<p>ML7291 Lamp 3893-1 6v 15W BA15D 125591</p>		
<p>ML7292 Lamp 3893-2 6v 18W BA15D</p>		
<p>ML7293 Lamp 3893-3 9v 18W BA15D</p>		
<p>ML7294 Lamp 3893-4 9v 20W BA15D</p>		
<p>ML7715 6V 15W BA15S Hosobuchi OP2118X, Lamp 3894/2</p>		
<p>ML7061 6V 15W BA15S Frosted Lamp 3894/0, Frosted</p>		
<p>ML7717 9V 18W BA15S, Lamp 3894/1, OLYMPUS, 8-M103 8000306, 031293410986. OL8M-103, OP2305.</p>		
<p>ML7075 12V 10W BA15S , 5008, Lamp 3894/3</p>		



<p>ML4009 NARVA 130676/230V/75W 230v 75w BA21D3 H38984 230V-75W BA21d3 60x26,2mm Silver-side</p>	<p>NARVA 130676/230V/75W BA21D3 H38984 230V 75W BA21d3 60x26,2mm Silver-side</p>	
<p>ML4181 PHILIPS 13244 24V 50W BA21-3S</p>	<p>PHILIPS 13244 24V 50W BA21 S. BA21S. BA21-3S I/24V/45W/BA21S3/G12.5/CL/13244 PHILIPS 60x40mm</p>	
<p>ML7017 24v 45w, 24v-45w, 24V/45W BA21D3d, Medidas 40x59mm I/24V/45W/BA21D3/G12.5/300H/H27 3411/24V/45W/BA21D3/G40MM/CC 8/26.5MMLCL/62MMMOL/600L/ O3411/24V/45W/BA21D3/G40MM/C C8/26.5MMLCL/62MMMOL/600</p>	<p>24v/45w, 24v 45w BA21D3d 24V/45W/BA21D3/130030, 0574/VC 130030-H273820-BA21D-40X-60-24V-45W 130030/24V/45W/BA21D3 NARVA 130030/H273820/BA21D-40X-60-24V-45W BA21D-40X60MM-24V-45W/H273820/130030 BBT-OPTOTAL/24V/45W/BA21D/G12.5 H-273820/BA21D-40X-60-24V-45W ORBI</p>	
<p>ML4541 Philips 6144-F 16v 8A BA21D3 Philips 6144-F 16v-8A BA21D3 Philips 6144-F 16v 128w BA21D3 Philips 6144-F 16v-128w BA21D3 6144F-16V-128W-BA21D3 PHILIPS. 21824 Dim: 63MMx40MM</p>	<p>16-128-06D MICROLITES 16V/8A/BA21D,6144-F DR-FISCHER 6144F/16V/8A/BA21D/G13/26.2MMLC L/65MMMOL/25H PHILIPS 112900-P6144F-BA21D-41X65-16V-8A 112900/P6144F/BA21D-41X65-16V-8A 112900, 370110U, 307843202</p>	
<p>ML7433 Trumpf Helion Lamp 12v 50w, 12v 120w, 12v 150w Lamp 12V 50W code 1388308 With DBGM 913 Sockets Trumph Helion</p>	<p>TRUMPF HELION</p> 	
<p>ML5614 Martin 89-907-02-04 Martin lamp set refª: 89-907-02-04 for Martin X8,</p>	<p>Martin lamp set 89-907-02-04 Martin lamp 89-907-02-04 Martin 89-907-02-04 Martin 899070204</p>	
<p>ML5052 Martin Berchtold CZ 908-22 (Martin OT-Light 86-907-01) (86-907-01-04). 24V 150W ML700, /ML701, C-570 Suportes DBGM 913 Incl. ML4409 Martin Berchtold CZ 908-24 22,8V 120W G6,35 ML7456 Sockets Martin ML500: 08-063-00-20</p>	<p>ML7000 Martin 86-905-01 : 24v 120w , 24v-120w Martin ML500, /ML501, C-450 Martin 86-907-01 : 24v 150w , 24v-150w Martin ML700, /ML701, C-570 Martin 86-910-02 : 24v 250w Martin ML1000, /ML1001, C-950 G, 24v- 250w With DBGM 913 Lamp holders ML4409 Berchtold Chromophare CZ 908-24 22,8V 120W G6,35</p>	
<p>ML5050 Lamp base socket holder Original DBGM 830. for Martin ML-701 to connect DBGM 913</p>  <p>ML5050: Only DBGM 830</p>	<p>ML5002 Holder Lamp, Socket Original DBGM 913 (DBGM913) holderlamp, lamp holder Lampholders, Sockets Socket for halogen lamp base G6,35 with screw contacts M3 unlooseable</p>	 <p>ML5002: Only DBGM 913</p>
<p>ML7105 Martin 86-905-01, 22,8V 120W 22,8V-120W Martin ML500, ML501, C-450 22,8V 120W G6,35 MT 905 (MT905),</p>	<p>ML7000 24V 120W G6,35 Trumph Helion L 24V-120W G6,35 Martin ML500, ML501, C-450 MT 905, MT-905, MT905</p>	 












<p>145RH Truphatek 512501 Truphatek Green Spec-2 CL954 Gold. 512717 Xenon Lamp</p>	<p>512501/6, 2,5V 0,7A, M3.5x0.35 Laryngoscope Replacement Bulb Truphatek 512501FO Green SpeX. The same lamp in LED (145RHL)</p>	
<p>ML7701 Truphatek 513501 For Rechargeable Handles</p>	<p>Xenon Lamp 3.5V 0.7A, 4.8 OD M3.5x0.35, Truphatek 513501 Laryngoscope Replacement Bulb</p>	
<p>196RH Truphatek Green Spec TRUPHATEK Green Spec</p>	<p>Truphatek Green Spec 2,5V 0,7A Laryngoscope Replacement BULB</p>	
<p>ML7809 Lâmpada para Laringoscópio FAZZINI Laryngoscope Laryngoscope Lamp FAZZINI</p>	<p>ML7136 Penlon 55094 Bulb 2.5V / 700mA halogen lamp. 3v 0,4A</p>	
<p>ML4091 Laryngoscope Martin 19-353-00. 1935300 MARTIN, MILLER, FOREGGER ML4091A (With Spring)</p>	<p>ML7199 Laryngoscope Martin 19-395-00-04. 193950004 MARTIN, MILLER, FOREGGER ML7199A (With Spring)</p>	
<p>ML5429 Replacement lamp Type 507. Nystagmus Spectacle with fixed lenses by Prof. Frenzel 501</p>	<p>The Nystagmus Spectacle Type 501 can be operated with a 4 Volt; 0.23 Ampere; 1 Watt replacement lamp. (Thread diameter: 3.5 x 0.35mm).</p>	
<p>NEITZ-L29 Ophthalmoscopes bulb NEITZ-L29 4V-2.5V (0.63A). Laryngoscope Replacement Bulbs</p>	<p>ML2849 Bailey Rectoscope 3.5v Morgan A24 2850387690, MFA99325, 8714681114346. Bailey MFA99325 - Rectoscope 3.8V Morgan A24</p>	
<p>ML8811 Lamps from Fiberscope Pentax REF^a BS-H2 Pentax lamp BS-H2 is used for the mini light source BS-LHS2</p>	<p>BS-H2. Pentax Lamp for BS-LH2 Light Source Lamps from Fiberscope Pentax REF^a BS-H2. HALOGEN LAMP. Diagnostic Lamp 2.5V 1.8W available only in a set of 5 lamps</p>	
<p>ML4912 Lamp. 2,7v p/Tonometer Perkins Tip. Haag-Streit UK/CC 1146005</p>	<p>PERKINS 1146005 PERKINS 5801006 PERKINS HAND HELD APPLANATION TONOMETR PERKINS-5801006 PERKINS MK1,MK2 PERKINS-5801006 Laryngoscope Replacement Bulbs</p>	
<p>ML5486 Laryngoscope Medicon lamp 45.01.15, Medicon 45.48.66</p>	<p>Lâmpada para Laringoscópio Medicon 45.01.15. Medicon 45.48.66 Laryngoscope Replacement Bulbs</p>	
<p>ML1610 Diagnostic Lamp 2.5V - Medicon 45.48.70 Crypton lamp</p>	<p>Medicon lamps original brand 45.48.70 LARYNGOSCOPE 2.5V 45.48.70 MEDICON Diagnostic Lamp 2.5V - Medicon 45.48.70</p>	
<p>ML1612 Diagnostic Lamp 2.5V - Medicon 45.48.72 Xenon lamp</p>	<p>Medicon lamps original brand 45.48.72 LARYNGOSCOPE 2.5V 45.48.72 MEDICON. Diagnostic Lamp 2.5V - Medicon 45.48.72</p>	
<p>ML7708 Gowllands nº 12 Gowllands 310 ML7709 GOWLLANDS-13 GOWLLANDS-13</p>	<p>ML7708 Gowllands nº 12 For otoscope Gowllands 310 ML7709 GOWLLANDS 13 Ampoules de rechange nº 13</p>	
<p>ML7177 W26 (ACMI) medidas 4x15mm, American Cystoscope (Wappler)</p>	<p>ML7177 Wappler, American Cystoscope W26 ACMI 4x15mm, ACMI W26 4x15mm</p>	













<p>ML5464 Lux lamp, KAWE 2.5V/H, E-28933, xenon bulb 2.5V, 12.75144.013, XL-otoscope, F.O. Laryngoscope, fiber optic diagnostic 28 933 light. Kawe Model Piccolight F.O.</p>	<p>ML3198 Proctoscope EFFNovac 426 lamp 3,5V. EFFNER -Proctoscope Protoscope 426 lamp 3,5V 55x 313mm</p>	 <p style="text-align: right;">ML3198</p>
<p>ML5498 Lamp 3.5V to a light source Karl Storz D-7200 3,8V lamp for Karl Storz D-7200 Model 102200B</p>	<p>ML5428 Riester Laryngoscope Lamp (used for blades sizes 2-4) GIMA KAWE TIMESCO AESCULAP</p>	 <p style="text-align: center;">ML5498 ML5428</p>
<p>ML4533 Hosobuchi M-1021 (M4.5 thread). 6V 0.2A. Keeler 1030-L-2085. M1021 HOSOBUCHI</p>	<p>Lâmpada Laringoscópio 6V (25x5,5mm) Bulbs for Laryngoscopes 6V (25x5,5mm) Laryngoscope Replacement Bulbs It's a fixation lamp for Topcon SL-3D</p>	
<p>ML8228 Lamp. Haag Streit para Laringoscópio Bulbs for Laryngoscopes</p>	<p>Lâmpada laringoscópio 8V (25x5,5mm) Bulbs for Laryngoscopes 8V (25x5,5mm)</p>	
<p>ML5905 12V FROSTED. 12-0.01-E3. Zeiss 38-61-07 Lamp.Fosca refª. 386107 HITEC, 12v</p>	<p>KRIPTON (4x21x 5.5mm) P/PONTA DE LUZ Z386107 12V T5.5MM METAL SLEEVE, Zeiss Slit Lamp Fixation Device Bulb</p>	
<p>ML9911 Lamps for glasses Type 507 and Type 502 DEHAG Brand, 3.5 x 0.35mm (Frenzel ref. 37077)</p>	<p>Lâmpada para Oculos Tipo 507 e Tipo 502, Marca DEHAG (Frenzel 37077) Nystagmus spectacles Dehag 507</p>	
<p>ML7699 Neitz L 29 4V-2.5W (0.63A) Halogen Ophthalmoscope Lamp. Neitz L-29, Neitz L29</p>	<p>ML7630 Neitz L 30 3V 2W T1.5, Neitz L30 Ophthalmoscope Lamp. Neitz L-30</p> <p>ML7631 Neitz L 51 6V 10W T2.25, Neitz L51 Ophthalmoscope Lamp. Neitz L-51</p>	
<p>ML1901 Wotan 391808 BL 7248</p>	<p>ML2511 Dimed DMV520 Lamp bulb XHL 2,5V from Otoscope Dimed DMV520</p>	 <p style="text-align: center;">ML1901</p>
<p>ML2145 Lamp 57310 6v 15w Sinoptoro F57000/58000 Synoptometer by Kueppers F57000/58000 6v/15w Falange base. Kueppers 58000. Kueppers 57000 Spare Bulb, centred 6v / 15w f. 57000 / 58000</p>	<p>6v 15w BA15S, 6v-15w BA15S 50x25mm, Falange Base 37x32mm Spare Bulb Precentred</p> 	
<p>ML4209 Bulb 12V 30W for perimeter HFA-II and HFA-730, [8059023] Halogen spare bulb for Zeiss/Humphrey perimeter HFA-II and HFA-730, 12V 30W, according. Zeiss 2660021106082 12V/30W</p>	<p>ML4209 Zeiss 266002-1148-110 Zeiss/Humphrey perimeter HFA-II and HFA-730, 12V/30W, Humphrey HFA-II Zeiss 266010-0030-323. 30323 Zeiss, Zeiss 266002-1106-082 12v 30W Perimeter Humphrey HFA-730 Zeiss 266002-1106-082/14078</p>	
<p>ML5055 Lamp for Colposcope 6v 30W BA20D (6v 5A) 60x35mm Lamp for Coloscópio Leisegang Black Side</p>	<p>ML5056 Lamp for Colposcope 6v 10W BA15D 46x25mm Lamp for Colposcope LEISEGANG Black Side</p>	 <p style="text-align: center;">ML5055 ML5056</p>



<p>ML8088 Lamp for Colposcope Lamp Colposcópico 12v 20w BA15D Silver Side Is the same like ML8089</p>	<p>ML8089 Lamp. for Colposcope 12v 20w BA15D 12v-20w BA15D Without Silver Side Is the same like ML8088</p>	
<p>ML4508 Philips 220-230v 115w E27 Special lamp E27 GLS Clear 230V-240V 115w Carbon filament</p>	<p>E27 GLS Clear 230V-240V 115w Carbon filament Antique Lamp. Vintage Lamp E27 GLS Clear 220V-230V 115w Carbon filament, Coal Filament Filamento de carvão</p>	
<p>ML3009 M01129 HIKARI 6v 10w SPECIAL-BASE. HIKARI M01129 6v 10w HIKARI M01129 6v-10w</p>	<p>5629-5/6V/10W 5629-5, M-01129 HIKARI M01129/6V/10W HIKARI M01129 HIKARI, Lamp 5629/5</p>	
<p>ML5501 1168 GHAMBS 6v 15W P45D 1168 GHAMBS 6v-15W P45D LRLS C 1000, LRLS C1000 LRLSC 1000, LRLS C1000</p>	<p>I/6V/15W/P45D/T6/FC6H/100H. 1168/6V/15W, 1168, I/6V/15W/P45D/T6/FC6H/100H/ I/6V/15W/P45D/T6/FC6H/100H</p>	
<p>ML4309 Roche 12v 100w Halogen Lamp Assembly 05070856001 Roche, Lamp Halogen 12V/100W Assy., 28051763001 ROCHE 2076896001</p>	<p>ML4309 Roche 12v-100w Halogen Lamp Assembly 05070856001 Roche 12V 100W, 12V-100W Roche 12v 100w Halogen Lamp (Man.ID 0507083001) Cobas Integra 400 Cobas Integra 700</p>	
<p>ML6606 6V 20W P21D 6V-20W P21D 177160-W, 177160W, 3251 177160 LEITZ/LEICA/WILD 10177160 LEITZ/LEICA</p>	<p>177160-W, 177160W, 3251 177160 LEITZ/LEICA/WILD 10177160 LEITZ/LEICA 302843384/6V/20W BT-70249-B, BT70249B LT70249, LT 70249, 125505 0070249 WILD/LEITZ 70249W/6V/20W, 70249W 125505/6V/20W/P21D</p>	
<p>ML4403 Roche 10641600001 12v 20W P15 (ASSY 705-0304) 12v-20w P15 (705-0840) The lamp used for Roche Modular P800 Hitachi 94K2136 Hitachi 94K 2136 for Roche Modular P800</p>	<p>11A/12V/20W/P15 HITACHI 12-020-20/M, 12V/20W/P15 705-0840/12V/20W HITACHI 705-0840 HIKARI, 5629/2/12V/20W 131923-OR1390-P15-28X38-12V-2- FILS-20W, 131923/12V/20W/P15-S 131923/12V/20W/P15 131923/12V/20W, 12v 20w, 12v-20w Hitachi/Roche 705-0840 for 700&900 series Chemistry Analyzer</p>	
<p>ML6004 OLYMPUS Halogen Lamp ASSY Model PG55134 12V 20W for AU400 AU640 Bio auto analyzer. Photometer lamp PG55134 12V-20W</p>	<p>OLYMPUS Halogen Lamp ASSY Model PG55134 12V 20W for AU400 AU640 Bio auto analyzer ANALIZADOR OLYMPUS 400 wholesale pg55134 12v 20w lamp bulb</p>	
<p>ML5599 Hitachi K9P1641 12v 50w for C501 device 12v-50w for C501 device</p>	<p>Hitachi K9P1641 12v 50w, 12v-50w</p>	













<p>ML4058 Lamp 6703/AX2 12V 20W G4 M-01131 HIKARI 12V-20W G4</p>	<p>Lamp 6703-AX2 12V 20W G4 M-01131 HIKARI 12V 20W G4 Base Alumínio: 19x10mm</p>	
<p>ML4054 Lamp 6703/AX, 1913 24V 50W G6,35</p>	<p>Lamp 6703/AX19, 1913 12V 20W G6,35</p>	
<p>ML4059 Lamp 6703/1, 6703-1 24V 50W G6,35 Base Alumínio: 19 x13mm</p>		
<p>ML4051 Lamp 6703/AX1, 197, M-01131 HIKARI, M01131 HIKARI, 12V 20W G4</p>		
<p>ML5891 H3 12V 55W PK22S Multi Task Halogen Dental Lamp</p>	<p>01007-BLISTER-PACK EIKO, 01007-BP EIKO 01007-BP USA, 01007-BULK EIKO, 01007- BULK USA, 01007-D EIKO, Guerra 6705/2 01007 EIKO, 01007 UPC. 3M 618 Microfilm QUANTIMATIC Microfilm. Ampco: CONTENDER Dental Lamp</p>	
<p>ML4099 48695 NARVA 24 65w PK22S Lamp 6705/3, Lamp 6705-3 24v 65W PK22S</p>	<p>NARVA 48695, Lamp 6705/3, 6705-3 122758-H192465-PK22S-12X42-24V-65W-H3, 122758/24V/65W/PK22S 122758/24V/65W, Guerra 6705/3 122758/H192465/PK22S-12X42-24V-65W-H3, 122758, 844409004562, 24/65w PK22S, 48695 NARVA 48695</p>	
<p>ML7694 24V 60W PK22S Black Top HBL-564, HY2460 Dental lamp 24V 60W PK22s 24V 60W PK22S Flat Contact Lamp suitable for Belmont Clesta 501, 502 & 505 dental light. H3 24v 60w Pk22s Black Top Dental Bulb - Belmont. HY 2460 H 3 F BLACK TOP</p>		
<p>ML8877 Kowa K9L9 FU3A lamp Kowa reference K9L9FU3(A), 300W Xenon flash bulb for fundus cameras Kowa FX- 50R, FX-500C, FX-500R, VX- 10, RC-XV, RC-XV2, RC- XV3,P PRO-1, PRO-2, PRO-3, Non-Mydriatic, VX-1, VX10</p>	<p>Flash Tube for Kowa fundus cameras P/Camera Fundus Kowa Pro III. Lamp FV3A for Kowa Mod. VX.10I. FU3 (A) for Kowa VX-10 Kowa K9L9- FU3A lamp. Lamp Angiografo KOWA RC-XV3. 99- K052. Kowa-Nr. K9L9FU3 (A). (GSI-128 (UCC/EAN-128) K 9 L 9 FU3. Lâmpada para Angiografo Kowa RC-XV3, Kowa VX-10 Angiografo Kowa VX.10I, Angiografo Kowa VX101, L9FU3A, K9L9FU3A Kowa Pro III</p>	
<p>ML9115 Lamp PF712E 230V 15W E27 PHILIPS. PF712E 230V-15W E27 Lâmpada Câmara Escura ML9116 Dr. Fischer 77841422 Photolamp 230v 15w E27</p>	<p>ML9115 Lamp. PF712E, Raio X 230V 15W E27 PHILIPS I/230V/15W/E27/A19/RED- SAFELIGHT/1000H/PF712E ML9116 Dr. Fischer 77841422 Photolamp 230v 15w E27 Dr. Fischer 230v-15w E27</p>	
<p>ML1737 Sockets, holderlamp (lamp holder) p/ Microscópio GZ6,35 ref. 865R (Refª 126RW4) Lamp Holder</p>	<p>BW 865R/986 (865R/3459)</p>	



<p>ML8721 Lamp Microscópio Zeiss 6V 10W Halogen bulb 386108-0000-000 5120202, Wotan 64223, 386108. 624, 624 OPTICAL, 700-065/6V/10W, 700-065 907-090 LEITZ/ZEISS, 1000-531 USHIO, 1000-857 USHIO6519/1 6v-10w</p>		 <p>ML8721</p>  <p>ML6877</p>
<p>ML6877 Lamp Microscope Zeiss 6V 10W MID 3600 p/espelho Frontal BFW (from portable headlight). BFW (ref. MID3600) is Ansi ESA. 3801-68 ZEISS, 3861-08 ZEISS, 3875-7880 ZEISS, 3878-0030 ZEISS, 3878-6570 ZEISS.</p>		 
<p>ML6878 Lamp Microscope Zeiss 6V 10W MID 3600 p/espelho Frontal BFW (from portable headlight). BFW (ref. MID3600)</p>	<p>Lamp 64410S 6V 10W G4 ZEISS, 3861-08 ZEISS, 3875-7880 ZEISS, 3878-0030 ZEISS, 3878-6570 ZEISS.. Ansi ESA. 3801-68 Filament Axial</p>	
<p>ML8909 Lamp Dr. Fischer 24V 45W Ba15s, 7669 - 45 watt - 1.88 amp - 24 volt - Single Contact Bayonet Incandescent Miniature 34x57mm, 36x58mm</p>	<p>Dr. Fischer 24V 45W BA15s Lumen (lm) 800-1200 Life (h) 100, 00144500, 7669, 24v-45w -1.88A, 24V 1,88A Osram 7669 Osram 24v 45w BA15s Single Contact Bayonet (BA15s) Base</p>	  <p>ML8909</p>
<p>ML7027 Lamp Nikon 220-240v/20w BA15D 25x55mm 220-240V 20W BA15D</p>	<p>Nikon 230-240v/20w BA15D Nikon 230 20w BA15D Nikon 240v 20w BA15D Nikon 220v 20w BA15D 4844/2 22x52mm</p>	
<p>ML5285 NARVA 131338 5,2V 0,85A E10 Halogen, Osram 6406530. 5.2V 4.42W/EY10/68001 NARVA 68001 NARVA 082006. 5,2V 0,85A E10, 5,2V 850mA E10</p>	<p>ML5285B Halogen Lamp 5,2V 850mA P13,5s. Merck photometer. Fotômetro Merck para análise de águas, efluentes e processos industriais. 109749. NOVA 30, NOVA 60 Spectroquant <u>Only lamp without plastic support</u></p>	 <p>ML5285</p>  <p>ML5285B</p>
<p>ML4721 Philips 12241 12V 45W BA15S, (ref. 7569) (1195) 12v-45w, 49174 NARVA 78191 HELLA- MARINE,, 122454 , 122455 ORBITEC, A-12066 , A7569 ORBITEC,</p>	<p>12V/45W/BA15S/49174 NARVA 0490, 546 ALTAN, 2066 EUROPEAN, 2066 OSRAM, 2066 USA, 3113/12V/45W/RP35MM/BA15S NAM-YUNG, 3113 NAM-YUNG, 3210 SUNG-WHEE,</p>	
<p>ML7782 CHY 230V 50W BA15d CHY 230V 50W 131403 Dim: 79MMx26MM Black Top</p>	<p>American Optical Special 1217LF Lamp 131403 CHY 230V 50W Lamp CHY BA15D T8 230v 50w CHY 230V 50W BA15D T8 Topo Negro, CHY Black Top</p>	 
<p>ML7781 CHY 230V 50W BA15d American Optical Corp. 11082 (Item No: R3110) PH/50T8/114DC 230V Lamp CHY 50W</p>	<p>CHY 230V 50W BA15d T8 GENERAL ELECTRIC CHY 50W T8 114-230VAC MOVIE/SLIDE PROJECTOR LAMP CHY 230V- 50W, dim:79MMx26MM Lamp CHY 230V 50W Without Black Top</p>	   
<p>ML4708 FUJI BXR 10V 5A BA15s FUJI BXR 10V 50W 50W 10V INC10V-5.0A. length: 74,4mm, Diameter: 25,3mm, Life: 100Hr length: 79,4mm, Diameter: 25,3mm, BXR/10V/50W/BA15S ANSI 6057N/10V/5W/BA15S/T8/40.5MML CL/78MMML/100H PHILIPS 23213-2 WESTINGHOUSE 23270, 1000103 USHIO 048777106860 UPC</p>	<p>FUJI BXR 10V-5A BA15s FUJI BXR 10V-50W BA15s Ushio BXR 50W 10V INC10V-5.0A BA15S Single Contact Bayonet Base M1-1156 PHOTOPHONE 35mm Projector TO-61 PHOTOPHONE M-1-1157 PHOTOPHONE 35mm Projector M-1-1158 PHOTOPHONE 35mm Projector M-1-1159 PHOTOPHONE 35mm Projector</p>	 












<p>3900/CH 220v 15w E27, 220v-15w E27 Clear, LM-6, LM-25</p> <p>3900/SM 220v 15w, 220v-15w E27 Frosted, used for Topcon Lensmeter LM-6, LM-10 Hosobuchi O-3637, Ophthalmology lamps</p> <p>3901/CH 220v 25w, 220v-25w E27 Clear, Takagi LM-10DX (O-3564)</p> <p>3901/CH-20 220v 20w, 220v-20w E27 Clear, Takagi LM-10DX (O-3623)</p> <p>3901/SM 220v 25w, 220v-25w E27 frosted, Lâmpadas Oftalmologia</p>	
<p>3905/SM 220v 15w E27 Frosted Lamps 220v-15w E27 Frosted Lamps</p> <p>3906/SM 120v 15w E27 Frosted Lamps 120v-15w E27 Frosted Lamps</p>	
<p>ML5442 220V 20W E16, 220V 17W E16 Microscópio TAKAGI 375, 220V-20W E16 230V-20W E16, 220V 17W E16, I/220V/17W/E16/T7/300H Microscope Takagi Mod 375</p> <p>ML5452 220V 20W E16 Frosted, 230V-20W E16 Frosted 220V 17W E16 Frosted, 220V-20W E16 Frosted 230V 20W E16 Frosted, 220V-17W E16 Frosted Microscope Takagi Mod 375</p>	
<p>ML5054 120v 30w BA15D 48x30mm, 3897-2/120v/30w BA15D CH-120VOLT-30WATT</p> <p>ML2055 220v 20w BA15D 48x30mm, 3897-1/220V/20W BA15D CH-220VOLT-20WATT</p>	
<p>ML7321 220V 15W E12 Frosted, Lamp TOC 220V 15W E12 Frosted 30x48cm For Takagi</p>	
<p>ML5448 Takagi MT 375. Microscope Takagi MT-375 220v 20w E16, 230v-20w E16, 230v 20w E16, 220v-20w E16 Microscópio Takagi MT375 (MT-375). SÉRIE 50. 841055. 220V 20W E16 (230V 20W E16). 3899/CH Takagi Frontofocómetro MT-375 82715 TAKAGI MT 375</p>	
<p>ML5584 220V 10W E16 Frosted 220V-10W E16 Frosted, 3899/SM</p>	
<p>ML2155 Lamps MA886 - Bulb 220V 20W (DC Bayonet Half Silvered) For Swift M3200 Microscope. MA886 220V 20W 220V 20W BA15D Lamp. Bulb 220/240V 20W 0-3601 P220.0020.3601 BAILEY. Olympus Microscope MA886 Lamps Halogen Half Silvered 220v 20w BA15D</p> <p>Frosted side Silver side, Half Silvered</p>	
<p>ML5422 3896-SP 12v 40w BA15S (Hosobuchi O-3403 12v-40w BA15S) lamp 30403, 3896-SP 12v 40w BA15S (Hosobuchi O-3403 12v-40w BA15S) lamp 30403</p> <p>Silver-Side</p>	
<p>ML5423 3896-SP2 6v 3A BA15S, (6v 18w BA15S) 3896-SP2 6v 3A BA15S (6v 18w BA15S) 6v-3A, 6v-18W BA15S</p> <p>Silver-Side, Half Silvered</p>	



<p>ML6551 6701/1 SH37, 24V 35W E11 ML6552 6701/2 SH42, 24V 40W E11 E11 24v 40w 130065, H422440 ML6553 6701/3 SH52, 24V 50W E11 ML6554 6701/4 SH62, 24V 60W E10</p>	<p>equipment: Skylux SH-42, SH42 Medical SKYTRON Operation Theatre A-1-010-18, A-1-010-28 Medical A-101018 UNIVERSAL Microscope A1-010-18, A1-010-28 Medical A101018 UNIVERSAL Microscope SH-42 Medical, SH42 Medical M01088 DONSBULBS 6701-1, 6701-2, 6701-3, 6701-4 Skylux-Light</p>	
<p>ML7054 Lamp 2954, 1740 (13347Z) 6V 15W P25D/4 6V 15W P15d25-4a 6V-15W P15d25-4a 3800-18-1740 ZEISS, 3801-77 ZEISS, 13347Z ZEISS, 41091 EIKO, 41091 UPC, 125519 ORBITEC. Philips 13347Z 6V 15W P25D 3800-18-1740 MICROLAMP 52x18mm</p>	<p>RA14 ZEISS, RA16 ZEISS, RA18 ZEISS, RP-14 ZEISS, RP-16 ZEISS, RP-18 ZEISS, SM-18-1740 USHIO, SM-3800-18-1740 USHIO, SM-3800-18-1740 ZEISS, 99-71174 EFFNER, 99-71177 EFFNER, 380-018-1740, 1740, 2954, 3800-18-1740/6V/15W ZEISS, 3800-18-1740/6VOLT/15WATT ZEISS, 3800- 18-1740/15W/6V/BA15D MICROLAMP, 3800-18-1740/ZEISS INTERLIGHT, 3800-18- 1740 BULB-DIRECT</p>	
<p>ML4179 ZEISS Fundus 110V 50W, 110V-50W Camera Viewing Bulb Lamp Zeiss Fundus Camera Viewing Bulb Double Filament 110v 50w (390413 Zeiss) Duplo Filamento 110v/50w</p>	<p>ZEISS Fundus 110V 50W LAMP FOR ZEISS 390413-Z - 110V 50W LAMP FOR ZEISS SM-39-04-13 ZEISS, 39-04-13 ZEISS, 8000178 USHIO, H406 HYBEC, SM-39-04-13 USHIO Double Filament 110v 50w</p>	
<p>ML2977 Guerra 2855 1730 (Philips 13347C) 6v 15w P22D Osram 8018 6v-15w Microscope Lamp 6V 15W P122D</p>		
<p>ML5588 Microscope lamp 120v 20w, lamp. 120v 20w BA15D Silver side, Side-Mirror 56x35mm</p>		
<p>ML5555 Olympus CH2 microscope bulb. 120v 30w BA 15D. 2191740 (8-B-152 CH2 KHC CHB CHT KH-2 120V-30W Olympus Device) 317428. 2-S-650-023 OLYMPUS 2S650-023-MICROSCOPE NIKON, 9UD/120V/30W NIKON, Olympus CH-2 microscope 30-S11-DC, Silver side, Side-Mirror --120VOLT, 56x35mm, 120v 30w SB</p>		
<p>ML7485 Lamp Life 6,3v 5,5A P30D Lamp 6,3v 35w P30D for device Shimadzu vertical filament with two contact wires Measure: 64x25mm</p>	<p>Lamp for device Shimadzu Lamp 6,3v 5,5A P30D</p>	
<p>ML2056 GE-USA 115/120v 30w BLC 115-125V/30W, 120V 30W OL8B-151 Microscope 120v 30w OL8B-152 Microscope 8B152/120V/30W OLYMPUS</p>	<p>GE-USA 115/120v 30w BLC 115-125V/30W, 120V 30W 60,3x35mm 120v/30W Tungsten OL8B-151 Microscope 120v-30w OL8B-152 Microscope Olympus 8-B152 OLYMPUS, 8-C-152,</p>	
<p>ML7116 14623 Projection Lamp 17V 95W G6,35 Dental Lamp Philps 14623 95W G6.35 17V 1CT/10X10F</p>	<p>ML7117 14623P Projection Lamp 17V 95W GZ9,5 This ceramic based capsule is a 17V 95W bulb with a GZ9.5 base which is used in certain medical machine lamps. It has a colour rendering index of 99 and a colour temperature of 2900K with a 2000hr rated life</p>	



























<p>ML1228 Merivaara Merilux 22,8V/40W F7 Code: 485761 A 084 7805 22,8V 40W IRC Germany</p>	<p>Surgical Light Merilux Merivaara 22,8v 40w Special IRC Merivaara Merilux 485761. Merivaara Merilux X1, X3, X5 A 084 7805, A0847805, A 0847805</p>	
<p>ML1297 NEW ORIGINAL PROJECTOR LAMP BULB FOR OSRAM P-VIP 180/0.8 E20.8 180 0.8 E20.8</p>		
<p>ML6441 Perkin Elmer Y1940, Cermax VQ Y1940 300W Xenon Lamp Endoscopy Medical Light Source</p>	<p>Y1940, Perkin Elmer Cermax VQ Y1940 300W Xenon Lamp Endoscopy Medical Light Source</p>	
<p>ML4719 Lamp Ushio FC72T4 Single Base Stocker & Yale. Ushio US3000120 FC72T4 Single Base Stocker & Yale Lamps Bulbs.</p>	<p>Ushio US3000120 Ushio Fluorescent Light Bulb #3000120 (Replaces Stocker & Yale 973-510) 11 watt - 33 volt - T4 - Miniature Bi-Pin (G5) Base - 5,100K</p>	
<p>ML4818 DR MACH 22,8/24v-50w BA20D DR MACH 67100103 24V 50W Base Ba20d, used for Soloflex Black Top</p>	<p>Dr. Mach 22,8/24V 50W Ba20d Dr.Mach 67100103 avec Top Noir Forme d'ampoule: T11 Diamètre: 11mm Longueur hors tout: 56 maxmm Centre de lumière: 29+/-1mm Dr. Mach 67100103, Finish: Black Top</p>	
<p>ML7110 DR MACH 22,8/24v 50w G6,35, DR MACH 67100201 D-85560, (22.8V 50W G6.35) used for Soloflex, Uniflex, Trigenflex, Quintaflex, M3, Mach 130, 380, 400, 500, 700, ref. 1660/91; 67100201</p>	<p>Dr. Mach 380F, Dr. Mach 380DF, Dr. Mach 400F, Dr. Mach 400DF, Dr. Mach 500F, Dr. Mach 500DF, Dr. Mach 700 F, Dr. Mach 700 DF. 6601/1. 99-611160. 130748. 64650. MA-050, 1660, 1691, 17100201. Dr. Mach 22,8/24V 50W G6,35</p>	
<p>ML7445 Hanalux 2000 24V 50W (64445U)</p>	<p>Hanalux 2003, 24V 50W Osram 64445U 24V 50W GY6,35 Osram 64445 U 24V 50W GY6,35</p>	
<p>ML9110 Dr. Mach 22,8V 50W Dr. Mach 22,8V-50W G6,35 Axial (Topo Negro) Dr. Mach Black Top</p>	<p>ML1900 Dr. Mach 22,8V 40W Dr. Mach 22,8V-40W G6,35 Axial (Topo Negro) Dr. Mach Black Top</p>	
<p>ML4599 Dr. Mach 22,8V 80W 22,8v-80w XIR G.6,35 64668 XIR. Dr. Mach M2F Serie Dr. Mach 67100105</p>	<p>Hanauchrome 22,8v-80w XIR G.6,35 IRC 23v-80w, IRC 23v 80w Dr. Mach 6710 0105 22.8V/80W G6.35 XIR 22.8V 80W G6.35 XIR</p>	
<p>ML5593 12V 50W GY6,35 lamps for medical equipment</p>	<p>12V 50W GY6,35 Lamps for medical and dental equipment 39x11mm</p>	
<p>ML7799 Massive 22331/92/50 (223319250) 12v 50w GY6,35</p>	<p>MASSIVE 50W GY6.35 12V 22331/92/50 MASSIVE 12V 50W MASSIVE 12V-50W</p>	



<p>ML2084 Welch Allyn M21E001 Welch Allyn Solarc Miniature Arc Lamp M21E001 21W-24W</p>	<p>HID/M21E001/12V/21W WELCH/ALLYN M/21/E-001 M21E-001/21W-METAL-DISCHARGE LAMPTECH M21E-001 LAMPTECH M21E001-12V-21W WELCH/ALLYN M21E001/12V/21W M21E001 WELCH/ALLYN</p>	
<p>ML2089 Welch Allyn M50E021. 50W HID Metal Halide Lamp Assembly. 52V-50W. WELCH ALLYN WA- M50E021, Miniature Metal halide Arc lamp Welch Allyn M50-E021</p>	<p>W.A. M50E021 HID Metal Halide Lamp Assembly. WAM50E021 WA-M50E021 50W METAL HALIDE lamps 310 009 020 short arc lamp 50 W (Xenon lamp 50W M50-E021 HG) For Xion Medical mod. EndoStrob DX</p>	
<p>ML9802 Welch Allyn M50E021 XION. HID Metal Halide Lamp Assembly. 52V 50W. Metal halide lamp for Xion EndoSTROB 50W, 800h, 5.700K.</p>	<p>ML9802 Ushio M50E021 XION EndoStrob DX XION MEDICAL XION 310 009 020, XION 310009020. Xion fixed an additional metal part to the front</p>	
<p>ML2082 Welch Allyn 09500 (Welch Allyn MFI & VDX 24w Solarc Lamp 09500U) 9500 SOLARC- 60v 24W -MR-11-METAL. USHIO MFI24V Solarc CE344. 49500, 46100. Ushio AL-0950, AL-1824</p>		
<p>ML2082A Welch Allyn M24E00S-001 18-24W MR11 METAL HALIDE USHIO AL-1824 USA DA-L400 DLS-400 SOLARC LAMP .40 AMP 60V</p>		
<p>ML2088 Welch Allyn M21E00S-001 Solarc MFI Replacement lamp</p>		
<p>ML2098 Welch Allyn 09800 Welch Allyn 21W Solarc Bulb, Model 09800</p>	<p>WA09800 WELCH/ALLYN 21w Solarc Welch Allyn Video Colposcope 88002A, 88004A, 89002A, 89004A. Welch Allyn 09800-U 21W Welch Allyn</p>	
<p>ML2099 Welch Allyn 50/60W - AL-5060 APLICAÇÕES: Fibra óptica Fonte de Luz, Vídeo, Microscópio Fototerapia, Foco Cirúrgico</p>	<p>USHIO AL-5060 WA AL5060 50/60W WELCH ALLYN Ø 52mm x 58mm APPLICATIONS: Fiber optics Light Source, Video, Microscope Phototherapy, Surgical Focus and endoscopes</p>	
<p>Welch Allyn M10P004 WA M10P004 - PARABOLICA - WELCH ALLYN. MR11</p>	<p>Socket B10R001 - WELCH ALLYN Socket para lâmpada parabólica WA M10P004</p> 	
<p>ML2110 Welch Allyn 90125 For CL300 REPLACEMENT LAMP. Product No 302990125 Description WA 90125 Cogent Light 290W, Brand Welch Allyn Remarks Metal Halide CL300 Watt (W) 290, Colour Temp. (K) 5600 Life (h) 1000</p>	<p>Part Number: 90133 CL100 Replacement lamps for Cogent- Light source 270W short arc lamp 90123. Part Number: 90125 CL300 REPLACEMENT LAMP</p>	
<p>Welch Allyn 01000 Bulb for Transiluminador 43200, 43000</p>		<p>Welch Allyn 01100, Bulb for 43100 Transilluminator</p> 



<p>ML7512 WELCH ALLYN HPX Welch Allyn 2.5v HPX Gold Lamp Cartridge Assembly, Lamp cartridge WA 608125-501</p>	<p>Lamp cartridge WA 608125-501 Gold Lamp Cartridge Assembly Socket for Laryngoscope HPX 608125-50 Soquete para Laringoscópio HPX 608125-501 HPX608125-50</p>	
<p>WELCH ALLYN 03100-LED</p>	<p>WELCH ALLYN 03100-LED, LED Lamp for Traditional Diagnostic Otosopes (20000, 20200, 21700, and 25020)</p>	
<p>WELCH ALLYN 03800-LED</p>	<p>WELCH ALLYN 03800-LED, LED Lamp for PanOptic™ Ophthalmoscopes (11810 and 11820) Laryngoscope Replacement BULB</p>	
<p>WELCH ALLYN 04900-LED</p>	<p>WELCH ALLYN 04900-LED, LED Lamp for Traditional Ophthalmoscopes (11720, 11730, and 11735)</p>	
<p>WELCH ALLYN 06500-LED</p>	<p>WELCH ALLYN 06500-LED, LED Lamp for MacroView™ Diagnostic Otosopes (23810, 23820, and 23920)</p>	
<p>WELCH ALLYN 00200</p> 	<p>WELCH ALLYN 00600</p> 	<p>WELCH ALLYN 00900</p> 
<p>WELCH ALLYN 01200</p> 	<p>WELCH ALLYN 01300</p> 	<p>WELCH ALLYN 01400</p> 
<p>WELCH ALLYN 01800</p> 	<p>WELCH ALLYN 02000</p> 	<p>WELCH ALLYN 02100</p> 
<p>WELCH ALLYN 02200</p> 	<p>WELCH ALLYN 02300</p> 	<p>WELCH ALLYN 02500</p> 
<p>WELCH ALLYN 02600</p> 	<p>WELCH ALLYN 02800 Lamp W. A. 02800</p> 	<p>WELCH ALLYN 02900</p> 
<p>WELCH ALLYN 03000</p> 	<p>WELCH ALLYN 03100</p> 	<p>WELCH ALLYN 03300</p> 












WELCH ALLYN 03400 	WELCH ALLYN 03500 	WELCH ALLYN 03600 
WELCH ALLYN 03700 	WELCH ALLYN 03800 	WELCH ALLYN 03900 
WELCH ALLYN 04100 	WELCH ALLYN 04200 	WELCH ALLYN 04400 
WELCH ALLYN 04430 	WELCH ALLYN 04450 	WELCH ALLYN 04500 
WELCH ALLYN 04600 	WELCH ALLYN 04700 	WELCH ALLYN 04800 
WELCH ALLYN 04900 	WELCH ALLYN 06000 	WELCH ALLYN 06100 
WELCH ALLYN 06200 	WELCH ALLYN 06300 	
WELCH ALLYN 06400 	WELCH ALLYN HPX 06500 	WELCH ALLYN 07200 
WELCH ALLYN 07600 	WELCH ALLYN 07800 	WELCH ALLYN 08100 
WELCH ALLYN 08200 	WELCH ALLYN 08300 	WELCH ALLYN 08400 
WELCH ALLYN 08500 	WELCH ALLYN 08800 	WELCH ALLYN 09600 











WELCH ALLYN 999079-15	WELCH ALLYN 999079-15. Welch Allyn specialty bulb. PC mount. Clear. Threaded base. 0.495" O.D. x 0.942"H.	
WELCH ALLYN 999080-5	WELCH ALLYN 999080-5. Welch Allyn specialty halogen bulb. 0.45"D x 1.085"L. Alt. P/N: 7013. NSN: 6240-00-297-3914.	
WELCH ALLYN 999079-6	WELCH ALLYN 999079-6. Welch Allyn specialty bulb. PC mount. Clear. Threaded base. 0.495" O.D. x 0.942"H.	
ML4174 Standard handle, Ø 30mm x Length 150mm Made in Germany 17 Stainless steel Laryngoscope handle		
ML1691 Welch Allyn 24620 MicroTymp3 Tip Box for MicroTymp 3 Tympanometric Instrument	Olivas Welch Allyn 24620 para microtimp3 Olivas 24620	
ML7128 Welch Allyn 45-22-189 Cuff Adult Cuff Adult Welch Allyn 45-22-189 Braçadeira AP tensão Arterial Ref. 45-22-189 Welch Allyn	Welch Allyn 4522189 Cuff Adult Cuff Adult for Primus Welch Allyn. HCM-SP-14965-45-22-189. Adult arm velcro cuff for Primus 68 and Primus Stabil 3. Cuff for Welch Allyn Iso-Stabil Welch Allyn Ref: 45-22-189 +H119045-22-1891	
Riester 10502 Ri-mini otoscope, head 10502	Riester 10490 Pen-scope otoscope, head 10490	
 ML5574	 ML5543	ML5574 Cabeça Welch Allyn refª 11720 p/Oftalm. Welch Allyn Welch Allyn 11720 Coaxial Ophthalmoscope Ophthalmoscope head ref. 11720 Welch Allyn ML5543 Cabeça Welch Allyn refª 11710 p/Oftalm. Welch Allyn Ophthalmoscope head ref. 11710 Welch Allyn
ML9304 Heine Mini 3000 D-001.70.210 (MINI3000 OTO.+HANDLE+1 X TIPS) Heine Otoscope Mini 3000		
ML8584 Cabeça Heine C-002.30.100, HEINE BETA 200 Ophthalmoscope 3.5V without handle	ML8585 Cabeça Heine C-001.30.100 HEINE BETA 200 Ophthalmoscope 2.5V without handle	 BETA 200
Heine Otoscope head BETA 200 F.O. B-001.11.500 in 2.5V and B-002.11.500 in 3.5V	Heine Otoscope head K180 F.O. B-001.11.550 in 2.5V and B-002.11.550 in 3.5V	 BETA 200 K180



<p>ML9488 Mirror for headlights from Aesculap AC 461</p>	<p>ML9488 Espelho para Espelho Frontal Aesculap AC 461 Aesculap AC-461 Aesculap AC461</p>	
<p>ML4909 Hanaulux Blue 30 Hanaulux 56053010 24V 50W Black Top</p>	<p>Hanaulux blue30 56053010 Hanaulux blue3056053010 Maquett Hanaulux Blue 30 24V 50W Black Top</p>	
<p>ML4077 Hanaulux Blue 80 Hanaulux 56053026 21,5V 130W</p>	<p>Maquett Hanaulux Blue 80 21,5V 130W ref. 56053026 Black Top</p>	
<p>ML5079 Hanaulux Blue 130/90 22,8V 90W Maquett Hanaulux Blue 130-90. Hanaulux 56053198 Blue 90/130</p>	<p>ML5080 Hanaulux Blue 100 Spec: 22,8V 113W Hanaulux 56053296 Hanaulux Blue 100</p>	 <p>ML5080</p>
<p>ML1597 L7414A Lamp 12v 35w for Microscope Leica DM750 P ILT 12v 35W Formerly known as GTL</p>	<p>Medical lamp L7414A 12V 35W Leica DM 750P, Leica DM750 P, Leica DM750-P. L7414A - Quartz Halogen Lamp (Bi-Pin) 4mm Base - formerly Gilway GW L7414A</p>	
<p>ML4029 Hanaulux 3000 567901947 22,8v 40w IRC 22,8v / 40W Hanaulux 22,8v-40w XIR001, XIR 64291 Hanaulux 567901947</p>	<p>Hanaulux Hanaachrome 22,8v 40w XIR 001 ALM567901948 (ARD 567901947) IRC Infra Red coated lamp with vertical filament and G6.35 base ALM567901948, 567901947 TUV-HESSEN-EV 22,8v/40w Hanaulux Infra Red Coated Lamp 567901947</p>	
<p>ML8588 Hanaulux 2000, M115282, Hanaulux 015282, 22.8V 50W Hanaulux H015282, 15282AX HIKARI, Lamp 6419/AX2, 55240/22.8V/50W NARVA. Hanaulux 56115282 axial 22.8V 50W</p>		
<p>ML4911 Lamp 6419/AX3, 24V 35W GY 6,35</p>		
<p>ML4916 Lamp 6419/AX4 24V 75W GY 6,35, LT-03026 Microscope</p>		
<p>ML8366 Hanaulux H018366 22,8V 75W Hanaulux H018366 Hanaulux H18366 22,8V 75W, Lamp 6709/1</p>		
<p>ML8421 6,6A 100W 10X58mm, Lamp 6709/2. HIKARI M-02002 J 6.6A 100W R7S</p>		
<p>ML9842 6,6A 200W 11X60mm, Lamp 6709/3</p>		
<p>ML7009 Hanaulux H018769 22,8v 40w ML7029 Hanaulux H018566 22,8v 50w ML7030 Hanaulux H018576 24v 50w ML5691 Lamp 5924/0 12v 30w Leica DM-1000 projector shin-nippon</p>	<p>6419/6 Lamp 6419/6 12v 35w 6419/7 Lamp 6419/7 12v 50w 6419-TS 1660, 6419-TS 24v 50w 6419/5 Lamp 6419/5 24v 50w 6419/3B Lamp 6419/3B 22,8v 110w ML8224 905-24, 6419/1 24V 100W</p>	



<p>Hanaulux H016170 110v 30w BA15D 0-3709B SILVER-TOP</p>	<p>016170 Hanaulux 16170 Hanaulux/SIEMENS 016170 MICROLAMP, 16170 USHIO, 130772-H111030-B15D-40X63-110V-CL-CAL- 30W, 130772/110V/30W NARVA 8000086 USHIO, 56016170/110V/30W/B15D 56016170 HANAULUX, 302936030 , B15D-40X63-110V-30W-CLAIRE- CALOTTE/H111030/130772 B15D-40X63-11, 0376/11</p>	
<p>Hanaulux H018250 220V 50W BX22D (Frosted) Hanaulux H018251 220V 50W BX22D Hanaulux H018252 240V 50W BX22D Hanaulux H018253 240V 50W BX22D Hanaulux H016544 24V 40W BX22D SILVER-TOP</p>	<p>Hanaulux H016578 125V 50W BX22D Hanaulux H016678 (0079/1) Lamp 00936009, Lamp P-9653393 24V 50W BX22D Hanaulux H019188 240V 50W BX22D Hanaulux H019189 240V 50W BX22D SILVER-TOP</p>	
<p>Hanaulux H016679 110V 50W BX22D. SILVER-TOP</p> <p>56-016679 SIEMENS 110-050-11R CANIMPEX 110-050-11R MICROLITES</p>	<p>H384 HYBEC H016679/110V-50W-1-2-SILVER-TOP H016679, H16679 HANALUX H016679 LAMPTECH H211050-B22D-50X80-110V- 50w H211050 H56016679-HANALUX/H384 HYBEC MH016.679 BAILEY OPERATING-ROOM-LIGHT/016679 SM-16679 USHIO, SM16679 USHIO</p>	
<p>Hanaulux H016344 Hanaulux 16344 120V 30W BA15D 125V 30W BA15D</p> <p>ML1799 230V 40W BA15D Guerra 0178/6 Top Silver</p>	<p>ML1298 Apallon A.5905 Apollon 24V 50W BA15D. 0178/4. M-04031 HIKARI Apallon 5905</p> <p>ML1294 Apallon A.5904 Apollon 24V 40W BA15D</p>	
<p>Hanaulux H016164 Hanaulux 56016164 24V 25W BA15D. ML7725, Hanaulux 016372, 0376/4 24V 40W BA15D. Hanaulux 016191 12V 25W BA15D Hanaulux 016170 110V 30W BA15D SILVER-TOP</p>	<p>ML7400 A.6004 24V 30W BA15D ML7401 R.16192, 0376/9 12V 35W BA15D ML1287 24V 50W BA15D ML7714 Dr Mach 1521 24V 60W BA15D Please see also 0574/5H</p>	
<p>ML7724 Hanaulux H0185500 Hanaulux 185500 24v 40w BA15D</p> <p>56028550 Hanaulux, CL70465 CARLEY-LAMPS , 8000121 USHIO, 0574/3 , 380363 ZEISS, 131820 NARVA SILVER-TOP</p>	<p>Hanaulux 0185500 H-018550 Hanaulux H-18550 Hanaulux 56028550 Hanaulux 56028550 HANAULUX CL70465 CARLEY-LAMPS 8000121 USHIO, 0574/3 380363 ZEISS, 131820 NARVA</p>	
<p>ML9223 Dr. Fischer 00936016 24V 60w BA20D Dr. Mach Soloflex 24V 60w BA20D Dr. Mach 00019489 Pantoff Dr. Mach D8017 0970/V1 (M1521)(M 1521) R.970-BA20D SILVER-TOP</p>	<p>M-04041 24V 60W Lamp 1150/24V/60W DR-FISCHER 936003 DR-FISCHER Dr. Mach Soloflex 1150/91 302936016 SPECIAL-INDUSTRIAL-LIGHTING CL70466 CARLEY-LAMPS DF936003/24V/60W DR-FISCHER Dr. Fischer 00936016 24V 60W Lamp Dr. Fischer P-9682299 24V 60W Lamp</p>	



<p>ML1814 MARTIN M8895001 24V 50W BA20D MARTIN 546 MM88.91.604 BAILEY 24V/50W/BA20D/130593 NARVA 130593, 0970/7 Lamp. M8891604 24V 50W BA20D SILVER-TOP</p>	<p>ML1892 MARTIN M8891602 24V 25W BA20D 131821 NARVA. 0970/3B Lamp M8891602 24V 25W BA20D 86-950-00/24V/25W/BA20D/G12.5/SB MARTIN / BERCHTOLD 88-916-04- MARTIN-SURGICAL CZ-916-04 BERCHTOLD. Lamp 651E SILVER-TOP</p>	
<p>RG1225 24V 50W B22D A22 SB 24V-50W B22D A22 SB 24V50/50W BA22D 50V 50W BA22D 524V 50W BA22D DOIS FILAMENTOS 70x117mm</p>	<p>24V 50/50W BA22D 8890-1/24V/50W 8891402 24V 50/50W I/24V/50W/50W/B22D/A22/SB I/24V/50W/50W/B22D/A22/SB/CL I/24V/50W/50W/B22D/A22/SILVER- BOWL/1000H, 24V 50/50W 24-50v 50W 70x117mm SILVER-TOP TWO FILAMENT 70x117mm</p>	
<p>ML8022 8022 12v 50W BA20D LWTF5-50W-12V BA20D 12v/50W, 12v 50W BA20D ML8023 M04025 6V 50W BA20D ML8024 M04024 6V 30W BA20D Osram 8025, Guerra 3749/1 ML8025 8029 15V 60W BA20D 8000274 USHIO</p>	<p>ML8022 12-050-05, 601/67, 3405, 3880/3, 41090 67612/12V/50W, 110002, 110005 NARVA 206677, 8000270, 302843171, 031293410900 4050300206677 OSRAM 12v-50w BA20D ML8023 M04025 6V 50W BA20D, 8022 ML8024 M04024 6V 30W BA20D, Osram 8025 ML8025 8029 15V 60W BA20D 8000274 USHIO</p>	
<p>ML3296 Dr. Mach Soloflex 00019489 24V 60w BA20D Dr. FISCHER 936003 Lamp Dr. Mach D8017 Guerra 0970/V1, M1521, M 1521 R.970-BA20D SILVER-TOP</p>	<p>936003 Dr. FISCHER 24v 60w BA20D Dr. Mach 24v 60w Ba20D Dr. Mach Soloflex 0970-V1, 1150-91 DR-MACH 302936016 SPECIAL-INDUSTRIAL-LIGHTING CL70466 CARLEY-LAMPS DF936003/24V/60W DR-FISCHER DF936003 BULB-DIRECT Dr. Fischer DF936003 DR-FISCHER</p>	
<p>ML6502 24v 50w RX 24v 50w BA20D RX 2058/3 , 555023 (55.5023) 24V-50W RX M04021 DONSBULBS, I/24V/50W/BA20D/G12.5/C L/500H</p>	<p>24v-50w BA20D RX 2058-3/24V/50W, 2058-3 4316/24V/50W/BA20D/B35MM/C2, 4316 SUNG-WHEE 4416 NAM-YUNG, 0178-FA 844409005040 HIKARI</p>	
<p>ML9902 7582, 1361 I/12V/25W/BA15S/S8/AX 5782, 7582/12V/25W 7582/12V/25W/BA15S 12325/12V/25W/BA15S/S25MM/C8 12325/BA15S, 12418, 17931 NARVA 7582, 49019 NARVA Medida: 49x25mm</p>	<p>Narva 12v-25W BA15S Narva 12v 25W BA15S 12-025-04 CANIMPEX 12V/25W/BA15S/17931 NARVA 12V/25W/BA15S/AX 050-535/12V/2.25A WILD/LEITZ 050-535/12V/25W LEITZ/LEICA 050-535 WILD/LEITZ 0483, 1361, 1566-4/12V/25W I/12V/25W/BA15S/S8/AX DONSBULBS</p>	
<p>ML5665 3706, 24V 40W BA15D 130771 NARVA, M04019 HIKARI, GUERRA 0178/L1 03706 HOSOBUCHI M04019 DONSBULBS (Top Silver)</p>	<p>ML1244 3701, 24V 50W BA15D HOSOBUCHI-0-3701 130769 NARVA GUERRA 0178/L5 M04012 DONSBULBS (Top Silver)</p>	












<p>ML1299 3721, 24V 40W BA15D GUERRA 0178/3F 03721 HOSOBUCHI M-04028 HIKARI M04028 DONSBULBS (Top Silver)</p>	<p>24V 40W BA15D O-3721 HOSOBUCHI 0178-3F GUERRA M-04028 HIKARI 3721 Apollon</p>	  
<p>ML1478 220V 50W E27 0079E/220V/50W/E27 GUERRA 0079/E I/220V/50W/E27/G16.5/SB 50x82mm Frosted glass Top Silver</p>		  
<p>0574/3H 24V 40W BA15D Hanalux 018550 Halogen 1210-91-H Dr-Mach</p>	<p>0574/TN 24V 40W BA15D (Top Black) 0574/3H 24V 40W BA15D 0574/4H 24V 50W BA15D 0574/5H 24V 60W BA15D Please see also ML7714</p>	  
<p>ML4101 Lamp Hitachi DT00701 Lamp Projector Hitachi CP-RS55 50MMX50MM-BARELAMP. Hitachi CP-RS56 · Hitachi CP-RS56 + · Hitachi CP-RS57 · Hitachi CP-RX60 · Hitachi CP-RX60Z. Hitachi HSCR165H11H</p>	<p>Lamp. Projector Hitachi CP-RS55 (HSCR 165W)(P/PROJECTOR Hitachi 165W refª HSCR165H11H), (70V-165W) 50MMX50MM-BARELAMP, Hitachi DT00671 CP-X340 CP-X345 X45 X55 S55. · Hitachi CP-RX61 · Hitachi CP-RX61 + · Hitachi PJ-LC7 · HITACHI CP-RX60Z 70V 165W HSCR 165H11H. HSCR165H11H</p>	 
<p>ML4102 Lamp Projector Panasonic model PTF-300NTE , ET-LAF100, PT-FW100NTU, PT-F100NTU, PT-F100NTEA, PT-FW100NT, PT-F100U, PT-F100NT, PT-F200NTU, PT- F200U, PT-F200, PT-FW300NTU, PT- FW300U, PT-F300NTU, PT-F300U, PT- F300E, PT-F300NTE, PT-FW300E Module for Panasonic PTF-300NTE</p>	<p>Lamp Projector Panasonic model PTF-300NTE , ET-LAF100, PT-FW100NTU, PT-F100NTU, PT-F100NTEA, PT-FW100NT, PT-F100U, PT-F100NT, PT-F200NTU, PT- F200U, PT-F200, PT-FW300NTU, PT- FW300U, PT-F300NTU, PT-F300U, PT- F300E, PT-F300NTE, PT-FW300E Module for Panasonic PTF-300NT</p>	  
<p>ML4804 Epson / UHE-170E3-C 170w 80v (80v-170W) projector Epson Lamp of supply Epson UHE-170E3, UHE-170 E3-C (V13H010L42 / ELPLP42)</p>	<p>LÂMPADA PARA VIDEO PROJECTOR EPSON UHE-170E3 para projetores Epson: S4 Powerlite Epson: EMPS4, EMPS5 Epson: EMPS42, EMP-S5</p>	 
<p>Projector XL-2U SP-LAMP-018 projector DATASHOW MARCA INFOCUS MOD X2</p>	<p>Lamp for Mitsubishi projector XL-2U Colorview (lamp VLT- XL2LP) XL2U SP-LAMP-018 projector lamp with housing for Infocus X2/LPX2/X3/LPX3 projector</p>	  <p>XL-2U SP-LAMP-018</p>
<p>ML4177 Lamp. P/datashow Benq mod. PB8220 UHP Lamp. Datashow UHP 250W 1,35 Refª. 75J4002001</p>	<p>datashow Benq 250w-1,35 ref. 65.J4002.001</p>	 




<p>ML5151 GEMINI 300 EZG GEMINI 300 (EZG) GE Gemini PAR20 - Optical Projection</p>	<p>300-EZG-GEMINI 300-GEMINI GENERAL ELECTRIC 300/16A GENERAL ELECTRIC 300WATT-16-A-GEMINI CANIMPEX 300WATT-16-A-GEMINI MICROLITES 300WATT MICROLITES 35V 300W</p>	
<p>ML7496 GENERAL ELECTRIC MARC-300/16/37.5V GENERAL ELECTRIC MARC 300/16 EZM</p>	<p>37,5V 300W FCB1011 ACMI, FCB1011 CIRCON FCB1011 CIRCON-ACMI GEM-MAR-300WATT, GEMI-MARC-300WATT Dim: 75MMx69MM</p>	
<p>ML4655 GENERAL ELECTRIC EZT-MARC-350/16T GENERAL ELECTRIC MARC 350/16T EZT</p>	<p>39936 GENERAL ELECTRIC GE PAR24 - Optical Projection 45V 350W</p>	
<p>ML2124 NARVA 131838 230v 50w BA20D 131838 NARVA, H-53050 ORBITEC, H53050 ORBITEC, 53050,131.838, 230v 50w BA20D dim:80MMx50MM</p>	<p>NARVA 131838, 230v 50w BA20D 131838 NARVA, H-53050 ORBITEC, H53050 ORBITEC, H 53050 Narva 131.838, 230v 50w BA20D dim:80MMx50MM</p>	
<p>ML2758 Nikon Lamp 6v 30W BA15D for Nikon Microscope. Olympus TB-1 is reference 8-C101. Hosobuchi OP2118 6v 30W BA15D, 25x56mm</p>		
<p>ML7750 Lamp 6v 30W BA15D, 6-8V/5A-TB1 OLYMPUS, 06-030-03N MICROLITES, 06-030-03V 6V/30W/BA15D/T8/FC6, 6V/30W TOSHIBA, 6V5A-TB1 OLYMPUS 6VOLT-5AMP-9-G NIKON, 6VOLT-5AMP-TB-1 OLYMPUS, 6VOLT-5AMP Hosobuchi 6VOLT-8VOLT-5AMP-T-B-1 OLYMPUS</p>		
<p>ML7749 Lamp. Microscope WILD 6V 15W BA15D 38-00-29-9020 Zeiss, 38-29-9020 Zeiss, 166-324/6V/15W LEITZ-LEICA 166-324/6VOLT/15WATT WILD, 166-324 USHIO, 166-324 WILD 500-012/6V/15W LEITZ-LEICA, 500-012 USHIO</p>		
<p>ML4303 WILD HLX3719 12v 30w, Microscope, Colposcope, Operating Microscope, Collimators , WILD HLX3719 12v 30w JC 12V 30W SP BASE</p>		
<p>ML4305 WILD HLX4643 12V 50W, Luca 10384643 12V 50W,12V-50W HLX4643 (13931) HLX384643, WILD 50 W/384643, Leica 38643, Leica M655. Leica M841. Bulb for LEITZ 384643, 10384643, 310-752, 690H, M600 SERIES, M655, M690, M691. Wild-50 12V 50W DN-28643</p>		
<p>ML4651 General Electric G.E. 150P25/2SB 120V 150W E26 "Half Silver" 120V-150W E26 "Half Silver" General Electric 150P25/2SB Spotlight 150P25/2SB Lamp used in Amsco 360 SILVER-TOP</p>	<p>120V 150W General Electric 120v-150w E2715014/150P25/2SB-120V 15014,15014 , 249383/150P25/2SB 8000071, 04613515014 SYLVANIA 048777322390 UPC P/150/P25/2SB, P/150W/P25/2SB P25-2SB HIKAR, IP25/2SB HIKARI P150P25-2SB/120V/150W AMSCO P150P25/2SB, SL15014 SYLVANIA SM-150P25-2SB, SM-150P25/2SB</p>	
<p>We have Discontinued Bulbs</p>	<p>Para qualquer tipo de Lâmpada Agradecemos nos contactem</p>	<p>For other Lamps, please contact us</p>



<p>ML4815 HAAG STREIT HS1200559 12v 30w G6,35, BC900/BD900 12V 30W G6.35 Base 20/18mm 11,5V-2,5 A, 11,5V/2,5 Amp Haag Streit BC 900/ BD900 Haag Streit BD 900/ BC900. Tonometro de Aplanação</p>	<p>HAAG STREIT 12V 30w Halogen SLIT LAMP BULB .SLIT LAMP MODELS BC-900; BD-900; 0727 1200559 and BA-904. Is used for BC900/BD900 12V-30W. HS1200559 12v 30w, 0757-82722818 (LOT) 9H, 11,5V 2,5 Amp, 11,5V 2,5 A Haag Streit BD900, Haag Streit BC900. HAAG-STREIT 12V / 30W Lamp for Aplation Tonometer</p>	
<p>ML9012 Kavo Dental Bulb Halogen Reflector Lamp KAVO 24V 150W FDV KIT /3 STAR5</p>	<p>Dental Bulb Halogen Reflector Lamp , 24V 150W , Model: FDV , Brand: STAR5 , Model: R7S</p>	
<p>ML2919 BSS (HYS) 1089 27,7V 180W, 27,7V 180W HYS1098</p>	<p>180W 22,7V PRIMA P1350 HILL- ROM TYPE Medical Lamp Hybec. HYS 1098 BSS(HYS)1089</p>	
<p>ML5099 Lamp. Vanity 120v 25w E10 (Diâmetro casquilho 11mm) 25G161/2/W7CD, 03753, G- 16 1/2</p>	<p>25G161/2/W7CD, 03753, G-16 ½ 50x70mm</p>	
<p>ML4794 B5-010-04 Instrument Type: SURGICAL LIGHT Instrument Model 1: HA60, B5-012-01 Instrument Model 2: ARGOS FOCUSING WAND Bulb Info: HALOGEN Base: 3-PIN BASE Average Life in Hours: 100 HOURS</p>	<p>Skytron - 54-010-04 Flash Lamp - B5- 010-04 Replacement Light Bulb. ARGOS STROBE BULB SKYTRON Overall length (mm"): 39mm SKYTRON WAND STROBE LAMP, 3-PIN BASE (STB501004) B5-010-04 Skytron Strobe (Argos) Skytron Strobe</p>	
<p>ML3098 DKK 24V 50W G8 DKK 24V-50W G8 DKK-50 / B1-010-28 24V 50W G8 - Medical Lamp. Application SKYTRON for Surgical Light</p>	<p>ML3098 DKK-50-24501, Skylux-SH101028, TUBE FC6 CERAMIC-TOP, 24-050-34 CANIMPEX, 24-050-34 MICROLITES, 99-62-H-80 ORBITEC, 031293413901 H, HYBEC 2300 B1-010-28, DKK-50, B1-010-28, 6702/1, DKK24501</p>	
<p>ML5528 22,8v 50w Zeiss 3800-79-9240 MARTIN 88-904-24</p>	<p>Higuchi M-0103, Hikari M-0103 18566, H018566, SIEMENS 56-115282, MARTIN-SURGICAL 88-904-24 Zeiss 3800-79-9240, H-018566, 6419/3A</p>	
<p>ML7455 Guerra 5629/3 6v 10w Base especial, H/6V/10W/WL- 31MM-PLASTIC- SPECIAL/T2.5/CL/1000H</p>	<p>ML7455 Lamp 5629/3 6v-10w Base especial, H/6V/10W/WL-31MM- PLASTIC- SPECIAL/T2.5/CL/1000H</p>	
<p>ML7499 Osram 64635 HLX 15v 150w GZ6.35 Gold Coating on Reflector Lâmpada 64635 HLX 150W 15V GZ6,35 OSRAM</p>	<p>Osram HLX 64635 Halogen Bulb 64635HLX MR16 ,Lamp Gold Coating on Reflector, Tungsten Halogen Infrared Lamp (IR) Osram 64635HLX Xenophot Bulb, US Code 54233, Low Voltage Bulb 15v 150w GZ6,35 MR16</p>	



<p>ML7437 Karl Storz ARC60 50W 81 1310 251, 81131021, Spare lamp 50w for use With Tele Pack X, Techno Arco 60. Techno Pack II, Techno Pack X/Xe and Portable Pack Digital</p>	<p>Lâmpada para Karl Storz ARC60 50W 81131021, 811310 251, 81131021. KARL STORZ TECHNO ARC 60</p>	
<p>ML5909 PST-2100BS (03133) FT/PST2100B GRASS FTS-107 F.T.S.-0.107 PST-2100-B GRASS</p>	<p>Lâmpada para estimulador GRASS Mod.: PS33 plus (bulb PST-2100 for Photic Stimulator Grass mod. PS33 Plus) PST-2100BS, PST-2100BS (03133) FT/PST2100B GRASS, FTS-107, F.T.S.-0.107, PST-2100-B GRASS</p>	
<p>ML3305 Lâmpada estroboscópio XS 201, XS201 92x 29mm 4 Pinos</p>	<p>Stroboscope lamp XS201 Stroboscope lamp XS 201 92x 29mm 4 Pin</p>	
<p>ML7042 ASTRALUX ULTRA – VIOLETS 220V 150W 220V-150W E27 Blue Glass Carbon filament Coal Filament</p>	<p>ML7041 ASTRALUX INFRA VERMELHOS 220V 150W ASTRALUX INFRA-RED 220V-150W E27 Black Glass Coal Filament, Carbon filament (Filamento de carvão)</p>	
<p>ML1259 Lamp E10 T10x24 S. LED Red 24-28V AC/DC. Lâmpada LED 4V/E10/5SMD cor vermelha</p>		 <p>E10 T10X24 S.LED Red 24-28V AC / DC</p>
<p>ML5544 ALM 23V 100W X514SP, ECL001. 23/24V 100W X514 ALM-Prismalix ECL-001 ALM SLIDE-50067 X514SP 23V-100W ALM ECL-0001 Prismalix / Cupola 4000 / 6000 / 8000</p>	<p>100W 23V LHX 001 for ALM Prismalix ALM-ECL0001, ALM ECL-0001, ALM50067, ALM 50067, SM-5067, X514 PRISMATIC OPTICAL 150 HR ALM50067, HIKARI M-01087, HYBEC 2362 ALM-50067, Higuchi ALM5067 M-01087. 45x13cm</p>	
<p>ML6669 ALM 24V 100W X514SP ALM 7400-8400, ALM ECA002 Medical. For new ALM 50067 PRISMATIC</p>	<p>23/24V 100W , 24v-100w ALM-ECL001/002 ALM7400-8400, Chick-730453, M-01057, M01057. The same as ML5544 45x13cm</p>	
<p>ML7700 24V 100W X514 ML6668 24V 120W X514 ALM H6950 24v 120w ALM (ECA-002) 24v 120w</p>	<p>ML7700 24V-100W X514 ML6668 24V-120W H6950 X514 Sylvania 0021921 (2329) 6023/2 Sylvania OO21921. L2329</p>	
<p>ML6688 24v-100w X514SP 24v 100w X514SP 42x13cm</p>	<p>ML6686 24v 140w X514SP 24v-140w X514SP 42x13cm</p>	
<p>ML6687 24V 140W X514 Dr. Fischer 00847137 L:42xD:13mm ICL: 19,5mm</p>	<p>ML6687 24V-140W X514 Dr. Fischer 00847137 Medical Lamp 24 V 140W, A26/14 L:42xD:13mm, ICL: 19,5mm</p>	



<p>ML4436 24v 70W X-511 (104316) (52174) (64175) 24V-70W X-511 (52174 X522) H2 13895. 64175-H2</p>	<p>ML4436 (104316) (52174) (64175) 24V 70W X511 (52174 X 522 H2) Foco 153092 ANGENIEUX AX4</p>	
<p>ML4437 24v 80W X-511 ML4438 24v 100W X-511, 153091 ML4439 24v 70W X-511, 14316 (Frosted)</p>	<p>ML4437 24v-80W X-511 ML4438 24v-100W X-511, 153091 ML4439 24v-70W X-511, 14316 (Frosted) (Fosca)</p>	
<p>ML4440 24V 80W X-513 (H 247240) 24V 80W X-513</p>	<p>ML4440 H/24V/80W/X513/TUBE/2000H H247240/24V/80W/X513 H247240/24V/80W</p>	
<p>ML5694 Suporte ECL 1460 ALM Foco Op. Ref.567077999, para Lâmpadas ALM H6950 24V 120W X-514 ECL1460 (lampholder)</p>	<p>lamp holder ECL 1460 Marquette 567077999 holderlamp ECL 1460 Socket ECL 1460 lampholder Marquette 567077999 (1x 24v 120w)</p>	
<p>ML5695 Suporte ECL 1470 ALM Foco Op. Ref.567078999, para Lâmpadas ALM H6950 24V 120W X-514 ECL (Maquet 567078999) ECL1470 2x 24v 120w</p>	<p>lamp holder ECL 1470 Maquet 567078999 holderlamp ECL 1470 Socket ECL 1470 lampholder 567078999 lamp holder ECL1470 Suporte ECL1470 ALM</p>	
<p>ML5696 Suporte ECL 1480 ALM Foco Op. Ref.567079999, para Lâmpadas ALM H6950 24V 120W X-514 ECL1480 (lampholder) (Maquet 567079999)</p>	<p>lamp holder ECL 1480 Maquet 567079999 holderlamp ECL 1480 Socket ECL 1480 lampholder 567079999 3x 24v 120w</p>	
<p>ML1717 Sockets (lampholder) suporte candeeiro marca Angenieux modelo AX14, AX4 / AX10 / AX14</p>	<p>ALM162798, lamp holder for AX14. BULB HOLDER 150 W OR 100 W FOR AX4, AX10 AND AX14. ARD397900555, ALM BULB HOLDERS ARD162798 Bulb holder 150W & 100W ARD 162798</p>	
<p>ML5533 Sockets 903 (lampholder) 6,35 (44x15x20mm) Para Retroprojector</p>	<p>lampholder 6,35 (44x15x20mm) Para Retroprojector</p>	
<p>ML1733 Sockets G6,35 lampholder G6,35 Sockets P16M 217 Lampholder P16/M</p>	<p>ML1734 Sockets G6,35/G4 lampholder G6,35/G4 Sockets MR16/MR11 Socket - GX5.3/GZ4</p>	
<p>ML7143 Sockets G6,35 Lampholder G6,35 Sockets P11/C Lampholder P11/C</p>	<p>MARTIN SURGICAL B680D Used two lamps 24V 100W G6,35 (905-24). See lamps ML8224 Used for ALM Angenieux (24V 100W, G6,35).</p>	













<p>Heine X-00.88.082,</p> 	<p>Heine X.04.88.068</p> 	<p>Heine X-04.88.104</p> 	<p>Heine X-01.88.032</p> 
<p>Heine X-01.88.034</p> 	<p>Heine X-01.88.035</p> 	<p>Heine X-04.88.111</p> 	<p>Heine X-01.88.037</p> 
<p>Heine X-01.88.038</p> 	<p>Heine X-01.88.039</p> 	<p>Heine X-01.88.040</p> 	
<p>Heine X-01.88.041</p> 	<p>Heine X-01.88.042</p> 	<p>Heine X-01.88.043</p> 	
<p>Heine X-02.88.044</p> 	<p>Heine X-01.88.056</p> 	<p>Heine X-01.88.057</p> 	<p>Heine X-01.88.058</p> 
<p>Heine X-01.88.059</p> 	<p>Heine X-01.88.069</p> 	<p>Heine X-02.88.070</p> 	<p>Heine X-01.88.077</p> 
<p>Heine X-02.88.078</p> 	<p>Heine X-01.88.084</p> 	<p>Heine X-01.88.087</p> 	<p>Heine X-02.88.049</p> 
<p>Heine X-01.88.098</p> 	<p>Heine X-02.88.089</p> 	<p>Heine X-02.88.090</p> 	<p>Heine X-02.88.062</p> 
<p>Heine X-02.88.091</p> 	<p>Heine X-02.88.099 Heine X-02.88.101</p> 	<p>Heine X-01.88.105</p> 	<p>Heine X-01.88.110</p> 















<p>ML4698 Heine Y-096.15.108 Projector XENON 1000/100 W</p>		<p>Projector XENON 1000/100 W Spare bulb for Xenon 1000 F.O. Projector Projektor XENON 1000/100 W</p>
<p>ML1155 Topcon 40524-18100. TP4052418100, 0016197. Topcon TRC-50VT/V Xenon Flash Tube. TRC-50-F/FT</p>	<p>Retinografo Topcon Flash 40524-18100. TP4052418100 , 0016197, Fundus Camera TRC50F/TRC50FT/TRC50V/TRC50VT, Ophthalmic Industry TRC50F/TRC50FC</p>	
<p>ML2251 Topcon 40526-18300 TRC-45AN Xenon Flash Bulb. 4052618300 TOPCON. TRC-45AN Xenon Flash Bulb. Xenon Flash Bulb (SKU B18300)</p>	<p>ML1449 CTS-10854, PerkinElmer Old B0160917 PerkinElmer, L6022728 PerkinElmer. Perkin Helmer B016-0917 / L6022728 PerkinElmer. Product Description: Deuterium Lamp (1500 hr). Model(s): 5, 7, 9, 10, 11, 12, 14, 15, 16, 17, 18, 19, 20, 25, 35, 40, 45, 800, 900, Bio 19, Lambda 2. For Perkin Elmer Lambda 25. Made by Heraeus Noblelight</p>	 <p>ML2251 ML1449</p>
<p>ML5005 DEUTERIUM 8453 Deuterium Lamp Type HP-404 for Agilent and Bios Detectors</p>		<p>Agilent (Hewlett Packard) - 1040 series I DAD (79880A), 1040 series II DAD (79883A), 1050 (79854A), 1050 DAD (G1306A), 1050 VWD (79853C), 1090 series I DAD (79880A), 1090 series II DAD (79883A), 1100 VWD (G1314A), 1100 DAD (G1315A), 1100 MWD (G1365A), 1200 DAD (G1315B), 1200 MWD (G1365B), 1200 VWD (G1314B), 1200 VWD-SL (G1314C), 8453A, CE G 1600 capillary electrophoresis system, HP 3D CE (G1600A CE) Bios Analytical - DAD 1100</p>
<p>ML1489 Deuterium lamp J15 80017500. Deuterium Lamp Type J15 (80017500) from Heraeus Noblelight. 80017500-AZ- 80017500-AZ</p>	<p>ML1451 Perkin Elmer B011-4620, Tungsten Lamp (Visible Region), for PerkinElmer Lambda Spectrometer Models. PerkinElmer B0114620, Tungsten Lamp 30W 10.8V G5.3 Adapter (P/N B0166354) may be required for older models of Lambda 5, 7 and 9. For Perkin Elmer Lambda 25</p>	 <p>ML1489 ML1451</p>
<p>ML4426 Lamp L-280 ML5004 DEUTERIUM D946 LH-51</p>	<p>ML4426 Hitachi 100-10, 100-20, 100-30, 100-40, 100-60, 100-80 ML5004 Hitachi U-1100, U-1500, U-1800, U-2000, U-2001, U-2010, U-2800, U-3000, U-3010, U-3110, U-3210, U-3310, U-3400, U-3410, U-3500, U-3501, U-4100</p>	 <p>ML4426 ML5004</p>
<p>ML5881 Deuterium Lamp Type SD 3651-02 J from Heraeus Noblelight Lâmpada Deuterium p/Espectrofotometro UNICAM UV2-100. Spectral distribution: 180 nm - 600 nm</p>	<p>Deuterium lamp suitable for: Thermo - Aquamate, Evolution 600, UV1, UV1-100, UV2, UV3, UV4 Unicam - Helios 2, Alpha, Beta, Delta, Epsilon, Gamma, UV1, UV2, UV3, UV4, UV300, UV500, UV550. Deuterium lamp for Unicam UV2 Deuterium lamp Type: 80013730 J 64 date code: 1314, J 64 - 800013730</p>	
<p>ML4679 Xenon Flash Bulb Topcon TRC- NW3/NW5/NW6/NW6S.</p>	<p>ML4677 Xenon Flash Bulb Topcon TRC-NW100/200/300 Xenon Lamp</p>	 <p>ML4679 ML4677</p>











<p>ML1246 Philips Master TL5 HE Xtra Eco 13=14W/840 1SL MASTER TL5 HE Xtra Eco 13=14W/840 1SL/40 Length: 563.2 (max) mm Diameter: 17 (max) mm</p>		
<p>ML4448 DR. Fischer 12V 35W-PR 305 Marine Lamps</p>	<p>DR .fischer 1 12 35W-PR 305 DR. Fischer 12V 35W P30S Lamp 304847085</p>	
<p>ML7754 Kowa SL-14, Kowa SL-15 Hand-Held Slit Lamp Bulb Kowa SL-15 and SL-15 12V 120W</p>	<p>12v 120W Bulb - Kowa SL-14 and SL-15 Hand Held Slit Lamp 12v 120W: Veatch Ophthalmic Instruments SLK004-500 and SLK005-500</p>	
<p>ML8007 Chromophare D510 D540/D660 CZ-909-22-04 22,8v-150w CZ 909-22-04 22,8v 150w CZ 909-22 / 22.8V/150W, D510, D540, D660</p>	<p>Chromophare D-510 22,8v 150w 22,8v-150w 22,8v/150w Lamp p/Berchtold OT light D530 Cod. p/Osram XIR lamp 99-612C1 XIR UVS 150W 22.8V G6.35</p>	
<p>ML8080 Chromophare D530 & D530 Plus CZ-908-22-04 22,8v-150w CZ 908-22-04 22,8v 150w</p>	<p>Chromophare D530 & D530 Plus 22,8v 150w Lamp p/ Berchtold OT light D-650 plus e p/ D-530 plus CZ 908-22 (908/22) 22815022</p>	
<p>ML8081 Chromophare C572 Chromophare (908/24) CZ-908-24 22,8/24V 150W Berchtold Chromophare C572 CZ 908-24 22,8/24V-150W ML4409 Chromophare CZ 908-24 22,8V 120W</p>	<p>ML8081 Chromophare C572 CZ-908-24 22,8/24V-150W 22,8V 150W, 22,8V-150W Angenieux 163508 CZ 906-24, CZ 906-24-04 ML4409 Berchtold Chromophare CZ 908-24 22,8v 120w 22,8V-120W</p>	
<p>ML8079 Chromophare D650 CZ-907-22-04 22.8V-250W CZ 907-22-04 22.8V 250W</p>	<p>Chromophare D650 CZ 907-22-04 22.8V 250W Lamp p/ Berchtold OT light D-650</p>	
<p>ML7786 Chromophare C450 CZ-905-22 22.8V 110W CZ 905-22 22.8V-110W</p>	<p>Berchtold Chromophare CZ-905-22 Berchtold Chromophare C452</p>	
<p>ML0555 Button on/off, Membrana interruptor para Berchtold Chromophare C452, Membrana interruptor 47091</p>	<p>Membrane switch for Berchtold Chromophare C 452 Membrane switch 47091</p>	



<p>ML5028 Handle for D650 Plus. Handle for Pantof D650. Chromophare D 650 plus Chromophare D650 plus</p>	<p>Punho para foco D 650 Plus. Pantof BERCHTOLD D650 Plus Handle CZ-499-06-04. Handle for D650 Plus CZ499-06</p>	
<p>ML5032 Handle S601 surgical light. Handle for surgical light Surgiris Model S601</p>	<p>S601 highlights Handle Surgiris surgical light S601 STERILIZABLE HANDLE FOR SURGIRIS</p>	
<p>ML0812 OT Dr. Mach Model M3SC 72000183 repair kit for OT light, consisting of one cluster electronic board, 3 new LEDs and the wiring of the OT light. Mach LED -3</p>		
<p>RL4469 Kavo 1601/653 (7404211) (07404211)</p>	<p>Kavo 1601/653 7404211 24V 150W medical replacement lamp K214</p>	
<p>JR532A Lamp DEL Microscope 250v 15w BA22D</p>	<p>Lâmpada para Microscópio, 250v-15w BA22D (medida: 60 x 30 mm). Lamp DEL</p>	
<p>ML7007 A47809 6V 20W (Electrical Spezial-lampen A47809) 6v-20w Falange 6v-20w G4 (A47809). M01130 6V 20W HIKARI , S.4140 Microscópio Sbisa Microscope Sbisa</p>	<p>6320-1/6V/20W, 6320-1 H/6V/20W/G4- SPECIAL/100H/M01130 HIKARI, H/6V/20W/G4-SPECIAL/CL/100H/M01130 HIKARI, H/6V/20W/G4-SPECIAL/CL/100H M01130 HIKARI, M-01130 HIKARI, M01130/6V/20W HIKARI, M01130 HIKARI, S- 0.4140, S-4140. Sbisa 3399</p>	
<p>ML7101 Lamp 000000-0234-793, Bulb 12V 35W Halogen - Axiovert 40/200 and Axio Observer</p>	<p>ML7097 6V 20W Halogen, Colposcope KSK 111, 302302-9090 Lamp. p/Colposcópico KSK111</p>	
<p>ML4151 NARVA SH5914 12v 23w G4 30x19,5mm HLL-23W-12V-G4 Donsbulbs, HLL-23W-12V-G4</p>	<p>Lamp Narva 55914 12v 23w G4 Narva SH5914 12v 23w G4 55914 NARVA, H/12V/23W/G4 H/12V/23W/G4/100H H/12V/23W/G4/CAPSULE JC-12V-23W/M</p>	
<p>ML9921 Bulb 6v 20w 02300 for Mediworks S350S Slit lamp Bulb for Tower typ slip lamp 6v 20w (Ceramic lamp) Slit Lamp Microscope S360S/S350S Slit Lamp Microscope , Only Bulbs</p>		
<p>ML4094 Lamp for Karl STORZ TELECAM SL II (TELECAM SLII)</p>	<p>KARL STORZ TELECAM® SL II Only Spare Lamp</p>	



<p>ML9008 HOSOBUCHI P-7002 P-7002 HOSOBUCHI 100v 100w BA15D 100v-100w BA15D 70x25mm</p>	<p>Hosobuchi P7002 Hosobuchi P-7002 Hosobuchi 100v 100w BA15D 100v-100w BA15D (Vidro Liso transparente)</p>	
<p>ML8029 T22x63 240V 20W E14 mat T22x63 240V-20W E14 mat 6323520. T22x63 240V 20W E14 frosted T22 x 63 240V 20W E14 mat FROSTED</p>	<p>T22x63 235V 20W E14 mat 6323520. T22x63 240V 20W E14 microscope KYOWO TOKYO - 220/240V 20W, 220/-40V 20W (Internal lamp) FROSTED</p>	
<p>ML9910 Zeiss 30-10-62-9901 6v 10w PG22 30-10-62-9901 ZEISS, 30-10- 80 ZEISS, 38-01-20-7030-D ZEISS, 38-01-20-7030-SLIT- LAMP ZEISS, 38-01-20-7030 38-01-207030 ZEISS, ML9906 Zeiss 0143-124-000 6v 20W PG22 64251HLX, 25198 048777332122 UPC 844409004203-M01007 30-13-87 ZEISS, 5628/WO</p>		<p>ML9910, ZEISS 30-10-62-9901 6v-10w PG22 130154-H164222-PG22-9X44-6V-3000K-10W 130154/6V/10W/PG22D, 5628/W1 130154/H164222/PG22-9X44-6V-3000K-10W 301080 ZEISS, 64222, 327273, 380120-7030-D ZEISS 910849-1010 ZEISS, 942202, 6470165 8000214 USHIO, 3010629901 ZEISS 3801207030-D, 3801207030 ZEISS Microscope Light Bulb, Halogen. Speciality Bulb, Zeiss 3801-20- 7030, Zeiss 3801-20-7030D, Zeiss 38-01-20-7030-D, Zeiss 10TS, Zeiss 10TF, Zeiss 10STG, Zeiss 10-S, Ushio SM64222, ML9906, Zeiss 0143-124-000 6v 20W PG22 20SL ZEISS, 30-13-87 ZEISS. 0143-124-000 ZEISS, 143-124 ALCON, 3013-87 ZEISS, 04613554021 SYLVANIA, 5628/WO microscope Slit Lamp 6V 20W PG22</p>
<p>ML9991 Zeiss 38-01-20-7040 12v 30W PG22 CP 1030, CP1030, Auto Chart P projector, 5628/W2 PG22, 12V 30W Microscope Slit Lamp. SL-30. 38-01-20-7040-D ZEISS · 1045 · S7 KOWA · 30TF ZEISS · KOWA · TS 1600 NIDEK · STG/7 KOWA · 1080/12V/30W MARCO · 34338-4 · TF1600 NIDEK · SL-1600 · TC2ZE1027 · ASL-7-A-26 Kowa · S2ZE1027 · YAG-1379 · ALCON 2500 LASER Microscope ALCON 3000LE Las Coherent 7970YAG</p>		<p>Zeiss 38-01-20-7040 12v 30W PG22 79761/12V/30W NIKON · 30-STG-M Zeiss · TS-2-ZE-1027 · 30/TF Zeiss · TF 7 KOWA · S/2/ZE/1027 · TS1600 NIDEK · TF 1600 NIDEK · MAIN/GXX/IV MARCO · 3801 20-7040 Zeiss · STG 7 KOWA · SL-2-ZE-1027 · ICM-405- SCOPE Zeiss · ASL7A26 KOWA · 380179/9540 Zeiss · STG1600 NIDEK · COMPACT S 7 KOWA · 30 TF M ZEISS · 380120/7040D ZEISS · STG 1600 NIDEK · 30 STG M ZEISS · COMPACTC7 KOWA · ZOOMPHOTOF2 · TF/1600 NIDEK · G-IV/SLIT-LAMP MARCO · YAG/1379 79761 OSRAM · NSP 1 NIKON · 1080 MICROLITES · COMPACT-SL-7 · FS 2 NIKON · SM/64260 USHIO · 30/TC/M ZEISS · FS/2 NIKON · 30SLM ZEISS · 30CVG ZEISS · TF7 KOWA · 30/S ZEISS · COMPACTTF7 KOWA · TC 7 KOWA · FS/3 NIKON · FS-2 NIKON · MR 1045 MARCO · 34338-4 · TF1600 NIDEK · SL-1600 NIDEK · ASL-7-A-26 KOWA · S2ZE1027 · YAG-1379 · TF2ZE1027 · MR1045 MARCO · 30/ CVG/M ZEISS · TS/1600 NIDEK · TC/7 KOWA · Projector Optotipos. 380120 7040, 12V 30W</p>
<p>ML5141 Zeiss 234755 ZEISS 12v 50W PG22 FF450 Zeiss, CL70098 CARLEY 64609 OSRAM, 64609 SYLVANIA 64609 USHIO, 5628-1 Guerra. Zeiss RC-XV2, 131184 NARVA, 8000217 USHIO, SM64609 USHIO</p>		<p>ML5148 Zeiss 38-00-75-1020 Zeiss 12v 100W PG22 38-00-75 Zeiss. 38-75-1020 Zeiss, 3800-75-1020 Zeiss, 3875- 1020 Zeiss, 5628-2 GUERRA, 5628/2 GUERRA, 54038 OSRAM, 54038 SYLVANIA, 64626-HLX OSRAM, 64626-HLX SYLVANIA, 64626-HLX WOTAN, AU500 Olympus, BW-EHE. CL70037 CARLEY-LAMPS</p>
<p>ML1472 EV2-000021 Light source, 300 W Xenon. Light Source 300W Xenon EV2-000021 module lamp. EV2-000021 300W XENON Light source For 3D System CF</p>	<p>Aesculap EinsteinVision Light Source 300W Xenon EV2-000021 module lamp. Lamp module EV-000039 for EinsteinVision Light source EV2- 000021. EV2-000021 300W XENON Light source For 3D System CF</p>	
<p>ML9482 Lamp Y 1958, lamp module for Haag-Streit Hi-R 700. Microscope Moller - Wedel (Haag-Streit Surgical) Model: XL2-23 Lamp module Y1958</p>		



<p>ML6677 Olympus CLV-160 Light Source Powerful 300W xenon lamp for Olympus color CCD videoscope system</p>	<p>Powerful 300W xenon lamp for Olympus color CCD videoscope system Only Xenon Spare lamp</p>	
<p>ML9880 300W xenon lamp for the GIMMI Alpha Light 300W Fonte de Luz Gimmi modelo Alpha Xenon Light 300 W, S.5110.05</p>	<p>Light source ALPHAXENON LIGHT 300 W ALPHA XENON LIGHT 300 W 300W xenon lamp for the Alphaxenon Light 300W Only Xenon Spare Lamp</p>	
<p>ML2008 Fonte Luz Storz 20133101-1 Cold Light Fountain XENON 300 20133101-1 201330 28. ArthrexSCB300 Only Xenon Spare Lamp</p>	<p>ML5500 300W Xenon lamp for Storz 201332 20, Storz 20133220 Xenon 300W para equipamento Storz 201332 20 (Storz 20133220) Only Xenon Spare Lamp</p>	
<p>ML6601 Storz Xenon Nova 300 STORZ 20134020 201340 20 Cold Light Storz Xenon Nova 300. STORZ 20131520 Only Xenon Spare lamp</p>	<p>ML4498 Storz Xenon Nova 175 Storz 20131520 (PN 20131580) (Xenon Nova 20131520) Xenon 175W p/ Storz 20131520 (Storz 201315 20) Only Xenon Spare lamp</p>	
<p>ML7756 Fujinon XL-4400 ref.^a LMP-002 Fujinon XL-4450 Only Xenon Spare Lamp</p>	<p>ML6788 Original module LMP-002 (complete with lamp). Fujinon XL-4400, Ref. LMP-002. Fujinon XL-4450. Replacement for FUJINON LMP-002 LAMP and Heat Sinks Light Bulb</p>	 <p style="text-align: center;">ML6788</p>
<p>ML5779 Fujinon XL-401 300W CERAMIC XENON LIGHT BULB LAMP LÂMPADA APARELHO LUZ FRIA Fujinon XL-401 Fujinon XL-402.</p>	<p>300W CERAMIC XENON LIGHT BULB LAMP Fujinon XL-401, Fujinon XL-402. XL-401/402 Light Source and VP-401/402 Processor modules Only Spare Lamp</p>	
<p>ML8161 Medical Lamps Light Source/Operating Microscope Olympus MAJ-1817. MD MAJ-1817 Olympus. EVIS EXERA III XENON LIGHT SOURCE OLYMPUS CLV-190. UXR-300SLY</p>	<p>Original Olympus MAJ-1817 300W Xenon Lamp for Olympus Light Sources CLV 190 and Evis Lucera Elite CLV-290SL. UXR-300 SLY Olympus MAJ-1817 OEM Xenon Bulb Light Source CLV-190. CL1854 Xenon Bulb Light Source CLV-190, CLV190, CLV 190</p>	
<p>ML6806 Pentax MSCR70-50E (MCR70K70H) MROL-M2 MSCR70-50E / OL-M2 70W METAL HALIDE LAMP ASSY NO TIMER, FOR PENTAX. EPK-700</p>	<p>Pentax MSCR70-50E/OL-M2, MODEL: PENTAX EPK-700, 70W METAL HALIDE LAMP ASSY NO TIMER, FOR PENTAX EPK-700 PENTAX EPK700, PENTAX EPK 700. Lamp Pentax MROL-M2 Only Spare Lamp</p>	
<p>ML9424 Pentax EPK-100P Pentax EPK-1000 Lamp ROL-X20 Pentax</p>	<p>ERLO/X20, ROL-X20-1, OL-X20 Pentax EPK100P, Pentax EPK 100P Pentax EPK1000, Pentax EPK 1000 00000ROL-X20 (21, 2995Y8 (30) 1) Only Spare Lamp (Complete kit)</p>	
<p>ML4197 KARL STORZ 20132026 175w Lamp. Xenon Nova 175 Storz Y1441, Storz Xenon nova 20131520</p>	<p>Storz Lamp Y1441 (20132026) Storz Xenon nova 20131520. Replacement for KARL STORZ Y1441 replacement light bulb lamp. KARL STORZ MODEL 20131520 / REF.Y1441</p>	



<p>ML7441 PE175BF Perkin Elmer Xenon CX-2175, ORC175F, ICL 175 PE175 BF (Xenon Storz Y1441 p/Storz Xenon nova Mod. 175 201315 20 (Ref. 640) (Ref. 9904) (ref.1106) Olympus XLS (175W) Y1176 or CLV-A (175W) 20131520 Nova, 20132020, 615B, 615C, 615XE, 615XES Storz. Power: 175 Watts, Current: 14 amps (DC), Operating Voltage: 12.5 volts (DC), Ignition Voltage: 23 kilovolts, Temperature: 150 C, Lifetime: 1000 hours (600 hours warranty).</p>		<p>610 XES STORZ, 610XES STORZ, 3110-32-9010 ZEISS, 5121 WOLF, 5121.012-AUTO-LP-LIGHT-SOURCE-BULB-175W WOLF, 5121.012 WOLF, 7000/175W LUXTEC, CORPORATION, 9100/LUXTEC INTERLIGHT, 9100 LUXTEC-CORPORATION, 9175-175 LUXTEC-CORPORATION, 9175-175 LUXTECCORPORATION, 9175-175W LUXTEC-CORPORATION, 9175-175W LUXTECCORPORATION, 9175-LIGHT-SOURCE-BULB-175W LUXTEC, CORPORATION, 9175-X, 9175-XE, 9175/175 LUXTEC-CORPORATION, 9175/175 LUXTECCORPORATION, 9175/175W LUXTEC-CORPORATION, 9175/175W LUXTECCORPORATION, 9175/LUXTEC INTERLIGHT, 9175/X, 9175/XE, 9175 175 LUXTEC-CORPORATION, 9175 175, LUXTEC-CORPORATION, 9175 175W LUXTECCORPORATION, 9175 LUXTEC-CORPORATION, 9175 X, 9175 XE, 9175X, 9175XE, 25180, 57015, 58122/175W OLYMPUS, 130835/175W ORBITEC. 130835 , 201315-20-XENON-NOVA STORZ</p>
<p>RL4466 Xenon CX-2300, Perkin Elmer PE300BF, PE 300BF, ORC300F. Philips Medical Ceramic Xenon 300w, Modelo nº CEX-300BF (Philips CEX-300BF, Olympus MD631, (ORC1064S) , CL300BF (USHIO UXR-300 BU/T) CLV-S30, CLV-S40, CLV-U20, CLV-U40, CLV-U60, CLVF10, CLVU40, CLVU20, CLV-S, CLVSBO Olympus, CLV 420-U40 CLV10, CLV-020, CLV-1, CLV-10R, CLV-A 300W, CLV-S20 (cod. 28148) CLV-180 Ceramic Xenon Lamp) Timer TM 1167, TM 1739 , TM2. Power: 300 Watts, Current: 21 amps (DC) Voltage: 14 volts (DC), Ignition Voltage: 23 kilovolts . Temp: 150 C Lifetime: 1000 hours (600 hours warranty) Visible output: 5000 Limens, Colour Temp:5600 Kelvin</p>		<p>PENTAX 3300 CL300BF PENTAX 500A CL300BF PENTAX 86192 CL300BF PENTAX 86194 CL300BF PENTAX 86357 CL300BF PENTAX 86514 CL300BF PENTAX 86515 CL300BF PENTAX 908506 CL300BF PENTAX EMP300 CL300BF PENTAX EMP3000 CL300BF PENTAX EPM 3000 CL300BF PENTAX EPM 3500 CL300BF PENTAX EPM-3300 CL300BF PENTAX EPM3500 CL300BF PENTAX LIGHT SOURCE #3000 300W: CL300BF PENTAX LX-3000 CL300BF PENTAX LX-500A CL300BF PENTAX LX300 CL300BF PENTAX LX3000 CL300BF PENTAX OL-4 CL300BF PENTAX OL-X17 CL300BF PENTAX OL-X4 CL300BF PENTAX OL-X6 CL300BF PENTAX OL-X8 CL300BF PENTAX OLX4 CL300BF PENTAX VIDEO PROCESSOR CL300BF PENTAX Y1089A, PENTAX Y1124 CL300BF PENTAX Y1251 CL300BF PERKIN ELMER LX4400 CL300BF PERKIN ELMER LXC4400 CL300BF PILLING 1300 CL300BF PILLING 300W XENON CL300BF PILLING 521393 CL300BF PILLING 521448 CL300BF PILLING 528300 CL300BF PRODUCTS FOR MEDICINE 846519, CL300BF PRODUCTS FOR MEDICINE XL300 LAMP ONLY: CL300BF QED 3010 300W LIGHT SOURCE CL300BF RUSCH 1300 CL300BF RUSCH 300W XENON, CLV-180 Olympus. 300W (14V 21A), ORC1052, MB-361 Y-1052S. 30765R-6 (P/CLVS20) LAMP H001- 300W H.L.S.- HIGH PRESSURE 3567658 Para Aparelho Marca GMB Modelo ML-L.. Olympus MD-631, Olympus MD 631. Diameter: 32.96mm/1.298", Nominal Length: 42.34mm/1.667", PERKIN ELMER / Excellitas / Luxtel. Olympus CLV-180, Olympus CLV180, Olympus CLV 180. Olympus CLV 520, Olympus CLV 520.Olympus MD 631</p>
<p>ML22446 PE300BF Perkin Elmer, Excellitas, Xenon CX-2300, ORC300F, MD 631 14v-300w. lamp Philips Medical Ceramic Xenon 300w, CEX-300BF (Philips CEX-300BF, MD631) (ORC1064S) CL300BF (USHIO UXR-300 BU/T) CLV-S30, CLV-S40, CLV-U20, CLV-U40, CLV-U60, CLVF10, CLVU40, CLVU20, CLV-S20 , CLVSBO Olympus,</p>	<p>CLV 420-U40 CLV10, CLV-020, CLV-1, CLV-10R, CLV-A 300W, CLV-S (cod. 28148) CLV-180 Olympus Ceramic Xenon Lamp. Timer TM 1167, TM 1739, TM2, LUXTEL CL 300BF. ILC Cermax. Olympus CV-180 Evis Exera II Endoscopy System includes CV-180 Video Processor, CLV-180 Xenon Light Source, Keyboard, Interconnect Cables, Mobile Cart, Monitor. Base: 5 Hole, Window Diameter: 25.4mm/1", Diameter: 32.96mm/1.298", Nominal Length: 42.34mm/1.667</p>	
<p>ML3145 Pentax ROL-X28 300w Xenon lamp for Pentax EPK-i5000 and Pentax EPK-i5010 Xenon lamp Y1911 ML3150 Lamp module ROL-X29 for Pentax EPK-i7000</p>	 <p>ML3145</p>	<p>ML3151 ME300BF 300w xenon lamp for use with Smith & Nephew 500XL light source & Leica M525, Leica F40, Leica F50, 7000i. For Pentax EPK - i7000, Surgical Microscope</p> 
<p>ML3151A Xenon 300W Complete module Lamp Y1907, for 500XL Xenon Light Source, brand Smiths&Nephew</p>	<p>REPLACEMENT BULB FOR SMITH & NEPHEW 500 XL LAMP & MODULE, Y1907 LAMP & MODULE. SMITH & NEPHEW 500 XL Endoscope Light Sources Original Module Y1907 300W xenon arc bulb heat sinker housing set</p>	



<p>ML2253 Lamp Cermax for Electronic Endoscope OL-X27 (ROL-X27) with heatsink for Pentax EPK-i5000 , i5010 (Cermax) / Luxtel (Ceralux).</p>	<p>Lamp OL-X27 (ROL-X27, OL-X28) for Pentax EPK-i5000 endoscopes. Included xenon lamp 300W and heatsink (cartridge). Fully assembled and ready for installation endoscopic equipment. Original lamps are made in USA under brands Excelitas Technologies (Cermax) / Luxtel (Cer</p>	
<p>ML3152 Xenon 300W Complete module Lamp Y1915, for microscope Leica M525 F40, Leica F40, ZEISS NC-4, Excelitas J-2018, BT J-2018 for Leica Y1915 300</p>	<p>Original Module Y1915 300W xenon arc bulb heat sinker housing set LEICA M525 F40, ZEISS NC-4 operation surgical microscope The ASL-Y1915 replacement xenon lamp module for Leica M525 F40. Replacement Lamp Module for Leica Y1915 Microscope Leica F40</p>	
<p>ML4503 Zeiss 3049539030 Superlux 300/301 Complete module Lamp CL1426, for microscope Leica MS2, Leica M720 OH5 Lamp Module CL 1426.</p>	<p>304953-9030, LX300/F, CL-1426 Microscópio Leica MS2 Leica OH5, OH4, OH3, OHS1, MS3, MS2, MS1, Leica F40, F50. Luxtel 1426 Superlux 300, Superlux 301 Linvatec/ConMed 8430, 8180, 8170. Luxtel module Lamp 304953-9030-000</p>	
<p>ML1207 STRYKER X6000 ENDOSCOPE NEW LIGHT BULB, PN 220-185-300</p>	<p>Complete module Lamp 220-185-300. STRYKER X6000 Complete module Lamp</p>	
<p>ML7181 300W Xenon Lamp Bulb Module for Stryker X7000 X8000 Endoscopy Light Source Lab</p>	<p>Bulb Module for Stryker X7000 X8000 Endoscopy Light Source Lab. Complete module Lamp</p>	
<p>ML7178 Stryker lamp Module X8000 Ref. 220-201-000. Xenon Lamp 300W Original Stryker Bulb Module 220-201-000</p>	<p>Stryker 220-201-000. Stryker Light Source Xenon Bulb Module X8000 220-201-000. iButton data logger found on the Stryker X8000 220-201-000. Stryker X8000 Complete Module 300W XENON Lamp, housing and with I-BUTTON. Stryker Endoscopy 220-201-000</p>	
<p>ML7179 Stryker Lamp 220-191-000, Stryker 220-191-000 300W Xenon Lamp Module X7000/X8000 Endoscopy Light Source Lab JV-ANT-C-STY220</p>	<p>ML7172 Y1830 300W Xenon Bulb X7000 Endoscopy Light Source. Y1830 220-190-300 PE300C-10F 300W Xenon Bulb for Stryker X7000. Only Lamp. PE300C-10F - 300W Elliptical. See Also ML7173</p>	 <p style="text-align: center;">ML7179 ML7172</p>
<p>ML2297 CL1576 Replacement Lamp for Pentax EPK-i, Y1882. Lamp for EPK-I, Lamp CL1576 Y1882-ASL-Y1882-0LX25. Luxtel CL1576 Pentax EPK-I. Pentax KBM P CL1576. Pentax Y1882. RLT2108. CL-1576, EPK I. ROL-X25-1 Only the lamp Y1882, CL1576</p>		
<p>ML2212 Lamp for EPK-I, lamp Y1882-ASL-Y1882-0LX25. Luxtel CL1576 Pentax EPK-I. Pentax KBM P CL1576. Pentax Y1882 CL1576. RLT2108. CL-1576, EPK I CL1576 Replacement Lamp for Pentax EPK-i, Y1882. ROL-X25-1 Complete lamp module ROL-X25-1 for Pentax EPK-i</p>	 <p style="text-align: center;">ML2212</p>	











<p>ML2290 CL1495 Replacement Stryker Xenon Lamp X8000. This replacement lamp is designed for Stryker 8000 in Luxtel's Assembly CL1567 only. This lamp will not fit into a Stryker Heat Sink, only into a Luxtel Heat Sink Assembly CL1567. CL1567: With CL1495 lamp</p>		 <p>CL1567</p>
<p>ML7173 Y2900 300w xenon lamp Stryker X7000/X8000 endoscopy cold light source endoscopic Y2900 VAC300-F13-B-MB. Only Lamp. See Also ML7170 and ML7172</p>	<p>If the light source is not first generation the lamp is a screw based type. Y2900 300W Lamp for Stryker X8000, X7000. But Please take note that as long as the iButton timer on the lamp housing is not replaced the light source will continue to ask you to replace the lamp. Old generation</p>	
<p>ML7170 Complete module including housing and iButton timer. Stryker Light Source Xenon Bulb Module X8000</p>	<p>Stryker Light Source Xenon Bulb Module X8000. iButton data logger found on the Stryker X8000. Source light Stryker X8000.</p>	
<p>ML6914 Y1152 304953-9030-000 Microscope Carl Zeiss NC-4, Zeiss NC 4. Only Bulb for module</p>	<p>ML4635 Original Zeiss module Zeiss module 304953-9030-000 Microscope Carl Zeiss NC4. Zeiss NC - 4 Microscope</p>	 <p>ML4635</p>
<p>ML1707 Fujinon/Hoya EPX-2200, EPX-2500, EPX2200, EPX 2200, EPX 2500 Fujinon EPX-2500</p>	<p>Pentax LMP-XE150 or LMP-001 (that is complete module) The lamp is PE150AF / CL150AF EPX-2500 processor/light source. CL150AF REPLACEMENT LAMP FOR FUJINON EPX2200, EPX2500</p>	
<p>ML7087 REPLACEMENT BULB FOR ZEISS 304949-9001-000. Zeiss OPMI Vario 700. 300W Xenon lamp for Zeiss Vario 700 Microscope ONLY LAMP</p>		<p>ML7087A, Complete module including lamp.</p> 
<p>ML7987 Carl Zeiss Xenon Lamp Module for S8 Ref 304977-9037 Surgical Microscope Lamp Module 304977-9037 Only replacement lamps</p>		
<p>ML4450 Lamp module Richard Wolf Mod. 5131.001, Mod.2431.101 ML4455 Only Xenon lamp</p>	<p>ML4450 High-Performance Light Source 300 W, Xenon, AUTO-LP 5131. incl. 300W xenon lamp module (Xenon lamp module 2431.101)</p>	
<p>ML4150 Lamp module Richard Wolf Mod. 2431.121. Endolight X300. Richard Wolf Ref. 2431.121</p>	<p>ML5105 Xenon 400W for Microscópio Leica. Leica M720 OH5 Lamp Module CL1585. Leica M720 F60, Leica M720F60</p>	 <p>ML5105</p>
<p>ML4140 Xenon Lamp Excelitas Cermex PE300BUV</p>	<p>Excelitas Cermex® PE300BUV 300W. Wattage (W): 300. Voltage (V):14th. Color temperature (K):5900. Lifetime (h):1000</p>	











<p>ML5255 Xenon CX-150AF Perkin Elmer PE150BF 150w 11,70V-150w 12,5A</p>	<p>ML1862 Storz Y1100-611 300W xenon lamp for 20133620 20133120 20133020 endoscope light source Karl Storz 300W xenon lamp</p>	 <p>ML1862</p>
<p>RT4457 Zeiss Superlux 40 Replacement Lamp Module Complete Assembly(310752)</p>	<p>Zeiss Superlux 40</p>	
<p>ML7090 Martin lamp 86-970-00-04 for KS Martin ML702 discontinued Not available any more</p>	<p>ML3949 300W Xenon lamp, Cermax V300-Y20A for Luxtel MLX300 Light Source</p>	 <p>ML3949</p>
<p>ML1929 Timer (500h) Part# 1739</p> 	<p>ML1928 Timer (500h) CL1529 Part# 1167</p> 	<p>ML1956 Timer (500h) CL1555, CL 1555, TM-2, Part# 1038 (80-0424-000)</p>  <p>300W xenon arc lamp 500 hours timer, for MD-631 Y1064S,CLV-S40,CLV-240,CLV-260,CLV-260SL time meter</p>
<p>ML1174 Martin Set lamp 86-907-02-04 for KS Martin ML702HX/ML702 HX lamp 869070204</p>	<p>Martin ML702 HX lamp set including 90V 70W + 24V 120W halogen lamps Original Set Martin 86-907-02-04</p>	
<p>ML7491 ROYAL 328, 6V 0,2A, 6v 200ma 16x6mm</p>	<p>Microscope lamps without lens Royal 6V 0,2A BASE MF-6</p>	
<p>ML7436 ROYAL 611(LNS-MF46) with lens, ROYAL 312N, 6V 0,2A , 312 N, MF6, TL-1 3/4, 16x6mm</p>	<p>Microscope lamps with lens Royal 6V 0,2A BASE MF-6 , MICROSC LABOPHOT NIKON ELEPSE E-400 Lamp Royal 611, Lamp Royal 312N NIKON HFX-II</p>	
<p>ML5017 Royal LNS-BP310-706K Pentax 15BS. 3v 0,5A 3v 1,5w G2,54 Tubular lens end H 3V 1.5W G2.54 TL1, 100H LNSBP310 ROYAL light source LNS BP310</p>	<p>Royal LNS BP310-706K for Bronchoscope Pentax 15BS Royal LNS-BP310-305C 3v 1,5w G2,54, 3v 0,5A G2,54 Fibrosópio Pentax FT-13RBS Fiberscope Pentax model FT-13RBS Royal LNS-BP310-711C, Royal LNS BP310</p>	
<p>ML8908 HCR145/3B, 24v-40w G.16. HCR-145-3B, 24v 40w G16. HCR - 145B - 3B. SKYTRON 24V 40W SEALED BEAM HSB-3S 8893/1</p>	<p>HCR 150/3S. HCR150/3S HCR-150-3S, H145B3B/24V/40W , H145B3B, SKYTRON, MSHCR14.53B, MSHCR15.03S, SH-42/HALOGEN SKYLUX, SH-42/ SKYTRON, SH42/24V/40W SH42/24V/ SKYTRON, SURGICAL-EXAM- LIGHT/A101018, LIGHT/A101028, 8893/2</p>	
<p>ML8911 Skytron A101018 24v 40w G.16 195x42mm</p>	<p>A-101018/HALOGEN SKYTRON A1-010-18/24V/40W/SKYTRON A101018/HALOGEN SKYTRON 130060 NARVA, 130060 ORBITEC</p>	
<p>MR2240 LEUCI, 24V 150W E27 SCIALITIQUE 24V-150W E27 K4168892/3CH, 8892/3CH H-295133, H295133, 130016, 130016/24V/150W NARVA, S7 BBT</p>	<p>24-150w, 24-25v/150w E27 medidas: 150x90mm 24 150W, 24V-150W E27 P/CANDEEIRO SCIALYTIQUE medidas: 150x80mm</p>	



<p>ML8118 NARVA 650442/ 65044 NARVA 65042 E40D 220V 150W - 12V 100W 110x220mm Frosted Lamp</p>	<p>12/220V-100/150W E40D (110x220mm) TGL 8795 Lamps for operating Tables Narva 12/220V-100/150W TGL 8795, Lamps for operating Tables, NARVA 65044 /65042 (Frosted Bulb)</p>	
<p>ML5081 13276 G, E40d, 8893/3. BBT 24V-250/250W E40. 954044. BBT H 307592 250W 24V. BBT 130004. 130004 H 130004 (130 004) HT307592</p>	<p>24/24v - 250/250W 24-24V/250-250W 24 24V - 250 250W E40D 307592 24-24 250-250 E40 80 X 170 OPTOTAL 4 130004 (130 004) HT307592, 80X170mm Antique Lamp. Vintage Lamp</p>	
<p>ML4044 LÂMPADA 8v 30w P30S LÂMPADA 8v-30w P30S Para Microscópio Lâmpada para Tonômetro de sopra For Tonometer Goldmann For Goldman Tonometer Goldan 56x25mm</p>	<p>Lamp 8v 30w P30S Lamp 8v-30w P30S For Microscope For Tonometer Goldmann Tonometer</p> 	
<p>ML6211 7582, 12v 25w P30S 7582, 12v-25w P30S</p>	<p>1762AX ,7582, 7582 USA 7582X/12V/25W 7582X EUROPEAN 7582X USA 122416-A7582X-P30S-26X52- 12V-25W-AXIAL 122416/12V/25W/P15S</p>	
<p>ML1151 1762/2, 6v 30w PD30D 6v-30w PD30D ML8151 1762-H 7v 35w PD30D 7v-35w, 7v 5A PD30D 7v-35W, 7V-5A PD30D</p>	<p>ML1151 17622/2 I/6V/30W/P15D30/S8/AX/100H/17 62-2 ML8151 1762H I/7V/35W/P15D30/S8/AX/100H/17 62H</p>	
<p>ML8016 GE1630 6v 2,75A GE1520 6,5v 18w ML8017 GE1460X 6,5V 18w ML8018 GE1468X 6V 4,5A ML8019 GE 1634 20v 1A GE 1634 20v-20w</p>	<p>ML8016 GE1630 6v-18w PD30D GE1520 6,5v-2,75A PD30D ML8017 GE1460X 6,5v 2,75A PD30D ML8018 GE1468X 6v-4,5A , 6v-27w ML8019 GE 1634 20v-1A, 20v-20w P30D GE 1634 20v 1A, 20v 20w P30D</p>	
<p>ML8084 GE 1096, 6v-27w , 6v-4,5A 6v 4,5A, 6v 27W, 6V 4.5A, P15D30S 6v 4,5A, 6v 27W, P15D30 double contact 15mm with flange 30mm 8000048 USHIO 031293404015 EIKO 50,8x27mm</p>	<p>GE 1096, 6v-27w, 6v-4,5A P15D30S 6v 4,5A , 6v 27W, 6V 4.5A, P15D30S 031293415424, 031293404015 EIKO 37169 GENERAL-ELECTRIC 40401 EIKO, 40401 UPC, 41542, 41542 UPC, 621700U MICROLAMP I/6V/27W/P15D30S/S8/27MMLCL/500H 6v 4,5A, 6v 27W, P15D30 double contact 15mm with flange 30mm</p>	






<p>ML7416 GE 88441 Q1000-T7-BP 120V 1000W G38 (190 x 130 x 230 (mm)) General Electric Q1000-T7-BP GE 1MT20/BP-120V 1M/T7BP</p>	<p>ARIWAY BEACON Q1000 T7BP -120V MG. BIP. Base Down C-13 FILAMENT 1000W T7 Bulb 1M/T20BP 1000W 120V, 1M/T7BP G38 Mogul Bi-post Base Down GE 88525 GE 22429 Mg Bip. Tungsten Halogen Airfield Lamp C-13 FILAMENT. MS-25015-2 1000W</p>	
<p>Chromophare CZ 9086004 E805 9085/C</p>	<p>Chromophare CZ 9085004 X65/X66/E655 9085/C</p>	
<p>RJ2256 P-129249.001 33V 235W GY9,5, P-129249.001 33V-230W GY9,5 Vidro Fosco, Frosted Lamp (P-129249.001 33V-230W Surgical) 33v 230w Amsco Surgical P2533235, P128249001, P-128249.001P129249001, P129249001 Q/235/T4/2PPF/FR/33V, Q/235/T4/3, Q/235/T4/GY9.5/33V, Q/235W/T4/2PPF/33V. Amsco Surgical P 129352 157. Amsco P129-249-001 Amsco -P129-249-001, Frosted Lamp</p>	<p>RJ2255 Amsco P-129362-228 22V 220W, P129362-228, P129362-228CL/22V/220W Amsco P-129362-228 22V 220W GY9,5 Clear Lamp (Vidro Transparente) Amsco P 129362 228 22V-220W GY9,5 Frosted Lamp (Vidro Fosco) P129362-228CL/220W/22V, P129362-228CL P2522220/125341, P129362228, P129362228 USHIO Q220T4GY9.5-22V300H, QUANTUM-300, Amsco P129362-228 P129362-228 22V 220W GY9.5 Amsco, Ushio 8000355, Amsco P129362228</p>	
<p>ML4944 USHIO 1000503 USHIO 24v 150w GY9,5 FDS 5974 FDS/DEZ Germany Sylvania FDS Ushio 1000-503 USHIO medidas 60x12mm</p>	<p>1000503 BLV, 1000503 USHIO, 1000-503 USHIO, 6436269 OSRAM, 41007850-EOC PHILIPS, 41007850 PHILIPS, 3129301910 WIKO, 3129303210 WIKO, 04613554277 SYLVANIA, 031293032102 EIKO, 048777132760 UPC, 69549121928, 69549121928UPC STANDARDPRO, 83869 IWASAKI, 27554 NANTONG- DONAR, 4074 EFOS, 5722/1</p>	
<p>ML4496 Philips 7389, 230V-500W- GY9.5 PHILIPS FRANCE Philips 7389, 230V 500W GY9.5 Philips France 75x22mm</p>	<p>41044 PHILIPS, A1-244-230V-500W LIF A1-244-230V-GY9.5 LIF, A1-244-230V A1-244-230V ATLAS, A1-244-230V THORN A1-244/7389-230V-PHILIPS A1/244/230V/500W/GY9.5 LIF A1/244/230V/500W LIF, 5968 A1/244/230V GENERAL-ELECTRIC HALOGEN SSTV LAMPS 75x22mm</p>	
<p>ML4111 DYP 120v 600w G9,5B 57.2MM x 39.7MM</p>	<p>78-8454-3481-6 3M 120-600-32 CANIMPEX 600WATT-QUARTZ-G7-TP 01820 EIKO, 3481 3M, 32071. 54545/DYP SYLVANIA 54545, 1000248 USHIO 3129301820 WIKO</p>	
<p>ML4114 DYS 120v 600w GY9,5 600WATT-QUARTZ-G7-TP-2 625-32-09 LUSCOMBE 63MM x 36.5MM</p>	<p>121 3M, 176-092, 221 3M 244-1018 ALLIED 600-120-DYS-A1/264-17000-BDTH- CC-6-12.7X6.4-3200-64-37-75-GZ9.5- 24-CM-32955-9 G E, 5825-2, 5825/2 820 PERKO-MARINE, 01850 EIKO 01850 UPC, 01850 WIKO, 5269 3M</p>	
<p>ML4011 DYR 220v 650w GZ9,5 63MM x 36.5MM</p>	<p>78-8000-5349 3M 01840 EIKO, 01840 UPC 01840 WIKO, 5825-1, 5825/1 25058, 54839, 1000249 USHIO 762160 BULBRITE, 942503 048777115961 UPC</p>	













<p>NARVA 55122, 55125 NARVA HLWS5 6v 10w PY16-1,25</p> <p>NARVA 55133 NARVA HLWS5 6v 20w PY16-1,25</p> <p>RL55168 NARVA 55168 NARVA HLWS5-A 6v 30w PY16-1,25 HLW5-A 6v 30w PY16-1,25 Narva 000000-0402-943 NARVA</p>	<p>NARVA 55127 NARVA HLWS5 12v 10w PY16-1,25</p> <p>NARVA 55140 NARVA HLWS5 12v 25w PY16-1,25</p> <p>NARVA 55171 NARVA HLWS5 12v 35w PY16-1,25</p> <p>NARVA 55147, 55149 NARVA HLWS5 6v 25w PY16-1,25</p>	
<p>Narva 55204 NARVA HLWS5 12v 50w PY24,25</p> <p>Narva 55205 NARVA HLWS5-A 12v 50w PY24,25</p> <p>Narva 55305 NARVA HLWS5 12v 100w PY24,25</p> <p>Narva 55307 NARVA HLWS5-A 12v 100w PY24,25</p>	<p>Narva 55204 NARVA HLWS5 12v-50w PY24,25</p> <p>Narva 55205 NARVA HLWS5-A 12v-50w PY24,25</p> <p>Narva 55305 NARVA HLWS5 12v-100w PY24,25</p> <p>Narva 55307 NARVA HLWS5-A 12v-100w PY24,25</p>	
<p>Narva 57055 NARVA HLW55-P 24v 200w PY24-1.5 NARVA 57055 HLWS5-P 24V 200W PY24-1.5, 226.847</p> <p>Narva 55606 NARVA HLWS5 36v 400w PY24-1,5 HLWS5 36v-400w PY24-1,5</p>	<p>Narva 55503 NARVA HLWS5 24v-250w PY24-1,5</p> <p>Narva 55507 NARVA HLWS5-A 24v-250w PY24-1,5</p> <p>Narva 57055 NARVA HLW55-P 24v-200w PY24-1.5</p>	
<p>ML7711 NARVA 67273 NARVA LWT-P5T, 6V 15W, Z16 (67273) 6v-15W MB16, -P5T 6V 15W Z16 *p.n. 26.1679/06 *art.n. 068.134 *art.n. 681.34/3 *art.n. 0227-986, Remat, LMA 1</p>	<p>Erhitzungsmikroskop, Auflichteinrichtung fur Laboval pol d / Laboval und Ergaval, Metaval (Durchlichteinrichtung), Beleuchtungseinrichtung U m. XBO, Ergaval, Amplival, Amplival pol u, Amplival photometrie, Telaval, Citoval, Epignost, Elypovist, Laboval pol d, Technival, Vertival, Peraval interphako, Epival interphako, Epival, Fluoval, Fluoval 2, Nadelziehgeraet</p>	
<p>ML5019 Narva 67471 NARVA LWT-P3 6v 30W Z16 P.N. 26.4859/06 art. N. 068.142</p>	<p>ML5019 Narva 67471 NARVA LWT-P3 6v 30W Z16 P.N. 26.4859/06 art. N. 068.142</p>	
<p>ML7711 Hosobuchi 6v 15W MB16, 6V 15W Z16, Filament C-8Z</p>	<p>Hosobuchi 6v-15W MB16, NARVA 67273 NARVA 6V 15W Z16, Filament C-8Z</p>	
<p>ML5022 Narva 67491 NARVA 6v 30w</p>	<p>Narva 67491 6v-30w</p>	
<p>MC8825 Light bulb 230V 40W B22d clear light lamp incandescent ball. Lamp 220v 25w BA22D 45x74mm</p>	<p>MC8840 Light bulb 230V 40W B22d clear light lamp incandescent ball. Lamp 220v 40w BA22D 45x74mm</p>	




<p>ML5149 Narva 6V 10W BA15S (16x36mm) (19x37,5mm) max</p>	<p>Narva 17301 Microscope Leitz 500-096 Osram 5001 6V 10W BA15S (TUNGSRAM P6V10WR10W E20 2A9.4.III)</p>	
<p>ML7533 17161 NARVA 17161 R5W 6V 5W BA15s</p>	<p>(K5627). 6v-5w BA15S ELBX205 1.47 inch (37.3MM) x 0.748 inch (19MM)</p>	
<p>ML1164 NARVA 55241 HLL 22.8V 50W G6.35</p>	<p>ML1178 Narva 55913 HLWS4 Projector W711 HOYA, Heine Y-096.14.103 6V 20W G4</p>	
<p>ML1165 NARVA 55242 HLL 22.8V 110W G6.35</p>	<p>NARVA 55242 HLL 22.8V 110W G6.35. NARVA 55242 HLL 22.8V-110W G6.35</p>	
<p>ML4088 NARVA 55939 HLWS4-J 12V 100W Collimator 062-65003-02</p>	<p>NARVA 55939 HLWS4-J 12V/100W Replacement for 55309 Old NARVA-No. 55939 Narva. Applications: collimator, medical.</p>	
<p>ML4454A Philips 12345SL 12V 20W G4. CL770 Shimadzu Spectrometer Shimadzu Spectrometer CL770</p>	<p>Lamp Philips 12345-SL 12V 20W G4 Product D (max) 8.5 mm, L (min) 19.25 mm, L (max) 19.75 mm, L 19.5 mm, C 30 mm, P (max) 0.75 mm , P 0.7 mm, P (min) 0.65 mm</p>	
<p>ML4454 Lamp for CS projector CP-2147 NARVA 55917, HLL 12V 20W 2-Pin (G4) Clear. Narva HLL 12v 20w, 12v-20w. Lamp for CS projector CP2147, CP 2147, C.P. 2147, C.P.2147.</p>	<p>ML4454 Lamp for CS projector CP-2147. 12V 20W G4 lamp Narva HLL 55919. Q20T3G4FC6-12V. WA01149. Lâmpada para projector CS Modelo CP-2147. Halogen lamp with heat resistant pins</p>	
<p>ML5144 Microscópio Nikon Eclipse 400 Nikon Eclipse 50i 6V-30W, 6V 30W Nikon Eclipse 50 Microscope</p>	<p>78021, 84246 Medical CS-1 SLIT-LAMP Medical CS1 SLIT-LAMP Medical ENDOTHELIAL-CAMERA Medical Microscope Medical Olympus OL-8-C-410, OL8-C410, OL8C410 Swift, MA-326, MA326 Medical</p>	
<p>ML1148 Microscópio Nikon Eclipse 55i 6V-30W, 6V 30W</p>	<p>Nikon Eclipse Microscope NIKON ECLIPSE 55i 8460-8663-7721 Nikon Eclipse 55i Microscope</p>	
<p>ML2279 12V 100W G6.35 AX-F28 12V 100W GZ6.35 04613554256 SYLVANIA 4050300242958 OSRAM CL70056 CARLEY-LAMPS GZ6.35-11X39-12V- 100W/H62138/132134</p>	<p>0021995 Capsule - High Pressure 100W GY 6.35 12V Sleeve 5410288219950 12V 100W GZ6.35. 2300 3000 C, Lamp I5A8 FL 5 12v 100w 100W 12V, 100W 12V GZ6.35 HALOGEN MEDICAL & MICROSCOPE Microscope Olympus BX50 d:11,5mm, l: 37mm, a: 27mm</p>	
<p>ML6077 Lamp for Nikon NP-3 6V 20W G4 11143, 1120, Wild 362-658, 11143,1120 , P7388, Ushio 1000532, M30/ESB, S7388, M/30, Noritsu I061012-00, 409867, M30, 25678-4, JC6V20W-G4, 54019-0, 012407, M30, M/30, M-30, M34, M-34</p>	<p>Microscope MA780, 6518/2, 7388, PNO20, 409867, M34, 64250. Optical Equipment Microscope 1120, Nikon Microscope 79099 , Olympus Microscope 8C405, FU-2852, 3040- 0280, 31-31-81-01, Zeiss 38-79/9690, 3801- 20-1730, 3801-20-1740D, 3801-43-1350, 380120-1730, 380143-1350, 3801431350, 3879-9690, 500-245, 55912I, 655C, MB-300</p>	
<p>ML6446 Replacement lamp for Spectrophotometer 12v 35w GY6,35 44x11mm 130352 NARVA, 130352 ORBITEC, 1000823 BLV, 1000823 USHIO</p>	<p>ML6446S Replacement lamp for Berchtold Chromophare D200 12v 35w GY6,35 44x11mm Axial filament 44x11mm</p>	
<p>ML4211 Lamp ELHV BA22D Medida 48x24mm</p>		



<p>ML7486 Led Chip for Faromed 08-795, LED headlight.</p> 	<p>ML7486 Led Chip E09-LED-1W-06, Capacete de otorrino Faromed - EFFNER 08-795/LED. Faromed 08-795, LED headlight.</p>	
<p>OSRAM HBO 50W/ACL1 HBO 50W AC L1 OSRAM, HBO 50W/ACL1 39V 50W SPEC. Mercury Fluorescence Bulb OSRAM HBO50W/ACL1</p>	<p>AOPHOT/78564/L1 NIKON BH-LSF OLYMPUS CL70084 CARLEY-LAMPS DIAPLAN-MICROSCOPE/H286 LEITZ/LEICA DIAPLAN/H286 LEITZ/LEICA EPI-FLUORESCENCE-FLUOPHOT/78564/L1 NIKON</p>	
<p>RL126104 OSRAM HBO 100W/2 HBO 100W/2 (20V-100W (5A)) HBO 100 W/2 OSRAM 69217 - 100W 20V DC Short Arc Mercury Lamp - Two Electrodes. Unidade de Fluorescência marca NIKON modelo Y-Fi</p>	<p>RL126104A USHIO 5000274 SUH-102DH USHIO 5000274 - 5A 20V 100W Mercury Short Arc Lamp Equal to Osram HBO 100W/2</p>	
<p>ML4117 OSRAM HBO 103 W/2 100W 23V DC 2 ELECTRODES , 69182 OSRAM BRAND , LONG LIFE VERSION 300+ HOUR , 90MM LONG OSRAM HBO 103W/2 , ALB UPC: 4050300382128</p>	<p>HBO 103W/2 Mercury short arc lamp, 100 watts, 23 volts, DC current, 90 mm in length, average rated life 300 hours, operating position s90, base anode SFa7.5-2, base cathode SFa9-2, product code 69182 253 65D 09 2G, 25365D092G</p>	
<p>RL126103 Lamps N1217800 Olympus USH-1030L Microscópio Olympus MXV10</p>	<p>MVX10, MVX-2RE, MVX-CA2X, SZX-STAD1, SZH-STAD1, SZX-STAD2, SZ2-FO, SZH-SG, SZH-SC. Olympus USH-1030L Mercury Lamp</p>	
<p>ML7419 Lamp LT03098 6v 20w Medical bulb 20W6V Special: Halogen. Biochemical Analyzer. Corning 560 20W6V Special. Medical bulb 20W6V Special: Halogen bulb 6V 20W</p>	<p>LT03098, 6v-20w Medical bulb. Corning 560, 6V 20W Halogen lamps Voltage: 6 Wattage: 20 Life:100 hours</p>	
<p>ML7901 43630 FLAX Lamp Flax BLV 28,3</p>	<p>43630 FLAX FOR CUPPERS OCULUS</p>	
<p>ML2486 Topcon lamp for Retinal Camara TRC-NW7SF MARKII</p>	<p>Topcon flash lamp for Retinal Camara TRC-NW7SF MARKII code 4042366510</p>	
<p>ML4005 Amasco 129359-934</p>	<p>22V 200W (24v-200w)</p>	



<p>ML4451 2756/464C 12v 25w P20S RODENSTOCK</p>	<p>2756/464C 12v-25w P20S RODENSTOCK 50x18mm</p>	
<p>ML7145 Dulux L 18w/31-830. DULUX-L 18W/31/830 4PIN</p>	<p>34500-9 PHILIPS, 41088 GENERAL-ELECTRIC, 109722 HALCO, 180312 OSRAM, 300276 OSRAM, 345009 PHILIPS, 504510 BULBRITE</p>	 
<p>ML7896 Dentsply Qhl 75w Replacement Lamp 644657</p>	<p>Lâmpada para fotopolimerizador de Estomatologia DENSPLY-COMPATIBLE WITH QHL 75 tm- CURING LIGHT- SPECTRUM 800- REORDER#603.17.507 REORDER#644657 For Curing Light of Stomatology</p>	
<p>ML4577 Ohmeda Biliblanket 20v 140w GX5,3 Lamp. For Phototherapy Device BiliBlanket</p>	<p>6600-0680-200 Ohmeda biliblanket 20v 140w GX5,3 Lamp. LAMP OHMEDA 66OO-0262-200 (6600-0680-200)(H368)(H 368). Lâmpada OHMEDA 6600-0262-200 (6600-0680-200) (H368)(H 368). Para Aparelho de Fototerapia Biliblanket</p>	
<p>ML4578 Datex-Ohmeda 6600- 0680-200 halogen lamp 20v 140w with cable and plug</p>	<p>Replacement lamp 6600-0680-200 for: OHMEDA BiliBlanket Plus High Output Phototherapy System. Datex-Ohmeda 6600-0680-200 20V 140W MR16 w/Lead Wires and Connector</p>	
<p>ML1411 Dr. Fischer 6V 5W BA15D For device photocell Guillotine Polar brand. Polar Cutter Model No. 92</p>	<p>ML4990 Philips 13910, 24v 3w BA9S 22x9mm Bulb PHILIPS 13910CP T3W 13910 24V BA9S We also have in 24v 2w BA9S</p>	  <p style="text-align: center;">ML1411 ML4990</p>
<p>ML8106 Tungsten Halogen 4,8V 500ma 4,8V 2,4A E10 (Ø x L) Max 9X30mm</p>	<p>ML5580 Tungsten Halogen 6v 10w E10 (Ø x L) Max 9,3X31mm</p>	
<p>ML2102 Framiral mod. stimulo Lamps. 2,25v 0,25A E10. 222 Lâmpada com Lente (Ø x L) 9x20mm</p>	<p>Watts:0.56W Volts:2.25V Amps:0.25A Base:E10 MINIATURE SCREW BASE 10MM DIAMETER Length:23.6 Millimeters Life:5 Bulb Shape:TL3 Bulb Finish:CLEAR, osram 3651 E10 MINIATURE SCREW BASE 10MM Glass TL2-9/10, Lamp with lens</p>	
<p>ML7114 Philips 13876 24v 3w E10 Base E10SK 24V 0.125A (Ø x L) 12x20mm</p>	<p>13876 24V 3W E10 RECESSED-GLOBE PHILIPS 313/24V/3W/E10/RECESSED- GLOBE NAM-YUNG, 313 NAM- YUNG, OSRAM 3767 OSRAM Philips 13876 Philips</p>	 
<p>ML1116 Lâmpada Tube 220-260V 5-7 W 026 mA base = E14 220/260V 5-7 W 026 mA (Ø x L) 16 mm x 35 mm</p>		
<p>ML1115 Lamp 220/260V 5-7 W 0.02A 0.03A Base = E14 Clear 220v 7w, 220v 7w, 260v 5w, 260v 7w (Ø x L) 16 mm x 36 mm</p>		 <p style="text-align: center;">ML1116 ML1115</p>



<p>ML3003 99-N027 4,8V 0,7A E10 4,8v 700mA krypton, 4,8V-0,7A. 4,5/5v 0,5-0,8A, 712. Bulbs 1512 4,8v 0,5A E10. 4,5v 0,5A E10, 4,5v 0,8A E10 5v 0,5A E10, 5v 0,8A E10 Cegaza lamps (Ø x L) 11,5x24mm</p>	<p>ML3298 6v 500mA lamp for Olympus BX40 Microscope BX40 (Ø x L) 11x24mm. Lamp 6v 3w (6v 0,5A) (6v 500mA) (Ø x L) 11x24mm, G11x24mm (Ø x L) 11,x24mm</p>	 <p>ML3003 ML3298</p>
<p>ML1034 Hess Screen 6890D, Ecran Hess 6890D 2,5V 200mA Haag-Streit 2,5V 0,2A, 2,5V-0,2A (Ø x L) 11x24mm</p>		
<p>ML1033 Ecran Hess, Hess Screen 2,5v 300mA HAAG STREIT 2,5v 0,3A, 2,5v-0,3A (Ø x L) 11x24mm</p>		
<p>ML9122 Olympus 6V 5W E10, 6V-5W E10 LAMP FOR Olympus 6V 5W E10, Olympus 8G101 (OP8G101) 5W 6V E10 / MINI SCREW G4 1/2 Incandescent. Olympus 8G101 OP 8G101 (6V-5W) Microscope Replacement Bulbs/Lamps for Olympus BX-U-DO; Olympus STM4; Olympus BX-U-MDOB; Olympus BX-U-SDOOLYMPUS 8G101; Olympus DO; Olympus STM4 5 WATT; OP 8G101 (6V-5W); Microscope Olympus 8-G-101, Lamp ML9122 14x27mm</p>		 <p>ML9122 ML2306</p>
<p>ML2306 Lamps 12V 6W E10, Lamps 12V 0,5A E10, Lamps 12V 500mA E10 15mm x 29mm</p>		
<p>ML5578 Headlights Lamps 6v-0,45A E10, 6v 2,7w Ø 17mm, 17x30mm 125235 M E 56T 8 8/10 3,5/5w Clair. Espelho Frontal 6v 0,45A E10 6v 2,7W Ø 17mm, 17x30mm, Headlights Lamps 6v 500mA, Ø 55mm</p>		
<p>ML5579 Headlights Lamps 6v-0,90A E10, 6v-5,4w Ø 17mm , Storz 102200 Espelho Frontal 6v 0,90A E10 6v 5,4w Ø 17mm, 17x30mm</p>		
<p>ML4431 Headlights Lamps 10V E10 EFFNER Ø 15mm, 29x15mm Espelho Frontal 10V E10. Philips bulb 7181D 8V/10v-50mA-E10 EFFNER, Ø 15mm, 29x15mm</p>		 <p>10V (ML4431)</p>
<p>ML5919 Lamp for Karl Storz 6v 2,4W Headlights Lamps. K16439K5, 3724. K16439 K5, 47643, 12910 (Ø x L) 15,5x28mm.</p>	<p>ML6545 Lamp Hosobuchi OP-2128 6V 15W BA15D. Hosobuchi OP 2128 Hosobuchi OP2128 For Hikari Op2114</p>	 <p>ML5919 ML6545</p>
<p>ML3110 Osram 3626 Osram 3,7v 0,3A E10 Headlights Lamps (Ø x L) 11x23mm</p>		
<p>ML3109 Osram 3632 Osram 6,2V 0,3A E10 (Philips 714D) Headlights Lamps (Ø x L) 11x23mm</p>		
<p>ML3122 Osram 3715/2 Osram 6v 3w E10 Headlights Lamps Lâmpada para espelho Frontal 6V-3W E10 (Ø x L) 15x27mm</p>		
<p>ML7471 6V 5W E10 Frosted Headlights Lamps Guerra 7114/SM Lâmpada para espelho Frontal 6V 5W E10 Fosca</p>		 <p>ML7471</p>
<p>ML6666 6v 5W BA9S Lamp Microscópio Leitz HM-Lux3 (HM-LUX 3) HM-LUX, SM-LUX Leitz MPV, Leitz SM Attachable Illuminator, Leitz SM-M, Leitz SM-POL, MPV</p>	<p>Lamp for HM-LUX3 6V 5W BA9S lamp (ref. 500-073) for microscope IT Essilor Posicentron, AIT Posicentron Essilor, Elipsoidal Leitz, Essilor Posicentron AIT., HMLux 3 Leitz, Nidek, Leitz, Leitz Elipsoidal, 3040-0240, 3901-55, 3925-22 ZEISS</p>	 <p>ML6666 C6</p>
















<p>ML6541 Olympus 8-S101 6V 2A P27S M8000307 USHIO 6v 2A P27S 1/6V/12W/P27S/G6/CL 8S101 L-0911-VO INAMI MICROSCOPE-ILLUMINATOR- OME/8S101 OLYMPUS SM-8-LL-101, SM-8S101</p>	<p>ML6541 Olympus 8-S101 6V-2A P27S 6VOLT-2AMP, 6v-12w P27S. 8-S-201 OLYMPUS, Hosobuchi OP2163 8-S101/6V/2A OLYMPUS 8S101/6V/12W INAMI 8S101/G5-6V-2A-SCB-RECT-FLG 8S101 OLYMPUS, 22-2001, 22-2004 111-C CHICAGO-MINIATURE</p>	
<p>ML6544 Lamp 010209 6V 15w P27S 6V 2A P27S for Colposcope MedGyn AL-102, AL-102S, 6V-2A P27S 6V-15w P27S 3895/2, 3895-2 OP-2128X, 6V 15W-SC Colposcópico MedGyn AL-102S, MedGyn AL 102</p>	<p>ML6544 Lamp 010209 6V 12w P27S 6V 2A P27S Colposcope MedGyn AL-102. 6V-2A P27S, 6V-15w P27S OP-2128X 6V 15w P27S OP 2128 X, OP2128X, OP-2128 X MedGyn Products Bulb F/Colposcope 6v 2amp F/al-102 Ea - Model 10209</p>	
<p>ML6543 Lamp Colposcope 9v-2A P27S, 9v 18w P27S INAMI-9V2A/G5-9V-2A-RECT-FLG- SC-BAY, L0610-VI/9V/2A INAMI, L0610-VI, MICROSCOPE/L0610-VI INAMI</p>	<p>ML6543 Lamp Colposcope 9v 2A. P27S Lamp Colposcope 9v-18w P27S 09-018-04M MICROLITES, 9V/2A/P27S 9V/2A INAMI, 9V/COLPOSCOPE 9V2A P27S INAMI, 3895-1/9V/18W COLPOSCOPE-9V, COLPOSCOPE-9VOLT</p>	
<p>ML7006 Ushio USH-450GS 450 W att Short Arc Lamp for ASML Stepper 2500/40 with PEP Upgrade, USHIO super high pressure mercury lamps are designed for maximum output in the wavelength area of 436, 405 and 365 nm</p>	<p>ML7006 MPN: USH-450GS ASIN: B0041M7OGG, Ushio ush 450-gs USH-450GS/XENON-GAP-SHORT-450W LAMPTECH, USH-450GS LAMPTECH, USH450GS/XENON-GAP-SHORT-450W, USH450GS LAMPTECH. Ushio USH 450GS Short Arc lamp</p>	
<p>ML7008 Ushio USH450BP, USH 450BP. Bulb GS 450W Super High Pressure UV Lamp USH- 450BP Basis Print UV STTER 57 Basis Print UVSTTER57 Graphic Arts</p>	<p>ML7008 Ushio USH-450BP, Ushio's super high- pressure UV lamps have been developed utilizing the three wavelengths of UV rays (436nm, 405nm and 365nm). The lamps are highly illuminant light sources that offer long life and stable radiation intensity</p>	
<p>ML7012 Ushio USH-1000DW. These lamps feature at a stable radiant power, and since their size closely approximates a point light source, the light they emit can be easily concentrated</p>	<p>Ushio USH-1000DW USHIO super-high mercury lamps are designed for use in the ultraviolet wavelength region, particularly at 436, 405 and 365nm.</p>	
<p>ML7010 Ushio UXL-S150M-0 (S150M-0 5000375) 17,5V 150w (8,5A) Olympus, OL-8-B-195 Microscope, OL8B-195 Microscope, OL8B195 Microscope</p>	<p>Ushio UXL-S150M-0 Long Life 150W 17.5V Projection UXL Series Long Life Short Arc Xenon Lamp , UXLS150M0</p>	
<p>ML7044 Lightness 5800000LX Xenon Bulb for Light Source,</p>	<p>350W Xenon Bulb for our Light Source, Lightness 5800000LX</p>	
<p>CL3634 Halogen Quartzline lamp for overhead projector. 36v 340w GX5.3. GE 36v 340w lamp. ERV 36v 340w</p>	<p>Color Temperature: 3300K, Average Rated Life: 75 Hours, Filament: CC-8, CC8 ERV 41874 GE 36v 340w Overhead Quartzline® Multi-Mirror® MR16 - Overhead Projection</p>	















<p>ML8899 EXZ (Q50MR16C/NFL25° EXZ) 12V 50W (code 20835) EXZ 25°/27°. GE Code 20835 Q50MR16C/NFL25°/EXZ AIR SHIELDS VIEMERS IIC 5-90/V210726. Light Phototherapy. Dräger MU03692 Draeger Medical Q50MR16C/NFL25° EXZ 12v 50W Only Medical Lamps pediatric phototherapy</p>	<p>GE #20835 Q50MR16C/NFL25 12v 50W GE EXZ 50 Watt Precise Constant Colour Halogen MR16 Without Cover Glass -GX5.3 Bi-Pin Base • 25° Narrow Flood Beam Spread • 12 Volt • M50 EXZ precise 330709 •Color Temp 3050K •Rated Avg. Life 6000Hrs, MOL 1.875in (47.6 mm) x Diam 2in • Ansi Code EXT. EXZ BLV 189852 12v 50w 24° Light Phototherapy, Light therapy</p>	
<p>ML7789 EXZ, General Electric 330709 M50-EXZ 12V 50W 24°. 048777127728 UPC,048777127933 UPC 807154171286, 807154174904 HALCO, 5021731307096 5994100094371 TUNGSRAM 5994100094432 TUNGSRAM</p>	<p>509437 TUNGSRAM, 509443 ,521015 GENERAL-ELECTRIC, 641250 BULBRITE 842091/12V/50W , 1000424 USHIO 1000427 USHIO, 04613554205 SYLVANIA 031293150462 EIKO, 043168274647 0043168398749, 0043168937603-045923019623, 046135585951 SATCO EXZ-12V-50W-GU5.3-24-DEG</p>	
<p>ML2211 EXN 50W 12V, GU5.3, MR16 36D 1CT (Diameter D: 50 (max) mm, Overall Length C: 46 (max) mm, Cod. 924047417122. USHIO 1000398 EXN USHIO JR12V-50W/FL36 EXN PQ50MR16/FL EXN, VISUAL EFFECTS GOBO SCANNER EXN AIR SHIELDS 78 SERIES INFANT WARM, EXAM: EXN, 44870 WFL</p>	<p>Lamp Halog. 12V 50W WFL MR16, GU5.3. 50MR16/FL35/EXN/C 12V (EXN) MR16 Halogen Light Bulb. LAMP, HALOGEN, 50MM, 12V-50W, 38DEG. WELCH ALLYN 04450 WELCH ALLYN 04450-U EXN, 44500 EXN WELCH ALLYN 44510 EXN, LS-150 EXN WELCH ALLYN LS150 EXN VISUAL EFFECTS V-8210-1 EXN, V-8210-4 EXN , AIR SHIELDS 7830 MAIN/EXAM BULB: EXN</p>	
<p>ML2591 G.E. 8V 50W GX5,3 EFM MEDICAL DENTAL DJ PROJECTOR BULB LAMP Lamp 64607, 8V 50W GX5,3 Lamp 5132/1, 8V 50W GX5,3</p>	<p>GE GENERAL ELECTRIC G.E 41251. 8V 50W EFM MR16 HALOGEN MEDICAL DENTAL DJ PROJECTOR BULB LAMP General Electric 41251 EFM. Dim: 50MM x 42MM</p>	
<p>ML5006 Dental Lamp 12v 75w para dental (5332/1). 403/OPTILUX, Osram 64617 Osram, para Fotopolimerizador. Flat Pins</p>		
<p>ML5029 Dental Lamp 12v 75w para dental. 403/OPTILUX, Osram 64617S 404-HIGH-Output Demetron, para Fotopolimerizador. Flat Pins</p>		
<p>ML4284 Dental Lamp 12v 80w Bulb For Demetron 21237 12V 80W Medical Dental Curing Optibulb Lamp. Fotopolimerizador Demetron 21237. Demetron Kerr Optilux 21237 Dental Optibulb Lamp. Flat Pins</p>		
<p>ML5007 Dental Lamp 12v 100w para dental (5332/2). JCR/M 12V 100W SD G5.3-4.8 - Dental Curing Lamp. Fotopolimerizador. Flat Pins</p>		
<p>ML1599 Dental Lamp 13V 100W GZ4 13v 100w G4 MR11 halogen Dental Lamp 14847 13V 100W 13v/100w GZ4 MR11 halogen lamp for dental microscopes 14847, Philips 14647 13v 100w</p>	<p>ML7466 Dental Lamp 6V 35W GZ4 Osram 64600, Dental 64600 WHITE 6V 35W G4 3M Litema Vivident Lamp Bulb (069). Lamp 5331/3 131465-H68635-GU4-35X40-6V-35W Dental lamp table L145</p>	
<p>ML4414 EJY Projector Bulb 19V 80W GX5.3 Base 25hr Life Fiber Optics. Equipment Use: Microscope ANSI Lamps, Proj/Fiber Opt. GE-EJY-32886-19V 80W, 10272 EIKO BRAND, CC-6 FILAMENT. UPC: 031293102720, P.N. 10272 Manufacturer UPC 031293102720</p>	<p>EJY/19V/80W/GX5.3/MR16D ANSI, 19V-80W GX5,3 have been found in equipment: American Optical Corporation, 1180 Medical 1182 Microscope Fairchild Camera and Instrument Corporation 07 SEVENTY 8mm Projector Movie 21 SEVENTY 8mm Projector Movie 27-P SEVENTY 8mm Projector.</p>	
<p>ML2337 Lamp. 14,5V 90W GX5.3 EPV for Microscope WILD 14.5-090-23R/EPV MICROLITES 580 SK, 800 SK, 900-G SK 1000-347 USHIO, JCR 14,5V 90W 2000/HIGHLIGHT OLYMPUS EPV/14,5v/90w/GX5,3/MR16D</p>	<p>Para Microscópio WILD Intralux Model MR 2-1. Lamp. 14,5V-90W GX5.3 EPV Microfilm Lamp - MR16 - 90 Watt Light Bulb - 14.5 Volt - GX5.3 Base - 3150K - EPV/14.5V/90W/GX5.3/MR16D ANSI. EPV 14,5v 90w GX5,3 MR16D, 13186 EPX/EPV Microscópio INTRALUX MR2-1 13186 EXP/EPV 14,5V 90W GX5.3 409805</p>	









<p>ML7672 H63144 EPN A1/230 12V 35W GX5,3 H63144 12V 35W GX5,3 EPN 12V-35W GX5.3 (mm): 44.5</p>	<p>12-035-23NR, 12-035-23NR, 99-LL-447 02660 EIKO, 02660 UPC, 30165 12V 35W, 02660 WIKO, Welch Allyn 04200, 9248-036- 17100, 9248-037-16600, 11377, 25118, 41232 EIKO 41232 UPC, 41430 48435, 48600, 48610, WA- 04200, 48625, 48635</p>	
<p>ML6531 LAMP 19V 80W GX5,3 DDM Projector Bulb 19V 80W DDM/19V/80W/GX5.3/MR16D, LAMP 5132/7 19V 80W, 04877111482 UPC, 69549120280, 69549120280UPC, STANDARD PRO 6240011487644</p>		
<p>ML7492 842039 Products for Medicine 24v 200w GX5,3 EJL, Projector Lamp 13164, 27980-7, 29150, 29305, 31508-5, 52240, 54730-EJL-24V-200W-GX-5.3, 54730- EJL-24V-200W-GX-5.3, 54730/EJL, 54730, 54730, Lamp 5132/12 24v 200w GX5,3 EJL</p>		
<p>ML7702 82v 360w GX5,3 MR16D ENX/82V/360W/GY5.3/MR16D 54984/ENX SYLVANIA 5132/14, 5132-14, 41705 G.E. ENX 82v-360w GX5,3 MR16D Flat Pins</p>	<p>ENX 82v-360w GX5,3 MR16D ENX/82V/360W/GY5.3/MR16D ANSI 9070 3M, 9080-PROJECTOR 3M 9080 3M, CL70009 CARLEY-LAMPS CL70075 CARLEY-LAMPS ENX 82v 360w GX5,3 MR16D Flat Pins</p>	
<p>ML4060 ESD WALLACH 150W 120V GY5.3 MR16. Lamp. 906117 Colposcópico ACG141C3 Res. #8114738. NIKON SMZ-1. 131242 120v-150w GY5.3 Flat Pins</p>	<p>120V 150W ESD, 120v-150w JCR GY5.3 Ushio BC2067 1000940 - JCR120V-150W/B MR16 Halogen Light Bulb PC2013060219 Wallach Zoom Scopl-Nikon (Atenção ver suporte, Pinos Rectangulares) Flat Pins</p>	
<p>ML8991 KAVO 120V 250W ETJ, 555-5001 120V 250W GY5,3, KAVO TYPE POLAR 2, KAVO 60513 ETJ 120V 250W, Polar II; Vicon DLS; Vicon LS; Vicon LS II (Older), DENTSPLY - 70353, DENTSPLY 70353, 131244, 5555001 KAVO-AMERICA. 2-pin; flat .7mm x 2mm</p>		
<p>ML9981 13095 ENH 120v 250w GY5,3 Japan 316216, 93506, 122007 1000-333 USHIO, 1439 BLV, 02570 EIKO GY5.3 2-pin; flat .7mm x 2mm; 5.3mm centered</p>		
<p>ML6116 Fujinon LMP-SD 15V 150W Fujinon EPX-201/ EPX-201H 15v 150w GZ6,35 LMP-SD FV7A PI-1(UCC/EAN-128) Lamp. Vidioprocessador XP / E200 Ref. LMP - SD. Fujinon 0802 FV 7A LN, NECK/15V/150W/GZ6.35/MR16DSN</p>		
<p>ML7633 15V 150W, Olympus CLE165, CLE 165, CLE-165, Evis Exera II. LH150P11 PENTAX, LH5000-ALGOTECH LITEMA, LT-64634, MD151 OLYMPUS, Olympus CLH-2, Pentax LH-150 PC, Pentax LH-150PC, Zeiss, 6423FO XHP 6423 FO Fiber Optic Lamp</p>		
<p>ML5412 ZEISS 12V 75W GZ.6,35 EFN, 3800-79-9030, 3800-79-9930, 3879-9030, 4170-80-9015 ZEISS, Lamp 5132/2</p>		
<p>ML5514 Hanaulux 18565, 12V 100W GZ6,35 EFP Hanaulux 018565 Lamp 5132/3, 3879-9040 ZEISS, 098-145 WILD LEITZ</p>		
<p>ML6117 Fujinon LMP-EXP for Vidioprocessador EPX-2200 up to EPX-4400EXP. The LMP-EXP is used for the newer light sources like EPX-2200 up to EPX-4400</p>		<p>No Picture</p>
<p>ML4134 424549 12v 20w ESX 10° Lamp. 12v 20w 10° ESX 12v- 20w OR 10° ESX 12v/20w OR 10° ESX MR16 Dicroic 20w 12v ESX 10° Beam Angle</p>	<p>ML8208 132491 NARVA 12V 12W GU4 FTA 7° Beam Angle MR11 Dicroic 12w 12v Sistema de VERTIGENS. Spinning system lamp Dim: 35.3MMx35MM.</p>	
<p>ML4133 Lamp. 12v 20w 12° ESX-P Lamp. 12v-20w 12° ESX-P 12v-20w OR 12° ESX-P GX5,3 12v/20w OR 12° ESX-P GX5,3</p>	<p>Lamp. 12v 20w OR 12° ESX-P 12v 20w GX5.3/GU5.12° ESX-P 12v 20w GX5.3 12° ESX-P MR16 Dicroic 20w 12v ESX-P / M268 12° Beam Angle Prolite MR16 Dicroic 20w 12v ESX-P/M268 12° Beam Angle</p>	



<p>ML4068 BAILEY HM301202030 FTD MR11 GU4 OPEN 12V 20W 30º Bailey, Diameter 35mm, Total leng 35 mm. EAN Code 8714681079454, Art. Nº 245691, Prod.N. 9000308 12V 20W GU4 FTD</p>	<p>12V-20W FTD Electra1400C/1600C MLA 900C 70307 12V-20W Lâmpada KIT LIGHT BLUE 12V-20W have been found in equipment: Micro Strip MICS-1 Studio / Stage, MICS-2 Studio / Stage MICS-3 Studio / Stage, MICS1 Studio / Stage MICS2 Studio / Stage, MICS3 Studio / Stage</p>	
<p>ML7788 Spare Bulbs for ADMECO- LUX 22,8V 50W GX5.3 art.40.0092.80 (40009280) 13938 22,8V 50W GX5,3 13938 22,8V-50W GX5,3 [100022] 170.52 RON. Lamp 00847160</p>	<p>ADMECO-LUX 22,8V 50W art.40.0092.80 (40009280) 13938 22,8V 50W GX5,3 13938 22,8V-50W GX5,3 [100022] 170.52 RON Mit Reflektor 00847160 22,8V 50W axial GX5.3. Philips 13938XHP 22.8 V 50 W GX5.3</p>	 <p>Medizinlampe 22,8V 50W GX 5,3 00847160 Made in Germany </p>
<p>ML6509 EKZ 10,8V 30W GX5,3 10,8V 30W EKZ GX5,3 EKZ/10.8V/30W/GX5.3/MR16D 63142, Leica CLS 50X</p>	<p>130301 NARVA, 130301, 239459 842028/10.8V/30W 842028 1000315 94000574 MILES/AMES, 3129302370, 031293023704 EIKO 0043168993388-EAN, 0043168993388</p>	
<p>ML6510 EPZ 13,8V 50W GX5,3 13,8V -50W EPZ GX5,3 . EPZ/13.8V/50W/GX5.3/MR16D</p>	<p>11-READER/PRINTER/13.8V/50W MICRON 13.8V-50W-MR16-GX5.3 USHIO</p>	
<p>ML4064 Philips 13528 6v 15w GZ4 (31504-4) Lamp Philips 13528 6v-15w GZ4 Lamp Bulb Replacement 13528 MR-11 MR11 MR 11 MR/11 15 watt (15w) 6 volt (6v) GZ4 (Miniature 2-Pin) C-6 BD BC1044 Order Codes: 13528, 409782 FocusLine 409782- 31504-4 315044 31504 4 31504/4, 5331/5, 5331-5</p>		
<p>ML7673 Philips 13529 6v 9w GZ4. 13529 Philips 6v-9w GZ4 para Microscópio Leica ATC2000 (Microscope Leica ATC 2000) 409911, 315077 PHILIPS</p>		
<p>ML2456 USHIO EKN 17,7V 120W GE EKN 17,5V 120W GX5,3. EKN17.7V/17.7V/120W/GX5.3 /MR16D ANSI.</p>	<p>USHIO EKN 17.7V 120W (1000311) Lamp, Ushio 1000311 - Voltage 17.7 Volt ANSI Code EKN Base Type GX5.3 Bulb Type MR16 Projector Bulb Color Temperature Kelvin 3200K Over all Length M.O.L 1.75 in. General Electric EKN 17,5V 120W</p>	
<p>ML9908 Philips 13298 10V 52W GZ4 Glass: MR11, 35 DICH, MR-11. Para aparelho fotopolimerizador Philips 13298 52w 10v (10v 52w) Light Bulb, replacement dental Lamp. 3M: XL 1500, XL1500, XL-1500 Spectrum 200R, Jetlite 4000 Plus Patterson: TCL-490, TCL490</p>	<p>Philips 13298 10V -52W GZ4 Philips fOCUSLINE dental lamp 13865, 14552, 13165, 13298. Cross Reference Order code Volts Watts Base Philips 13298 LT 05015 10 52 GZ4 Philips 14552 LT 05026 12 75 GZ4 Philips 13865 LT 05027 12 75 G5.3- Philips 13165 LT 05036 14 35 GZ4</p>	
<p>ML5508 DNF 21v 150w DNF/21V/150W/GX7.9/MR18D ANSI, Fonte de Luz APPIED FIBEROPTICS 21V 150W</p>	<p>DNF 21v-150w GY7,9 4411, 54411, 61351 STAR-DENTAL, 63154 86186. 86303 PENTAX, 93631 /21V 150W GX7.9 , 1500-LIGHT-SOURCE-BULB-150W- 21V APPLIED-FIBEROPTICS</p>	
<p>ML5509 DNE 120v 150w DNE/120V/150W/GY7.9/M R18D , NIKON SMZ-1</p>	<p>DNE 120v-150w GY7,9 DNE WIKO, DNE120V DNE120V/54409 SYLVANIA DNE120V/04613554409 SYLVANIAGTL-OS- 54409, WALLACH ZOOM SCOPE</p>	
<p>ML4084 Leica Original HLX 64606. Art.Nr. 10446809, 21V 150W EKE Lamp for Microscope Leica M500/M300 Leica 10446763, HLX64608, HLX64606, 21V-150W EKE. 10446809. For Leica 10446763, 10446809, 313-050, 313062, 902-058. Osram 64606 Osram HLX64606, Leica M300/ Leica M500, 21V 150W GX5,3. Leica Original Lamp Osram 64606. Art.10446809 lamp. Fibber Optic Lamp</p>		 <p>2.0 in. (50.8mm) 1.75 in. (44mm) OSRAM HLX 64606 for Leica M500 / M300 Art.Nr. 10446809 21V / 150W Made in Germany</p>
<p>ML7429 Lamp 21v 150w GX5,3 MR16, 04050 54753, EJA/EJV Bulb Ushio 1000297, 32882, 54753, 44142-8, Lamp 5132/11/18,. Fibber Optic Lamp</p>		



INCUBADORAS, MATERIAL DE SUBSTITUIÇÃO REPLACEMENT PARTS FOR INCUBATORS OTHER REPLACEMENT PARTS		
CÓDIGO DO ARTIGO	DESCRIÇÃO	FIGURA
ML5349 Transformateur électronique 230v/12v TLP 60.230/12 STR Verre et Quartz 15x3x2,8cm	Electronic transformer TLP 60.230/12 STR Electronic transformer 230v/12v Verre et Quartz Transformateur électronique TLP 60.230/12 STR 15x3x2,8cm	
ML7522 Motor 240v, para Incubadora Isolette Mod. C86 e C86H	Motor replacement part number 68 905 73 Isolette Infant Incubators Models C86, C86H, C100 and C200	
ML7523 Motor 240v, para Incubadora Isolette Mod. C86 e C86H	Motor replacement part number 68 905 73 Isolette Infant Incubators Models C86, C86H, C100 and C200	
ML7524 Sensor incubadora Isolette Mod. C86 e C86H	Isolette infant warmer incubator Sensor Mod. C86, C86H C100 and C200	
ML7525 Humidificador Incubadora Isolette Mod. C86 e C86H, C100 E C200	Isolette Infant Incubators Models C86, C86H, C100 and C200, Humidity Fill Spout replacement Kit part nº 68 902 70 Air-Shields a healthdyne Company	
ML1712 Mercury redistilled, 99.9+% MU5005 Glass bottle for solids 500g. Mercúrio Puro a 99,9% frasco com 500g. Mercurio Puro		
ML7527 Motor fan Incubator Isolette Mod. C86 e C86H, C100 E C200 (86x47mm) furo de 6mm. Diameter of hole 6mm	ML7526 Motor fan Incubator Isolette Mod. C86, C86H, C100 and C200 (67x25mm) Diameter of hole 6mm	
ML7104 Ohio humidificador Humidificador Ohio Diameter 170mmx125mm	Ohio fan humidifier. Diameter 170mmx125mm	
ML2603 Fujitsu Computer Board EZ-0963HUE-C. Fujitsu 9703523014 PCB Control EZ-0000HUE-C. 0963HUE-C W91109-CCP-3400 FUJITSU GENERAL AIR CONDITIONER (OUTDOOR UNIT Model AOG25RZAL)		