






LÂMPADAS ESPECIAIS PARA EQUIPAMENTO HOSPITALAR		
CÓDIGO DO ARTIGO	DESCRIÇÃO	FIGURA
<p>ML7488 F4T5/BL-B 7G Lamp. Fluorescente UV, Luz Negra, Black Ligth Blue 6 INCH 4W, F4T5/BLB, HR4T5BLB</p>	<p>ML7489 Philips TUV 4 Watt G4T5 ultravioleta Tube / Bulb, TUVG4T5, Germicidal lamps G4 T5 (TUV4W) 15x134mm</p>	 <p>ML7489</p>
<p>Fuji Medical Part Number 123N0007D lamps</p>	<p>Ref: 123N0007D (length 480mm diameter 15mm) for device digitalisator Fuji photo film Model: CR-IR 346RU S/N: 36526210B replacement erase lamps for Fuji computed radiography readers</p>	
<p>ML4717 Lamp TUV PL-L 95W/4P HO 1CT/25 TUV PL-L 95W/4P HO ref. 88829740</p>	<p>Lamp TUV PL-L 95W/4P HO 1CT/25, Lamp 88829740 Emitem radiação UV de onda curta com um pico de 353.7 nm(UV-C) para acção germicida; -O vidro da lâmpada filtra a linha de formação de ozônio de 185 nm</p>	
<p>Waldmann F85/100W-PUVA F85/100W PUVA, F85/100W- UV21, F85/100W-TL01. PUVA 100, 1000, 3001, 3002, 3003, 6001, 6002, 701, UV 1000KL, UV 1000K, UV 3002K, UV 3003K, UV 7001K, UV 8001K, UV 7002, UV 5040 KL, UV 5001 KL, UV 1000 KL</p>	<p>Indications: -Psoriasis, atopic dermatitis -Vitiligo, Lichen ruber planus -Mycosis Fungoides (MF) -Scleroderma, acne vulgaris -Uraemic pruritus -Polymorphous light dermatosis</p>	
<p>TL100W/01, TL 100W 01 Dermatological phototherapy lamp Dim: 1760mm x 38mm medical lamp nm:300 nm:315 phototherapy lamp therapy lamp ultra-voilet- B/01. Emission: between 305 to 315 nm with a peak of 311 nm</p>	<p>Lâmpada dermatologia fototerapia dim: 1760MMx38MM Medical lamp, phototherapy lamp therapy lampultra-voilet-B/01. (F/72INCH/T12/R17D/100W/300N) TL100W 01. 007873/122V/100W/R17D, 007873 (639745),8711500614254 PHILIPS F/69INCH/T12/R17D/100W/310NM/12 2V F/70INCH/T12/R17D/100W/UVB-01 F/72INCH/T12/R17D/100W/01 F/72INCH/T12/R17D/100W/01-UVB</p>	
<p>Electronic Starters N61347-2-1, 4-125w wadmann OKO6 - the most popular starter to work in temperatures from -25 for C to 100 of C, the fluorescent 4-125W, OKO2 - for serial connections, for temperatures of -25 with C to 100 of C, for 4-22W fluorescent lamps and compact fluorescent lamps and wheeled the dim power range. Electronic.</p>	<p>Arrancadores electrónicos FRIGOSTART - to work in temperaturachod -40 for C to 100 of C, for 4-125W fluorescent lamps (used in freezers) TURBO - for applications where an instant start fluorescent lamps (-5 for C to 100 of C) for fluorescent 4-65W SOLARSTART - to tanning salons (4-180W). Source: company materials S12 Sunstarter, 115w - 140w Tubes Cleo Body Tone Starte 100-180W Cosmostarter/E (15w-225w Sunbeds)</p>	



<p>UVA7001K, UV7001K Phototherapy for psoriasis Ultraviolet light therapy for atopic dermatitis. 97612638 UV lamps for device UV7001K, UVA7001K, UVA 7001K, UV 7001K, UV 7001 K , UV7001 K Lamp. dim:1760MMx38MM Phototherapy lamp, therapy lamp ultra-voilet-B/01. TL100W/01, F85/100W-UV01 (F85/100W-TL01) UVA GOLD 100W (TL85/100W/09) (WALDMANN F85/100W-PUVA), TL100W-10R. When replacing the fluorescent lamp it is recommended to also replace the starter Waldmann electronic ballast</p>	<p>Combined therapy: 27 pcs of F85/100W-PUVA +13 pcs of F85/100W-TL01 or 27 pcs of F85/100W-PUVA +13 pcs of F85/100W-UV21 or 27 pcs of F85/100W-PUVA +13 pcs of TL100W-01 or 20 pcs of F85/100W-PUVA +20 pcs of F85/100W-TL01 or 20 pcs of F85/100W-PUVA +20 pcs of F85/100W-UV6</p> <p>Waldmann Starter (for 100W tubes) Waldmann ballast (for 100W tubes)</p>	
<p>UVA7001K, UV7001K Phototherapy for psoriasis Ultraviolet light therapy for atopic dermatitis. 97612638 UV lamps for device UV7001K, UVA7001K, UVA 7001K , UV 7001K , UV 7001 K, UV7001 K Lamp. dim:1760MMx38MM Phototherapy lamp, therapy lamp ultra-voilet-B/01. TL100W/01, F85/100W-UV01. UVA GOLD 100W (TL85/100W/09) F85/100W-01, (WALDMANN F85/100W-PUVA), TL100W-10R.) When replacing the fluorescent lamp it is recommended to also replace the starter Waldmann electronic ballast</p>	<p>Monotherapy: 40 pcs of F85/100W-PUVA or 40 pcs of F85/100W-TL01 or 40 pcs of F85/100W-UV6 or 40 pcs of F85/100W-TL10R or 40 pcs of F85/100W-UV21 or 40 pcs of F85/100W-UV01 or 40 pcs of TL 100W-01</p> <p>Waldmann Starter (for 100W tubes) Waldmann ballast (for 100W tubes)</p>	
<p>UV 7002 UV7002, UVA 7002, UVA7002 - Therapy System. Radiation (UVA): Lamp waldmann: F79/120W-PUVA 320-410nm. Radiation (UVB): Lamp waldmann: F79/120W-TL01 310-315nm. When replacing the fluorescent lamp it is recommended to also replace the starter Waldmann electronic ballast</p>	<p>Combined therapy: 21 pcs F79/120W- PUVA + 21 pcs F79/120W -TL01</p> <p>Monotherapy: 42 pcs F79/120W - PUVA or 42 pcs F79/120W - TL01</p> <p>Waldmann Starter (for 120W tubes) Waldmann ballast (for 120W tubes)</p>	
<p>UV 5040 KL (UV 5040KL) UV5040KL UV-B (UV6) 285-350nm UV-B (UV21) 285-350nm</p> <p>UV 1000 KL UV1000KL, UV 1000KL, UV1000 KL</p> <p>Monotherapy: 26 F85/100W-PUVA 26 F85/100W-UV6 26 F85/100W-UV01(TL01)</p> <p>Combined therapy: 20 F85/100W-PUVA + 6 F85/100W 16 F85/100W-PUVA + 10 F85/100W- UV01</p>	<p>UV 5040 KL Combined therapy: 26 F85/100W-PUVA + 16 F85/100w-UV21 or 24F85/100-PUVA + 16 F75/100W-TL01</p> <p>Monotherapy: 40 F85/100W-PUVA or 40 F85/100W-UV6 or 40 F85/100W-TL01</p> <p>Waldmann Starter (for 100W tubes)</p>	



<p>PHILIPS TLK 40W/03 40W T12 24 LONG 47V VHO SUPER ACTINIC 308072 LIGHT EMISSIONS TECHNOLOGY TLK 40W/03 Philips BL TL-K 40W/10-R With Reflector</p>	<p>Super Actinic PHILIPS TLK 40W/03 for treatment of hyperbilirubinie Philips Actinic BL TL-K 40W/10-R Reflector Philips Actinic BL 40W with reflector</p>	
<p>F20T12/BB 31745 Luz Azul TL20w/52. TL20w/03T (ref.9511780661). Fototerapia Blue Vickers ref. 80030 BL4. AMEDA ZUG mod. FTL. (TL20W/03T) TL 20W/52 Philips LÂMPADA TL 20W/52 RS G13 Phototherapy of Hyperbilirubine</p>	<p>Philips F20T12/BB 31745 Luz Azul TL20w/52. TL20w/03T (ref.9511780661). Fototerapia Blue Vickers ref. 80030 BL4. AMEDA ZUG mod. FTL. TL20W/01RS. UV-B Medical. (TL20W/03T) Phototherapy of Psoriasis. TL 20W/52 SLV. F20T12/BB 20 Watt T12 Fluorescent Bulb Bilirubin Blue. Philips TL20/52 AIR-Shilds mod. PT.533-1E Phototherapy of Bilirubin</p>	
<p>TL20W/01 RS. UV-B Medical. This is the most efficacious waveband for the treatment of psoriasis, and is therefore much more focused. What's more, because the overall dosage of this narrowband radiation can be closely controlled, these lamps are suitable for home therapy.</p>	<p>Lamps have a very narrow waveband emission: between 305 to 315 nm with a peak of 311 nm. The UVB Narrowband (01) low pressure fluorescent Phototherapy lamps are proven to be the most effective lamps for treatment of skin diseases. Can be used for treatment of: Psoriasis and Parapsoriasis Vitiligo Atopic Dermatititis Mycosis fungoides Other Skin diseases Prescription Required!!!</p>	
<p>F20T12 C-50. F20T12-C50, TL20W/33RS (F20T12/D). AMEDA ZUG mod. FTL. F20T12/D / ECO "Daylight". Average Life Hours 9000, , Length 24 In, Pre Heat, 1 1/2 In (38mm) Diameter</p>	<p>For treatment of hyperbilirubinie F20T12/C50; F20T12/C50/ECO, Fluorescent Lamp, T12 Lamp Shape, 20 Watts, F20T12/C50/ECO, Base Medium BiPin (G13), Color Temperature 5000, Lamp Tone Enhanced Daylight, CRI 90, Rated</p>	
<p>ML9229 UV Lamp UV083245 UV LAMP K-14 5000 HOURS Brand: Olympus Manufacturer code: KAT14L(UV083245) LPC0005,</p>	<p>KAT14L(UV083245) UV LAMP K-14 5000 HOURS Original Lamp Lâmpada UV 083245 Lamp UV Lamp UV 083245</p>	
<p>ML1490 Sylvania F14w/133 G13 Coll White 4300K. 5410288000398. Dimensões Prod. Comprimento (mm) 375.4 Diâmetro 26</p>	<p>ML1490 Used for many applications such as frige lighting and under cabinet lighting. As used in various juke boxes. 14w T8 Fluorescent tube. Manufacturer EAN: 5410288000398 F14W/133-T8, Lamp 14W, 220V, G13, 0,395A. 361mm</p>	
<p>GE F17T8/GO/ECO/CVG Fluorescent Clean Room Shatter Resistant Light Bulb Case of 24</p>	<p>General Electric F17T8/GO/ECO/CVG Gold UV Blocking - GE Super Gold CVG lamps, provides maximum UV blockage and outstanding shatter protection which is critical in semiconductor fabs. The exclusive gold sleeve blocks all wave lengths below 520 nm.</p>	
<p>GE F32T8/GO/ECOCVG Fluorescent Shatter- Resistant Gold UV Blocking Light Bulb Case of 36 length 48.000 in (1219.2 mm) Diameter (DIA): 1.000 in (25.4 mm)</p>	<p>General Electric F32T8/GO/ECOCVG GE Super Gold CVG Lamps Provides maximum UV blockage and outstanding shatter protection which is critical in semiconductor fabs. Gold sleeve blocks all wave lengths below 520 nm.</p>	
<p>GE F40/GO/CVG CovRguard Fluorescent T12 Shatter-Resistant Gold UV Blocking Light Bulb Case of 30</p>	<p>General Electric F40/GO/CVG - GE covRguard® T12 - Shatter-Resistant; Gold UV Blocking - GE Super Gold CVG lamps Provides maximum UV blockage and outstanding shatter protection which is critical in semiconductor fabs. The exclusive gold sleeve blocks all wave lengths below 520 nm.</p>	



<p>ML3124 MEDESTIME-Bilicrystal System Duo 2 PLS11w/25000K/2G7 DZS 11/75W, DZS11/75 DZC 11/75W, BC101/8 Only Medical Lamps</p>	<p>ML3124 MEDESTIME PLS11w/25000K/2G7 MEDESTIME- Bilicrystal G11W-4P/BLUE DZS 11/75W DZC 11/75W, BC101/8 for Bilicrystal phototherapy devices</p>	
<p>ML6102 Dulux L Blue 18 W/71 2G11 Osram Dulux L Blue 4Pin. Dulux L 18w/71.</p>	<p>ML6102 Danger 2M21010. RADIUM Ralutec Long 18w/71 2G11. Fototerapia Luz Azul Only Medical Lamps</p>	
<p>ML4494 WEDECO KATADYN 8013927 Lamp. U.V. Wedeco Katdyn K-36L. WEDECO KATADYN UV Lamp, 36W. 8013927. WEDECO KATADYN UV Lamp K-36L Part No.8013927 germicidal lamp</p>	<p>WEDECO KATADYN UV Lamp, 36W.8013927 (Lâmpada para equipamento de esterilização de ortoscópios) K36L</p>	
<p>ML9019 AQUAFINE 3084 ALR-2-3084 C161</p>	<p>AQUAFINE 3084 ALR-23084 C161</p> 	
<p>Sylvania 21625 - F15T8/BLB Fluorescent Tube Black Light</p>	<p>Product Features Color/Finish: Black Light Blue Average Lifetime: 7500 hours F15T8/BLB 15 W MINIATURE FLUORESCENT UV LAMP 18 INCH G13 BASE BLACK LIGHT</p>	
<p>ML7454 G36T5/4P/SE 39W UV Output (watts): 12W Fluorescent Lamp,, G36T5/4P/SE G36T54P/SE, G36T54PSE</p>	<p>Fluorescent Lamp, T5 Lamp Shape, 4-Pin (G10q) Base, Length 34 In., 24 Watts, Lamp Designation G36T5/4P/SE, Clear, UVC source, 5/8 In. Dia., Average Life 9K hr., Cool Lamp Tone, Rapid, Standards Meets 2012 Lighting Legislation Energy Standards. Only shipped in quantities of 24.</p>	
<p>ML4415 Sankyo Denki G8T5E / 8W Lamp UV-B G8T5E / 8W - Sankyo Denki 8W T5 11.3IN GERMICIDAL SANKYO UVB UV-B 306 NM, Length: 290mm, Diameter: 15.5mm, Life: 3000 Hours</p>	<p>ML4415 USHIO G8T5E / 8W Lamp UV-B G8T5E / 8W - USHIO 8W T5 11.3IN GERMICIDAL USHIO UVB UV-B 306 NM, USHIO 3000318 Length: 290mm, Diameter: 15.5mm, Life: 3000 Hours</p>	
<p>ML9006 TUV-6W G6T5 PHILIPS TUV 6W (G6T5 TUV 6W). XG1068/6W, XG1068/42V/6W/G5, TUV-6W-T5- G5-D16-X-226/3 PHILIPS, XG-1068, 34-0013-01-UVP, 4DA16 TUV-6T5-6W-250-PK PHILIPS, LT6W-UV-FILTERED-C NARVA, LT6W/UV-C, T5/GERMICID/6W/G5, 020010, 23375, 7856, 8596.76,</p>	<p>G6T5L2P DONSBULBS, G6T5OF GERMICID/T5/6W/6000H, 6032GL, GERMICIDAL-6-G6T5-1.7-28-PERCENT- 6000-24-15873 GENERAL-ELECTRIC , GL-G6T5 HIKARI, GL/G6T5 HIKARI, GLG6T5 HIKARI, GRM0060 USHIO, LT6W/UV-C, T5/GERMICID/6W/G5 TUV-6/6W/G5/G6T5, 10809, 0002216. 218X16MM</p>	



<p>ML9030 Sankyo Denki 30W G30T8 Ultraviolet 99V UVC T8. 30 watt 99 volt T8 Purification Germicidal UVC Fluorescent Lamp</p>	<p>30 watt 99 volt T8 Purification Germicidal UVC Fluorescent Lamp</p>	
<p>ML9016 Sankyo Denki G16T5 16W Germicidal UV-C Tube G.5 XG1068/6W, 34-0013-01-UVP, XG1068/42V/6W/G5, TUV-6W-T5-G5-D16-X-226/3</p>	<p>16 watt 11" T5 Miniature Bi-Pin (G5) UVC Ultraviolet Clear Linear Straight Fluorescent, G16T5 UV replacement fluorescent light bulbs are germ killing fluorescent UV lamps</p>	
<p>ML7006 Ushio USH-450GS 450 Watt Short Arc Lamp for ASML Stepper 2500/40 with PEP Upgrade, USHIO super high pressure mercury lamps are designed for maximum output in the wavelength area of 436, 405 and 365 nm</p>	<p>ML7006 MPN: USH-450GS ASIN: B0041M70GG, Ushio ush 450-gs USH-450GS/XENON-GAP-SHORT-450W LAMPTECH, USH450GS/XENON-GAP-SHORT-450W, USH450GS LAMPTECH. Ushio USH 450GS Short Arc lamp</p>	
<p>ML7008 Ushio USH450BP, USH 450BP. Bulb GS 450W Super High Pressure UV Lamp USH-450BP Basis Print UV STTER 57 Basis Print UVSTTER57</p>	<p>ML7008 USH-450BP, Ushio's super high-pressure UV lamps have been developed utilizing the three wavelengths of UV rays (436nm, 405nm and 365nm). The lamps are highly illuminant light sources that offer long life and stable radiation intensity</p>	
<p>ML5018 Ushio Q 402 Z7, (Q402 Z7/Q406 Z7) (23030103) multi-spectrum, code: multi-spectrum, diameter: 0.64, graphic arts lamp, mercury vapor, metal halide gas, nm: 300-450. Drager 800 Blue lamp Q402 Z7.</p>	<p>HID/MERCURY-REPROGRAPHIC/450W/MHL450. HID / Mercury Vapor - medium pressure - long arc. volt: 140 amp: 3.21 watt: 450 base: R7s 23030103, 23060101, MHL-450, MHL450/450W/MULTI-SPECTRUM Ushio Q 402 Z7, model 56033065 Spectral Range 300-450 Graphic Arts</p>	
<p>ML7042 ASTRALUX ULTRA-VIOLETAS 220V 150W 220V-150W E27 Blue Glass Carbon filament Coal Filament</p>	<p>ML7041 ASTRALUX INFRA VERMELHOS 220V 150W ASTRALUX INFRA-RED 220V-150W E27 Black Glass Coal Filament, Carbon filament (Filamento de carvão)</p>	
<p>ML9115 Lamp PF712E 230V 15W E27 PHILIPS. PF712E 230V-15W E27 Lâmpada Câmara Escura ML9116 Dr. Fischer 77841422 Photolamp 230v 15w E27</p>	<p>ML9115 Lamp. PF712E, Raio X 230V 15W E27 PHILIPS I/230V/15W/E27/A19/RED-SAFELIGHT/1000H/PF712E ML9116 Dr. Fischer 77841422 Photolamp 230v 15w E27 Dr. Fischer 230v-15w E27</p>	
<p>ML8899 EXZ (Q50MR16C/NFL25° EXZ) 12V 50W (code 20835) EXZ 25°/27°. AIR SHIELDS VIEMERS IIC 5-90/V210726. Light Phototherapy Drager MU03692 Draeger Medical Q50MR16C/NFL25° EXZ 12v 50W Only Medical Lamps</p>	<p>GE #20835 GE EXZ 50 Watt Precise Constant Colour Halogen MR16 Without Cover Glass •GX5.3 Bi-Pin Base • 25° Narrow Flood Beam Spread • 12 Volt • M50 EXZ precise 330709 •Color Temp 3050K •Rated Avg. Life 6000Hrs, MOL 1.875in (47.6 mm) x Diam 2in • Ansi Code EXT. EXZ BLV 189852 12v 50w 24° Light Phototherapy, Light therapy</p>	
<p>ML4577 Ohmeda Biliblanket 20v 140w GX5,3 Lamp. For Phototherapy Device Biliblanket</p>	<p>6600-0680-200 Ohmeda biliblanket 20v 140w GX5,3 Lamp. LAMP OHMEDA 6600-0262-200 (6600-0680-200)(H368)(H 368). Lâmpada OHMEDA 6600-0262-200 (6600-0680-200) (H368)(H 368). Para Aparelho de Fototerapia Biliblanket</p>	